



**ASPIRA-aXA**  
XL-OPTIC - VISION WITHOUT LIMITS

# ASPIRA-aXA

THE EXTENDED 7.0 mm XL-OPTIC



## INNOVATIVE IOL DESIGN

**For standard cataract surgery and especially for patients with large pupils, traumatic mydriasis or iris defects**

- Aspira-aXA as the promising approach to avoid dysphotopsia<sup>1</sup> and deliver a fast adaption to the patient's visual habits
- An investment in the future

**For patients at increased risk of retinal diseases or need for combined vitreoretinal surgery**

- Enlarged panoramic view into the outermost fundus periphery during diagnosis, progression control and therapy<sup>5,8</sup>
- Offers stability even in highly myopic eyes or extreme situations such as intravitreal injections or combined vitrectomies<sup>5,7</sup>



## EASY TO INTEGRATE INTO THE ROUTINE

- Astigmatism-neutral implantation conveniently through small incisions
- Precise and reliable IOL supply by preloaded SAFELOADER<sup>®</sup> autoloading system
- Intuitive, easy handling for a quick and efficient surgical routine

# PREVENT PHOTIC PHENOMENA – FOR AN UNTROUBLED VISUAL OUTCOME

**Pseudophakic dysphotopsia is the most important dissatisfier for patients after successful cataract surgery.<sup>2</sup>**

## BENEFIT FROM A SMOOTH ADAPTION PHASE

- Reduced dysphotopsia from the early postoperative course leads to a fast adaption to the patient's visual habits and results in high patient satisfaction<sup>1</sup>
- Overlap of pupil and IOL optic prevents light passing by the IOL. Incident light is safely directed through the XL-optic
- Interfering edge effects get diminished, frequency and extent of dysphotopsia are reduced in comparison to conventional 6.0 mm IOL designs

## PREVENT PERSISTENT DYSPHOTOPSIA


Bag-to-bag IOL exchange with the 7.0 mm Aspira-aXA showed complete resolution of dysphotopsia in almost all patients. An IOL exchange with a wide optic diameter IOL seems a promising surgical treatment for dysphotopsia.<sup>3</sup>

## SIMULATED BEAM GUIDANCE WITH CONVENTIONAL IOL



Schematic illustration according to Holladay<sup>10</sup>

## SIMULATED BEAM GUIDANCE WITH ASPIRA-aXA



Reduces positive and negative dysphotopsia effectively<sup>1</sup>

### AREA A:


Retinal image is generated by light rays refracted by the IOL optic

### AREA B:

Retinal image is generated by light rays projected next to the IOL optic or refracted at the IOL edge

The dark gap between the two areas A & B can lead to the perception of negative dysphotopsia

# THE PANORAMIC IOL – THINKING ONE STEP AHEAD




## THE SOLUTION FOR VITREORETINAL SURGERY

- 7.0 mm XL-optic ideal for vitreoretinal visualization, facilitation of preventive examinations and follow-up surgeries<sup>4</sup>
- Possibility to perform wide capsulorhexis contributes to undisturbed wide view<sup>5</sup>
- XL diopter range from **-10.0 D** appointed to highly myopic patients

## AN INVESTMENT IN THE FUTURE

- Early detection and management of high-risk groups is key to preventing visual impairment
- Patients suffering from central visual loss that adapt to using unaffected retinal areas for peripheral vision might benefit from the XL panoramic optic.<sup>5,6</sup>

## EXCELLENT INTRAOPERATIVE VIEW



Opening of the anterior capsule membrane with a rhexis diameter of up to 6.5 mm

- ① **Edge of the XL optic**
- ② **Theoretical optic edge of a 6.0 mm IOL**

Extended fundus view for convenient assessment of the tissue structures of the posterior segment of the eye

Facilitates therapeutic measures in the presence of peripheral retinal diseases

Images: Courtesy of Professor M. Bolz, Linz, Austria



*"We have been using Aspira-aXA for several years - especially in combination with vitreous surgery. The particular advantage for me is the excellent stability of the IOL position and the vision, especially when using gas and oil tamponades. Even after almost two hundred implantations, no abnormalities, such as calcification or deposits of the hydrophilic material, occurred."*

Dr. J. Schrecker, Glauchau, Germany

# SAFE POSITIONING – STABLE REFRACTION



## QUATTRO POINT CONTACT ZONE

Large contact area provides stability in the capsular bag




## DOUBLE POSITION MARKING

For correct and safe placement




## CUT-OUT HAPTICS

As absorption elements for embedding the IOL in the capsular bag



- Excellent overall capsular stability<sup>7</sup>
- Even in myopic eyes<sup>5</sup>
- Withstands increased pressure conditions providing stability also under extreme situations such as intravitreal injections or combined vitrectomies<sup>5</sup>



*"Aspira-aXA combines tilt-free fit, excellent A-constant fidelity with optical precision, and its size makes it the number one monofocal IOL for eyes with a WTW greater than 12.0 mm and/or large pupils."*

Dr. E. Becker, Oranienburg, Germany

# CLINICALLY CONVINCING



## STABLE RESULTS<sup>8</sup>

<b>Refraction</b>	Postoperative CDVA stable over time: After 1 month, median 0.00 logMAR; <b>-0.10 to 0.22</b> After 1.5 years, median 0.00 logMAR; <b>-0.10 to 0.10</b>
<b>Mean Tilt</b>	<b>&lt;5.5°</b> , stable in the follow-up period
<b>Mean Decentration</b>	<b>&lt;0.2 mm</b> , stable in the follow-up period
<b>Median Rotation</b>	<b>1.8°</b> one week postoperatively No significant changes between surgery and 1.5-year follow-up

## AN IOL TO WIN YOU OVER

### Effective Reduction of Dysphotopsia

Incidence, frequency and extent of dysphotopsia are lower for patients treated with Aspira-aXA compared to patients treated with a 6.0 mm optic<sup>1</sup>

### Excellent Stability


Confirmed long-term stability with no significant IOL displacement or shift, hence no change in refraction<sup>5, 7, 8, 9</sup>

### Highest Level of Patient Satisfaction

97 % of patients were highly satisfied or satisfied with their vision 1.5 years after surgery<sup>8</sup>

## REFERENCES:

- <sup>1</sup> Bonsemeyer, M.K., et al. (2022). Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. *J Cataract Refract Surg.* 48(1):75-82.
- <sup>2</sup> Welch, N.R., et al. (2010). Satisfaction and dysphotopsia in the pseudophakic patient. *Canadian Journal of Ophthalmology.* 45(2):140-143.
- <sup>3</sup> Rozendal, L.R.W., et al. (2021). Bag-to-bag intraocular lens exchange with a wide optic intraocular lens for treatment of positive and negative dysphotopsia, in 39th Congress of the ESCRS.
- <sup>4</sup> Lee, A.C. and S. Fekrat (2008). Intraocular Lens Implantation From the Vitreoretinal Perspective, in *Retinal Physician.*
- <sup>5</sup> Borkenstein, A.F. and E.M. Borkenstein (2022). Efficacy of Large Optic Intraocular Lenses in Myopic Eyes with Posterior Segment Pathology, in *Ophthalmol Ther.* 11(1):443-452.
- <sup>6</sup> Borkenstein, A.F. and E.M. Borkenstein (2019). Creating Hybrid Monovision with 7.0 mm XL Optic and High-Add AMD Intraocular Lenses (XL-MAGS) in a Patient with Retinitis Pigmentosa, in *Case Rep Ophthalmol.* 10(3):304-311.
- <sup>7</sup> Pilger, D., et al. (2021). Postoperative Lens Rotation of a 7.0 Mm Optic IOL with Plate Haptics, in *J Ophthalmic Vis Res.* 5(1):014-020.
- <sup>8</sup> Schrecker, J., et al. (2022). Performance einer neuen 7-mm-Intraokularlinse mit Nachbeobachtung über 1,5 Jahre [Performance of a new 7 mm intraocular lens with follow-up over 1.5 years], in *Ophthalmologe.* 119(4):367-373.
- <sup>9</sup> Wendelstein, J., et al. (2021). Rotational Stability, Tilt and Decentration of a New IOL with a 7.0 mm Optic, in *Curr Eye Res.* 46(11):1673-1680.
- <sup>10</sup> Holladay, J.T. et al. (2017). Negative dysphotopsia: Causes and rationale for prevention and treatment. *J Cataract Refract Surg.* 43:263-275



## XL DIOPTER RANGE

ASPIRA-aXA/-aXAY PRELOADED IN SAFELOADER®	ASPIRA-aXA IN COMPACT LINE
10.0 to 30.0 in 0.5 D steps	-10.0 to 9.0 in 1.0 D steps 10.0 to 30.0 in 0.5 D steps

## RECOMMENDED IOL CONSTANTS

For laser interferometry and immersion ultrasound biometry

Haigis	Hoffer Q (pACD)	Holladay (surgeon factor)	Holladay 2	SRK/T	Barett LF/DF
$a_0 = -0.521$ $a_1 = 0.269$ $a_2 = 0.215$	5.924	sf = 2.17	5.724	119.568	2.18/-

## PLEASE NOTE

IOL implantation in combination with vitrectomy may influence the predicted lens position. Contact our customer service for details and adjusted power calculation. For the latest constant optimizations, please refer to [www.humanoptics.com](http://www.humanoptics.com) and [IOLcon.org](http://IOLcon.org).



E-MAIL

[customerservice@humanoptics.com](mailto:customerservice@humanoptics.com)

Direct dial to our Customer Service

+49 (0) 9131 50 66 5 44

Follow us on our social media channels!

YouTube LinkedIn

