INTRAOCULAR
IRIS IMPLANT

CUSTOMFLEX® ARTIFICIAL IRIS

de	Gebrauchsanweisung
en	Instructions for use
es	Instrucciones de uso
fr	Mode d'emploi
it	Istruzioni per l'uso
ko	사용설명서
nl	Gebruiksaanwijzing
pl	Instrukcja użytkowania
pt	Instruções de utilização
sl	Navodila za uporabo
sv	Bruksanvisning
tr	Kullanım talimatları
zh	使用说明

V5.1/2025-07_CF_AI_ROW

HumanOptics Holding AG
Spardorfer Str. 150
91054 Erlangen
Germany

Tel.: +49 (0) 9131 50665-0
Fax: +49 (0) 9131 50665-90
mail@humanoptics.com
www.humanoptics.com

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CUSTOMFLEX® ARTIFICIALIRIS

de

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CUSTOMFLEX® ARTIFICIALIRIS

en

These instructions for use are for the following models:

MODEL	DESIGN
CUSTOMFLEX® ARTIFICIALIRIS With Fiber	Silicone elastomer with meshwork
CUSTOMFLEX® ARTIFICIALIRIS Fiber Free	Silicone elastomer without meshwork

1. Description

The CUSTOMFLEX® ARTIFICIALIRIS device is a foldable iris prosthesis that is custom-made for each individual patient and is used for the treatment of iris defects in eyes which are pseudophakic, aphakic or requiring cataract extraction. The CUSTOMFLEX® ARTIFICIALIRIS prosthesis is manufactured from a commercially available ophthalmic silicone and is delivered sterile (sterilized using steam) in a blister filled with isotonic saline solution. Colored silicone paste is applied by hand in a pattern to match the color of the natural iris using a photograph of the existing iris or, in the case of aniridia, the color of the photograph selected by the patient. This custom color-match provides a cosmetically acceptable aesthetic restoration with high patient satisfaction.

The CUSTOMFLEX® ARTIFICIALIRIS device is manufactured as a full 360° iris prosthesis with an overall diameter of 12.8 mm, which can be trephined as needed to custom-fit the device for placement in the posterior chamber (ciliary sulcus or capsular bag).

The device is available in two different variants: With Fiber or Fiber Free. The two models are identical in every respect, except that the With Fiber model has an embedded polyester meshwork layer to provide adequate strength to avoid tearing when suturing. However, the With Fiber model is stiffer and more difficult to fold than the Fiber Free one.


Fig. 1: Anterior view of a CUSTOMFLEX® ARTIFICIALIRIS (left), cross section with dimensions of the models With Fiber (middle) and Fiber Free (right)

For further information regarding the specification of the CUSTOMFLEX® ARTIFICIALIRIS please visit www.humanoptics.com.

2. Mode of action / operating principle

The CUSTOMFLEX® ARTIFICIALIRIS functions as an iris prosthesis. It has a fixed aperture of 3.35 mm, with an opaque perimeter and a black posterior surface to absorb light completely, reducing photic phenomena. The device closely mimics the appearance of the natural iris and at the same time reduces the symptoms associated with aniridia. The small central aperture might increase visual acuity, depth of field, and contrast sensitivity (pinhole effect).

3. Material

The CUSTOMFLEX® ARTIFICIALIRIS is made of a medical-grade hydrophobic silicone elastomer consisting of diphenyl co dimethyl polysiloxane and silicone reinforcing resin. For colorization inorganic pigments are used. The model With Fiber is addition-

ally reinforced by an embedded polymer fiber meshwork consisting of polyethylene terephthalate (PET).

4. Intended use

4a. Indications

The CUSTOMFLEX® ARTIFICIALIRIS is indicated for the treatment of full or partial aniridia resulting from congenital aniridia, acquired defects, or other conditions associated with full or partial aniridia.

4b. Intended purpose

The CUSTOMFLEX® ARTIFICIALIRIS is intended for use as an iris prosthesis for the treatment of iris defects in eyes which are pseudophakic, aphakic or requiring cataract extraction. The device is intended for implantation in the posterior chamber (ciliary sulcus or capsular bag).

4c. Intended patient population

The CUSTOMFLEX® ARTIFICIALIRIS is indicated for use in adults and children for the treatment of full or partial aniridia resulting from congenital aniridia, acquired defects, or other conditions associated with full or partial aniridia.

4d. Intended users

The CUSTOMFLEX® ARTIFICIALIRIS must be handled by health professionals and implanted by specially trained and certified ophthalmic surgeons (see point 20. Further requirements for usage).

4e. Contraindications

The CUSTOMFLEX® ARTIFICIALIRIS device is contraindicated in eyes with any of the following conditions:

- Uncontrolled ocular inflammation (e.g. uveitis)
- Severe chronic uveitis
- Microphthalmus
- Untreated retinal detachment
- Untreated chronic glaucoma
- Rubella cataract
- Rubeosis of the iris
- Proliferative diabetic retinopathy
- Stargardt's retinopathy
- Pregnancy
- Intraocular infections

5. Caution

Implantation of the CUSTOMFLEX® ARTIFICIALIRIS is not recommended in patients with the following conditions and situations:

- Children who are less than 6 years of age because their eyes are still in a stage of major growth development that would be disrupted by ocular surgery
- Preoperative intraocular pressure (IOP) above 21 mmHg that does not respond to pressure-lowering medication, unless the IOP above 21 mmHg is due to a known underlying condition that is wellcontrolled with glaucoma treatment, such as ocular hypertension or open-angle glaucoma
- Patients with severe endothelial corneal dystrophy, because the surgical procedure to implant the CUSTOMFLEX® ARTIFICIALIRIS may damage the cornea sufficiently so that the potential benefits of implantation do not outweigh the risks
- The CUSTOMFLEX® ARTIFICIALIRIS is not designed for placement in the anterior chamber

- No useful vision or visual potential in the fellow eye, unless the patient has debilitating visual symptoms, so that the potential benefits of CUSTOMFLEX® ARTIFICIAL/IRIS implantation clearly outweigh the risks
- Presence of a condition or finding in the fellow eye that would make it unsafe to implant a CUSTOMFLEX® ARTIFICIAL/IRIS prosthesis in the eye to be treated
- Allergy to any of the planned postoperative antibiotic or anti-inflammatory medications, unless a suitable alternative medication can be prescribed
- Implantation for cosmetic color changes of the iris
- Post-partum women who are nursing or lactating and for whom postoperative medications are contraindicated
- Patients with gastric ulcers or diabetes mellitus in whom high doses of orally administered systemic steroids are required postoperatively
- Any other condition that would interfere with the planned surgical procedure to implant the iris device

6. Warnings

The CUSTOMFLEX® ARTIFICIAL/IRIS should be used with caution in the following situations:

- A clear natural crystalline lens – Do not implant the CUSTOMFLEX® ARTIFICIAL/IRIS in phakic eyes
- The visual potential of the fellow eye cannot be evaluated preoperatively (e.g., poor visual acuity due to cataract)
- Preoperative IOP > 21 mmHg that is known to be stable and well-controlled with glaucoma treatment (e.g., medication, tubes or shunts)
- Presence of any other medical condition that might be expected to make the patient an unsuitable candidate for CUSTOMFLEX® ARTIFICIAL/IRIS implantation
- Anticipated complexity of the planned surgical procedure that might increase the potential for complications
- Implantation in the fellow eye before stabilization of the first implanted eye (typically 1 month or more)

The aperture of the CUSTOMFLEX® ARTIFICIAL/IRIS pupil is fixed at 3.35 mm. In the event that a larger pupillary opening is required for posterior segment surgery, the CUSTOMFLEX® ARTIFICIAL/IRIS can be explanted and a new CUSTOMFLEX® ARTIFICIAL/IRIS may be implanted after completion of the posterior segment surgery. The technique for a secondary implantation would be determined in the same manner as for a primary implantation.

7. Potential complications and undesirable side effects

Below is a list of the potential complications and undesirable side effects associated with the use of the device, surgical procedure, or IOL.

Device-related complications associated with the CUSTOMFLEX® ARTIFICIAL/IRIS may include but may not be limited to:

- Worsening of photosensitivity and vision
- Elevated intraocular pressure
- Decrease in uncorrected distance visual acuity
- Decrease in best-corrected distance visual acuity
- Eye infection/inflammation
- Incorrect device positioning, dislocation, and decentration
- Secondary (additional) surgical intervention

Surgical repositioning, replacement or removal of the device may be necessary to correct device dislocations. Device defects can occur if the device is not handled properly.

Surgery-related adverse events may include, but may not be limited to:

- Cystoid macular edema
- Hypopyon
- Endophthalmitis
- Device migration
- Pupillary block
- Retinal detachment
- Secondary surgical intervention (unplanned)
- Corneal edema, persistent at 3 months or later
- Chronic iritis/anterior segment inflammation persistent at 3 months or later

8. Clinical benefits

The CUSTOMFLEX® ARTIFICIAL/IRIS implant has the potential to improve the quality of life by reducing visual symptoms and improving the cosmetic appearance of the eye.

9. Safety and clinical performance

For products registered under Regulation (EU) 2017/745, the summary of safety and clinical performance (SSCP) will be published in EUDAMED, the European Database on Medical Devices, under the URL (<https://ec.europa.eu/tools/eudamed>).

10. Handling

Before use, check the package for the correct model and the expiration date. The implant should not be implanted after the indicated expiration date.

Before use, check the integrity of the sterile barrier system. The device is sterile only if the sterile pouch is undamaged. The implant container may only be opened under sterile conditions.

To remove the CUSTOMFLEX® ARTIFICIAL/IRIS, hold the flap of the sealed foil lid of the container and pull it off. After removing the implant from the container, ensure that the device surface is free of any adhering particles or any other defects.

Please note that the color of the CUSTOMFLEX® ARTIFICIAL/IRIS in air versus in saline solution differs. The actual color of the device in aqueous humor may vary due to the cornea.

In the event of a malfunction of the device or a change in its performance, please return the affected product including all available documentation (e.g. labels, packaging) to your local distributor or the manufacturer. Contact to the manufacturer via email: complaint@humanoptics.com.

11. Surgical use instructions

11a. Model selection and preparation

The CUSTOMFLEX® ARTIFICIAL/IRIS is a foldable device that is custom-made for each individual patient. The CUSTOMFLEX® ARTIFICIAL/IRIS is available in two models: With Fiber or Fiber Free. The two models are identical in every respect, except that the With Fiber model has an embedded polyester meshwork layer to provide adequate strength to avoid tearing when suturing. The selection of the surgical technique should be dictated by the preoperative iris and anterior segment anatomy and pathology. The With Fiber model is generally used when suture fixation is planned, whereas the Fiber Free model is designed for sutureless implantation techniques. For placement in the capsular bag, only the Fiber Free version should be used. See section "General Surgical Procedure" for further considerations and a description of implantation techniques.

Important considerations for preparing the device:

- **Capsular Bag Trephining:** When implantation in the capsular bag is planned, the CUSTOMFLEX® ARTIFICIAL IRLS should be trephined to an appropriate size. For an adult eye with an average-sized natural lens, the appropriate diameter is typically 10.0 mm, though it can vary from patient to patient. The capsular bag diameter should be estimated based on the size of the evacuated capsular bag once a capsular tension ring has been placed, especially in smaller eyes, pediatric eyes, larger myopic eyes, or megalophthalmic eyes.
- **Ciliary Sulcus Trephining:** For passive placement into an anatomically suitable ciliary sulcus, the iris device should be trephined to the estimated smallest diameter of the ciliary sulcus.
- Always use sharp, sterile instruments for cutting or trephining the CUSTOMFLEX® ARTIFICIAL IRLS.
- The CUSTOMFLEX® ARTIFICIAL IRLS can be delivered via forceps or injector. Suitable injection systems may be inquired about or purchased from the manufacturer.
- When suturing the CUSTOMFLEX® ARTIFICIAL IRLS With Fiber model, the suture pass should be at least 1 mm away from the edge in order to guarantee the stability of the device after suturing and thereby minimizing the risk of device decentration or dislocation.

12. General surgical procedure

12a. Anterior segment preparation

The anterior segment should be appropriately prepared by cataract removal and IOL placement, and/or vitrectomy, as dictated by the preoperative anterior segment anatomy and pathology, in preparation for implanting the device by one of the surgical methods described below. The limbal-corneal wound should be of adequate size for the selected delivery method. Typically a 2.75 mm wound is required to insert the device, and at least a 4 mm wound is required if forceps will be used to insert the device.

12b. Capsular bag placement

The Fiber Free model should be used for placement of the iris device within the capsular bag. The anterior segment should be appropriately prepared, as described in the "Anterior Segment Preparation" section above. The anterior capsule should be stained with trypan blue or indocyanine green at the beginning of the surgical procedure. A capsular tension ring is inserted into the capsular bag to prevent post-operative retraction of the capsule with subsequent tilt and decentration of the artificial iris device. The selected IOL is then implanted in the capsular bag. The limbal-corneal wound should be enlarged, if necessary, for implantation of the iris device. During the procedure, the anterior chamber should be deepened as much as possible by a cohesive ophthalmic viscosurgical device (OVD) to allow adequate space for the iris device to unfold, minimizing contact with other intraocular structures. If the initial dye has faded, additional trypan blue or indocyanine green may be applied or instilled along the anterior capsule margin just prior to iris device implantation.

The iris device is either folded for implantation with forceps or rolled and placed in the injection cartridge with the colored side outward. The leading edge of the folded device should be placed under the distal capsule margin, visualized by noting the trypan blue or indocyanine green over the iris device, before the iris device is unfolded. It should be allowed to unfold with the edges of the implant oriented posteriorly (curled with the colored side outward), so that contact with the corneal endothelium is minimized. A spatula can guide the unfolding or injecting process. Once the iris device is unfolded, the edges can be completely tucked into the capsular bag, with care being taken to avoid undue pressure on the bag margins, especially in patients with congenital aniridia. The iris device can be manipulated either by hooks or micrograsping forceps to facilitate

positioning. If the iris device does not go into the bag easily, grasping it with intraocular microforceps at the pseudopupil margin and folding it can facilitate implantation. Once the IOL and iris device are centered and stability is confirmed, the OVD can be removed. If the chamber shallows, the iris device may escape from the capsular bag and require repositioning. Removal of the OVD using a bimanual approach may help to maintain a deep chamber and avoid dislocation of the artificial iris. The incision should be sealed and secured according to surgeon's preference. Instillation of intraocular carbachol is recommended to reduce the risk of post-operative pressure elevation.

12c. Passive sulcus placement

Either the With Fiber or Fiber Free model of the iris device can be used for implantation in the ciliary sulcus without suture fixation. If suture fixation is necessary, the With Fiber model should be used. The anterior segment should be appropriately prepared, as described in the "Anterior Segment Preparation" section above. The iris device should be trephined to the estimated sulcus size, as measured preoperatively by ultrasound or intraoperatively by direct measurement of the pressurized globe. The limbal-corneal incision should be of adequate size. The anterior chamber should be deepened as much as possible with a cohesive OVD to allow adequate space for the iris device to unfold, minimizing contact with intraocular structures.

The iris device should be folded for implantation with forceps or rolled and placed in the injection cartridge with the colored side outward. The leading edge of the folded device should be placed in the ciliary sulcus and allowed to unfold with the edges of the implant oriented posteriorly, so that contact with the corneal endothelium is minimized. The iris device can be manipulated either with hooks or with a micro-grasping small-gauge intraocular forceps to facilitate the positioning. A snag fit should be confirmed. If the iris device appears to buckle or fits too tightly, it should be removed, trephined to a smaller size, and then reinserted. If the iris device is freely mobile in the sulcus due to the device being trephined too small, it can be removed and replaced with the standby device after it is trephined to a larger diameter. Alternatively, gently placed and carefully tightened suspension sutures can be placed through the scleral wall at the ciliary sulcus to prevent movement of the device. The sutures should be tied with only enough tension to prevent movement and achieve centration. Overtightening the sutures can tear the device if a Fiber Free device is utilized. Once acceptable centration and stability are confirmed, the OVD can be removed. Removal of the OVD using a bimanual approach may help to maintain a deep chamber and avoid dislocation of the artificial iris. The incision should be sealed and secured according to surgeon preference. Instillation of intraocular carbachol is advised to reduce the risk of postoperative pressure elevation. Patch graft material can be placed over fixation sutures, as deemed necessary by the operating surgeon.

12d. Sulcus placement with suture fixation to the scleral wall

The With Fiber model should be used for fixation in the ciliary sulcus with sutures. The anterior segment should be appropriately prepared, as described in the "Anterior Segment Preparation" section above. The iris device should be trephined to at least 1 mm less than the estimated sulcus size, as measured preoperatively by ultrasound or intraoperatively by direct measurement of the pressurized globe.

The limbal-corneal incision should be of adequate size. The anterior chamber should be deepened as much as possible by a cohesive OVD to allow adequate space for the iris device to unfold, minimizing contact with intraocular structures.

After the preparation ex-vivo of the sutures for the fixation to the scleral wall, the iris device should be folded for implantation with forceps with the colored side outward. The leading edge of the folded device should be placed in the ciliary sulcus and allowed to unfold with the edges of the implant oriented posteriorly, so that contact with the corneal endothelium is minimized. The iris device can be manipulated either with hooks

or with a micro-grasping small-gauge intraocular forceps to facilitate the positioning. The iris device should be placed within the ciliary sulcus, and an adequately snug fit should be confirmed. If the iris device appears to buckle or fits too tightly, it should be removed, trephined to a smaller size, and then reinserted. If the iris device is freely mobile, then the sutures should be passed and tightened to achieve good centration of the device. Overtightening of the sutures could lead to ovalization of the pupil, distortion of the device. Once acceptable centration and stability are confirmed, the OVD can be removed. Removal of the OVD using a bimanual approach may help to maintain a deep chamber and avoid dislocation of the artificial iris. The incision should be sealed and secured according to surgeon preference. Instillation of intraocular carbachol is advised to reduce the risk of postoperative pressure elevation. Patch graft material can be placed over fixation sutures, as deemed necessary by the operating surgeon.

12e. Sulcus placement of the iris device and Posterior-chamber IOL (PCIOR) with suture fixation to the scleral wall

Only the CUSTOMFLEX® ARTIFICIALIRIS With Fiber model can be used for this surgical technique.

The suture fixation of both a PCIOR and iris device can be achieved using one of three methods:

- 1) Fixation of the iris device to the PCIOR ex vivo on the surgical field, then the PCIOR-iris device complex is affixed by using non-absorbable sutures passing through the scleral wall with these sutures affixed to the IOL portion of the complex;
- 2) Fixation of the iris device to the PCIOR ex vivo on the surgical field, then the PCIOR-iris device complex is affixed by using non-absorbable sutures passing through the scleral wall with these sutures affixed to the iris device portion of the complex; or,
- 3) The PCIOR and iris device can be independently fixated to the scleral wall using non-absorbable sutures, either placed through same scleral wall openings or separate scleral wall openings.

Note: Gluing the CUSTOMFLEX® ARTIFICIALIRIS to the IOL is not a recommended method to achieve fixation.

13. MRI safety status

The CUSTOMFLEX® ARTIFICIALIRIS devices containing the maximal amount of magnetic pigments were tested according to ASTM Standards F2052-15:2015, F2119-07:2013, F2182-11a:2011 and F2213-17:2017.

Non-clinical testing has demonstrated the CUSTOMFLEX® ARTIFICIALIRIS is MR Conditional. A patient with this device can be safely scanned in an MR system meeting the following conditions:

- Static magnetic field of 1.5 Tesla, 3 Tesla and 7 Tesla
- Maximum spatial field gradient of 200 G/cm (2 T/m)
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 2 W/kg (normal Operating Mode)

Care has to be taken when moving the patient inside and outside the imaging area, because spatial field gradients for almost all scanners are higher than the limits defined here. However, these higher spatial field gradients could affect the patient only for a very short period of time (a few seconds) and for a very short distance (within a few centimeters).

Under the scan conditions defined above, the CUSTOMFLEX® ARTIFICIALIRIS device is expected to produce a maximum temperature rise of less than 2.0°C after 15 minutes of continuous scanning.

In non-clinical testing, the image artifact caused by the device extends approximately 24.2 mm from the CUSTOMFLEX® ARTIFICIALIRIS when imaged with a gradient echo pulse sequence and a 7 Tesla MR system.

If the MRI scan can only be performed under different conditions than the above-mentioned ones, it is recommended that the patient sees the treating ophthalmologist after the examination.

This information is also available on the website www.humanoptics.com/mri.

14. Reprocessing

The CUSTOMFLEX® ARTIFICIALIRIS is for single-use only. Reprocessing or re-sterilization of the CUSTOMFLEX® ARTIFICIALIRIS is strictly prohibited, and may compromise device performance, which could cause serious harm to the patient's health and safety.

15. Disposal in accordance with national and local regulations

Discarded CUSTOMFLEX® ARTIFICIALIRIS devices (used or unused) are classified as medical or clinical waste due to their potentially infectious nature and must be disposed of accordingly in accordance with national and local regulations.

16. Patient information

A patient card is included in the package of every product. Enter the patient data on the patient card and apply the self-adhesive label containing the product identification to the designated space on the card. Instruct the patient to keep this card as a permanent record and to show it to any eye care professional consulted in the future. For further patient information please visit www.humanoptics.com/patient-information.

17. Lifetime of the CUSTOMFLEX® ARTIFICIALIRIS

Once implanted, the CUSTOMFLEX® ARTIFICIALIRIS implant remains permanently in the patient's eye. Due to the mechanical, optical and biocompatible material properties of the CUSTOMFLEX® ARTIFICIALIRIS, the device is expected to be stable indefinitely over the lifetime of the patient.

18. Reporting

Serious incidents should be reported to HumanOptics and to the relevant competent authorities.

19. Disclaimer

The manufacturer is not liable either for the implantation method or the operative technique used by the physician performing the procedure or for the selection of the CUSTOMFLEX® ARTIFICIALIRIS in relation to the patient or his/her condition.

Furthermore, the manufacturer is not liable for a postoperative difference in color between the natural iris tissue and the iris implant.

20. Further requirements for usage

A high level of surgical skill and experience in the field of anterior segment surgery is required for CUSTOMFLEX® ARTIFICIALIRIS implantation. Before the first implantation, the surgeon must have successfully completed the OCC (Online Certification Course) by HumanOptics. Each participant receives a certificate number, which is required to proceed with the product order.

EU law restricts this device to practitioners who have been trained and have experience in the surgical management and treatment of aniridia and iris defects.

The CUSTOMFLEX® ARTIFICIALIRIS is restricted to sale by or on the order of a physician or any other health entity.

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