INTRAOCULAR LENS ACRYL

monofocal / multifocal / toric

| | |
|-----------|---------------------------|
| de | Gebrauchsanweisung |
| en | Instructions for Use |
| bg | Инструкции за употреба |
| cs | Návod k použití |
| el | Οδηγίες χρήσης |
| es | Instrucciones de uso |
| fr | Mode d'emploi |
| hu | Használati útmutató |
| id | Petunjuk Penggunaan |
| it | Istruzioni per l'uso |
| ko | 사용설명서 |
| lt | Naudojimo instrukcija |
| nl | Gebbruiksaanwijzing |
| no | Bruksanvisning |
| pl | Instrukcja użycia |
| pt | Instruções de utilização |
| ro | Instrucțiuni de utilizare |
| ru | Инструкция по применению |
| tr | Kullanım Talimatları |
| vi | Hướng dẫn sử dụng |

MODEL

AS
MC 611 MI
Aspira-aA
Aspira-aAY
Aspira-aQA
Aspira-aXA
Aspira-aXAY
MC X11 ASP
Aspira3P-aVA
Torica-aA
Torica-aAY
Diff-aA
Diff-aAY
Triva-aA
Triva-aAY
Triva-aXAY
ToricaDiff-aA
ToricaDiff-aAY

HumanOptics Holding AG
Spardorfer Str. 150
91054 Erlangen
Germany

V10.0/2021-07

Tel.: +49 (0) 9131 50665-0
Fax: +49 (0) 9131 50665-90
mail@humanoptics.com
www.humanoptics.com

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Acryl-Intraokularlinse

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Instructions for Use

Acrylic Intraocular Lens

These instructions for use are for the following intraocular lenses (abbreviated "IOL"):

| MODEL | DESIGN |
|--|---|
| Monofocal | |
| AS / MC 611 MI | One-piece, spherical |
| Aspira-aA / -aAY / -aQA / -aXA / -aXAY | One-piece, aspheric aberration-free, 360° lens epithelial cell barrier on the posterior edge |
| MCX11 ASP | One-piece, aspheric aberration-correcting |
| Aspira3P-aVA | Three-piece, aspheric aberration-free |
| Monofocal toric | |
| Torica-aA / -aAY | One-piece, aspheric aberration-free, toric, 360° lens epithelial cell barrier on the posterior optic edge |
| Multifocal | |
| Diff-aA / -aAY | One-piece, aspheric aberration-free, multifocal (bifocal) diffractive, 360° lens epithelial cell barrier on the posterior optic edge |
| Triva-aA / -aAY / -aXAY | One-piece, aspheric aberration-free, multifocal (trifocal) diffractive, 360° lens epithelial cell barrier on the posterior optic edge |
| Multifocal toric | |
| ToricaDiff-aA / -aAY | One-piece, aspheric aberration-free, toric, multifocal (bifocal) diffractive, 360° lens epithelial cell barrier on the posterior optic edge |

1. Description

A sterile, foldable, UV-absorbing acrylic posterior chamber IOL in isotonic saline solution for implantation into the capsular bag after phacoemulsification. The models labelled with "Y" additionally contain a blue-light filter.

For further information regarding the specification of the IOLs listed in the table above, please visit www.humanoptics.com.

There are two different packaging versions. Not all products are available in both versions:

- Compact line: The IOL is packaged in a flat container for manual loading into a conventional cartridge injector.
- Safeloader™: The Safeloader® Autoloading System is composed of an ACCUJECT™ injector with integrated cartridge (manufacturer: Medical AG) and an autoloading container with a preloaded acrylic IOL.

A list of suitable injection systems can be found under www.humanoptics.com.


Fig. 1: Transmittance spectra of the hydrophilic acrylic IOLs

10% cut-off wavelength:

Curve (1): the cut-off wavelength at 10% transmission is 375 nm

Curve (2): the cut-off wavelength at 10% transmission is 380 nm

* Source: Boettner EA, Wolter JR. Transmission of Ocular Media, Investigative Ophthalmology, 1962; 1:776–783

2. Material

One-piece acrylic IOLs:

Hydrophilic acrylic copolymer consisting of 2-hydroxyethyl methacrylate (HEMA) and methyl methacrylate (MMA) with a UV absorber polymerized into the polymer matrix. Water content 26%, refractive index 1.46, Nd:YAG laser compatible.

Models followed by "Y" have a yellow coloring agent polymerized into the polymer matrix to absorb the high-energy portion of the blue light.

Three-piece acrylic IOLs:

Optic material: Refer to "One-piece acrylic IOLs". Nd:YAG laser compatible.

Haptic material: Highly flexible, blue colored polyethersulfone.

3. Intended purpose

3a. Indications

All IOLs covered by these instructions for use are indicated for the correction of aphakia after surgical removal of the natural lens.

In addition, toric IOLs are indicated to correct pre-existing regular, stable corneal astigmatism.

Multifocal IOLs are indicated for patients who desire near, and/or intermediate, and distance vision with increased spectacle independence. However, accommodation will not be restored. Toric-multifocal IOLs can further correct pre-existing regular, stable corneal astigmatism.

3b. Intended use

All these IOLs are intended for implantation in the capsular bag to replace the natural lens.

3c. Target patient group

Aphakic adult patients.

3d. Target user group

IOLs must be handled by health professionals and implanted by ophthalmic surgeons.

3e. Contraindications

Outside of general contraindications of ocular surgery, there are no specific contraindications for the IOLs covered by these instructions.

4. Caution

Careful preoperative evaluation and sound clinical judgment should be used by the surgeon to decide the benefit/risk ratio before implanting a lens in a patient with one or more of these conditions:

- Uveitis
- Proliferative diabetic retinopathy
- Uncontrolled chronic glaucoma
- Corneal endothelial dystrophy
- Microphthalmos or macrophthalmos
- Suspected ocular infection
- Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule adhesion that may impact negatively visual acuity
- Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules)
- Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss)
- Circumstances that would result in damaging the endothelium during implantation

For multifocal IOLs, caution should also be exercised in the following conditions:

- Amblyopia
- Optic nerve atrophy
- Significant corneal aberrations
- Medical conditions that may impair vision
- Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction

5. Warnings

- Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can theoretically be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL.
- Multifocal effects can cause changes on the surface of the IOL postoperatively (e.g. elevated concentration of calcium and/or phosphate in the aqueous humour, caused by dysfunction of the blood-aqueous humour barrier).

6. Special considerations prior to multifocal IOL implantation

- It is recommended to target emmetropia.
- Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes.
- In cases of corneal astigmatism > 1.0 D, implantation of a toric-multifocal IOL is recommended.
- Lens tilt and decentration may negatively affect the quality of vision.
- As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in poor light conditions.

• Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions.

7. Special considerations prior to toric IOL implantation

- The dioptric power of toric IOLs can be either indicated as "sphere (SPH) and cylinder (CYL)" or "spherical equivalent (SE) and cylinder (CYL)". Please check the label carefully.
- Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism.
- Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
- Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation.
- If possible, posterior corneal astigmatism (tomography) should be considered, especially with astigmatism < 2.0 D.

8. Potential complications

- As with any surgical procedure, cataract surgery with IOL implantation presents risks, which your surgeon must evaluate. Potential complications may include but are not limited to:
 - Corneal endothelial decompensation, corneal edema, retinal detachment, cystoid macular edema, increased intraocular pressure, inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis), iris trauma, posterior capsular and zonular rupture, secondary surgical intervention (e.g. repositioning, removal, or exchange).

9. Calculation of dioptric power

Accurate biometry is essential to successful visual outcomes. Preoperative calculation of required lens power for the IOL should be determined by the surgeon's experience, preference, and intended lens placement. In this context, the incision location and the surgeon's estimated surgically induced corneal astigmatism should be considered, particularly in the case of toric IOLs. As a starting point for IOL power calculations, please use the constants for the respective formula and product as recommended by the manufacturer (www.humanoptics.com). Lens constants must be "personalized" to address differences in instrumentation, measurement and surgical techniques, and IOL power calculation methods.

Physicians requiring additional information for dioptric power calculation should contact the manufacturer (application@humanoptics.com).

9a. Orientation of toric IOLs

The axis of the plus cylinder (meridian of lowest dioptric power, r1) is indicated by two opposite indentations (markings) at the edge of the optic (Fig. 2).


Fig. 2: Markings on toric intraocular lenses, C-loop design

Correct the corneal astigmatism by aligning the markings of the IOL with the steep axis of the cornea.

10. Handling

- Store the lens between 10°C and 30°C.
- Before use, check the lens package for the correct lens model, dioptric power and expiration date. The lens should not be implanted after the indicated expiration date.
- Before use, check the integrity of the sterile barrier system. The IOL is sterile only if the sterile pouch is undamaged. The lens container may only be opened under sterile conditions.
- Before use, the IOL should be warmed to temperatures between 18°C (operating room) and 36°C (intraocular temperature) to avoid any risk of damage to the IOL during implantation.
- To remove the IOL, hold the flap of the sealed foil lid of the container and pull it off. After removing the IOL from the container, ensure that the IOL surface is free of any adhering particles or any other defects.
- For Safeloder® products, the integrity of the IOL must be checked after successful loading of the loading chamber.

Attention: The IOL should not dehydrate! Hydrophilic acrylic IOLs may only be wetted with sterile isotonic saline solution.

In the event of a malfunction of the device or a change in its performance, please return the affected product including all available documentation (e.g. labels, packaging) to your local distributor or the manufacturer.

11. Preparation of the patient before implantation of toric IOLs

If manual marking is performed, please consider the following steps:

- With the patient sitting upright, mark the horizontal axis (0°) or the vertical axis (90°) on the cornea, as the reference axis. Positioning the patient in an upright sitting position is important to prevent ocular cyclotorsion.
- Then, mark the steepest axis of the cornea using the reference axis (0° or 90°).

12. Implantation

The capsulorhexis size should be about 0.5 mm smaller than the optic diameter of the IOL.

To ensure smooth and safe IOL implantation, fill the anterior chamber and capsular bag with sufficient viscoelastic material.

During implanting, ensure correct anterior/posterior orientation of the IOL: the extremities of the C-loop haptics point in a counterclockwise direction (anterior view). The models with other haptic geometries are correctly oriented when one marking appears on the top right and the other marking on the bottom left (see graphics).

Implantation of the foldable acrylic IOLs can be performed using either forceps or an injection system. A list of suitable systems can be found under www.humanoptics.com.

When using Safeloder® products, please refer to the enclosed Safeloder®-Instruction for use.

When using an injector system for implantation, please refer to the specific instructions for use provided with the injector system to ensure that it is suitable to use with the IOL to be implanted.

Special recommendations for toric IOLs:

- In order to achieve the intended position, the IOL can be rotated (clockwise in the case of the C-loop) with a push-pull hook placed at the optic-haptic junction. It might be useful to position the IOL 10° to 20° shy of the desired position, remove the viscoelastic material, then rotate the IOL to its final position.
- At the end of surgery, it is important to completely remove all viscoelastic material from behind the implant.
- After removing the viscoelastic material, recheck the correct IOL positioning.
- Patients should be kept at rest immediately after surgery in order to keep the IOL stable.

13. MR compatibility

The implant is MR compatible, and in MRI examinations it does not display any increase in temperature, image artifacts or changes in position. All tests to investigate MR safety were performed at 7 Tesla.

14. Reprocessing

Reprocessing of the implant is strictly prohibited, since material changes, for example, may cause serious complications and can be fatal.

15. Disposal in accordance with national and local regulations

Discarded IOLs (used or unused) are classified as medical or clinical waste due to their potentially infectious nature and must be disposed of in accordance with national and local regulations.

16. Patient information

A patient card is included in the package of every product. Enter the patient and clinic data on the patient card and apply the self-adhesive label containing the identification information of the product to the designated space on the card. Instruct the patient to keep this card as a permanent record of the implant and to show it to any eye care professional consulted in the future. The patient card also refers to the HumanOptics website for further patient information.

17. Reporting

Serious incidents should be reported to HumanOptics and to the relevant competent authorities.

18. Disclaimer

The manufacturer is not liable for the implantation method or the operative technique used by the physician performing the procedure or for the selection of the IOL in relation to the patient or his condition.

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