Summary of safety and clinical performance (SSCP) for the Artificial Iris

SSCP reference number: P036

Revision history

Rev number	Revision date	Description of change	Revision validated by the Notified Body
Rev.01	2025-03-11	Initial creation	☒ Yes Validation language: English ☐ No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB)
Rev.02	2025-08-01	Adaptation according to Change C-2025-008 based on decision by case conference TÜV NORD (see e-mail 2025-07-07): - Vulnerable populations: part 1 chapter 2 and part 2 chapter 2 amended - Product lifetime: part 1 chapter 5.4 and part 2 chapter 5.1 amended	☒ Yes Validation language: English ☐ No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB)

Table of contents

	Page
Part 1 – SSCP intended for healthcare professionals	3
1. Device identification and general information	3
2. Intended use of the device	4
3. Device description	5
4. Risks and warnings	6
5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)	9
5.1. Summary of clinical data related to equivalent device	9
5.2. Summary of clinical data from conducted investigations of the device before the CE-marking	9
5.3. Summary of clinical data from other sources	9
5.4. Overall summary of the clinical performance and safety	11
5.5. Ongoing or planned post-market clinical follow-up	16
6. Possible diagnostic or therapeutic alternatives	16
7. Suggested profile and training for users	17
8. Reference to any harmonised standards and CS applied	17
9. References	18
Part 2 – SSCP intended for patients	20
1. Device identification and general information	20
2. Intended use of the device	20
3. Device description	21
4. Risks and warnings	24
5. Summary of clinical evaluation and post-market clinical follow-up	25
5.1. Clinical background of the device	25
5.2. The clinical evidence for the CE-marking	25
5.3. Safety	26
6. Possible diagnostic or therapeutic alternatives	26
7. Suggested training for users	26

Part 1 – SSCP intended for healthcare professionals

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The SSCP is not intended to replace the Instructions For Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients. This SSCP is reviewed at least annually and will be updated if needed to ensure that any clinical and/or safety information in the SSCP remains correct and complete.

The following information is intended for users/healthcare professionals. Following this information there is a summary intended for patients.

1. Device identification and general information

Manufacturer's name and address	HumanOptics Holding AG Spardorfer Straße 150 91054 Erlangen Germany
Manufacturer's single registration number (SRN)	DE-MF-000017892
Authorised representative if applicable; name and SRN	N/A
Notified body name and organization number	TÜV NORD CERT GmbH NB 0044


Device trade name(s)	Artificial Iris • Variant Artificial Iris with Fiber (wF) • Variant Artificial Iris Fiber Free (FF) Different color categories ("CustomFlex®", "Select" and "Standard") are available ¹ for both the wF and FF variants.
Basic UDI-DI(s)	04049154_AI_M1_H1_O1_X7
CND/EMDN code(s):	P0399 (Y063021, PROSTHESES, OCULISTIC – OTHERS (OCULAR PROSTHESES))
Class of device	☐ I ☐ IIa ☒ IIb ☐ III
Year when the first certificate (CE) was issued covering the device	2011

¹ Not all color categories are available for sale in all countries.

2. Intended use of the device

Intended Use	
Intended purpose	The device is intended for use as an iris prosthesis for the treatment of iris defects in eyes which are pseudophakic, aphakic or require cataract extraction. The device is intended for implantation in the posterior chamber (Fiber Free model: ciliary sulcus or capsular bag; with Fiber model: ciliary sulcus).
Intended users	The Artificial Iris must be handled by health professionals and implanted by specially trained and certified ophthalmic surgeons.
Use environment	Ophthalmic operating room
Indications	The device is indicated for the treatment of full or partial aniridia resulting from congenital aniridia, acquired defects, or other conditions associated with full or partial aniridia.
Medical condition(s) to be treated	Indicated for the treatment of full or partial aniridia resulting from congenital, aniridia, acquired defects, or other conditions associated with full or partial aniridia. Aniridia is defined as the complete or partial absence of iris tissue and includes iris defects like missing iris segments, absent or deficient iris pigment epithelium, or complete iris absence. Intended purpose includes iris defects in eyes which are pseudophakic, aphakic or requiring cataract extraction.
Intended patient population	The device is indicated for use in adults and children from 6 years for the treatment of full or partial aniridia resulting from congenital aniridia, acquired defects, or other conditions associated with full or partial aniridia. The manufacturer does not hold clinical data relating to breastfeeding women or immunocompromised patient populations.
Intended anatomical location	Ciliary sulcus ("Fiber Free" model or "with Fiber" model) or capsular bag ("Fiber Free" model only) of the human eye
Contraindications	The Artificial Iris is not to be used for purely cosmetic reasons. The device is contraindicated in eyes with any of the following conditions: • Children who are less than 6 years of age because their eyes are still in a stage of major growth development that would be disrupted by ocular surgery • Uncontrolled ocular inflammation (e.g. uveitis) • Severe chronic uveitis • Microphthalmos • Untreated retinal detachment • Untreated chronic glaucoma • Rubella cataract • Rubeosis of the iris • Proliferative diabetic retinopathy • Stargardt's retinopathy • Pregnancy • Intraocular infections

3. Device description

Short description of the devices	Foldable iris prosthesis for implantation in the ciliary sulcus or capsular bag
Description of the devices	The Artificial Iris is a foldable iris prosthesis that is used for the treatment of iris defects. It is manufactured from a commercially available, biocompatible ophthalmic silicone, with a 3D-structured front surface and a smooth back surface. In the case of the “with Fiber” variant, a fiber meshwork layer is embedded into the silicone. The Artificial Iris is manufactured as a full 360° iris prosthesis with an overall diameter of 12.8 mm, which can be trephined as needed to custom-fit the device for placement in the posterior chamber (ciliary sulcus or capsular bag). The device can also be cut into a suitable size and shape and sutured to an iris remnant for treatment of partial aniridia or traumatic iris injury. The device is sterile and designed for single-use and can be implanted with a forceps or injector (the latter is only possible for the “Fiber Free” model). During the manufacturing, colored silicone paste is applied by hand in a pattern to match the appearance of a natural iris. Three different color categories are offered, “CustomFlex®”, “Select” and “Standard”. The “CustomFlex®” Artificial Iris is matched for each individual patient (patient-matched). During its manufacturing a photograph of either the existing iris or, in the case of aniridia, a photograph selected by the patient serve as template for painting. In contrast, the “Standard” and “Select” versions of the Artificial Iris are offered ready-made with fixed colors but are produced with an identical technical design and following the same processes. While the “Standard” versions have a unicolor front surface (brown or black), the “Select” versions are painted with predefined multicolor variations (different shades of brown, blue and green available). The posterior surface of each device variant is black to absorb light completely.  Schematic frontal and lateral view of a customized Artificial Iris. The dimensions are given in mm.
Mode of action / operating principle	The device has a fixed aperture of 3.35 mm, with an opaque perimeter and a black posterior surface to absorb light completely, thus reducing photic phenomena and other symptoms associated with aniridia. The device closely mimics the appearance of the natural iris in order to achieve a good cosmetic effect. The small central aperture might increase visual acuity, increase depth of field, and increase contrast sensitivity (pinhole effect).

Materials or substances in contact with the patient's tissues	The Artificial Iris is made of a medical-grade hydrophobic silicone elastomer consisting of diphenyl co dimethyl polysiloxane and silicone reinforcing resin. For colorization inorganic pigments are used. The model With Fiber is additionally reinforced by an embedded polymer fiber meshwork consisting of polyethylene terephthalate (PET).
Method of sterilisation and allowable number of uses	The devices are sterile (steam sterilised) and designed for single use. Reprocessing of an implant is strictly prohibited.
Previous generations of the device(s)	None.
Variants and/or configurations covered by this document	Variants with regards to embedded polyester meshwork layer: • "Fiber Free" (FF) variant: developed for sutureless implantation in the ciliary sulcus or capsular bag • "with Fiber" (wF) variant: developed for implantation in the ciliary sulcus, where suturing is indicated Variants with regards to the colorization of the implant: • CUSTOMFLEX® version: a photo of the existing iris (remnants) or of the iris of the contralateral eye will serve as a template for the manufacturing of the Artificial Iris. Aniridia patients who have no native iris available will use the photo of another person's eye as the template photo. • "Standard" and "Select" versions: they are offered ready-made with fixed colors (technical design and production process are identical to the CUSTOMFLEX® version)
Accessories which are intended to be used in combination with the device	None
Other devices and products which are intended to be used in combination with the device	The Artificial Iris can be folded and implanted via forceps. Only the Fiber Free model can also be implanted using commercially available injector systems. Injectors tested in this context are the Mediceal ViscoJect™ 2.2-1P BIO (Mediceal AG, Switzerland) and the Lenstec Metal screw-type Injector I-9012FS with Cartridge CART 45S (Lenstec Inc. USA). Note: Implantation of a with Fiber model through an injection system can lead to permanent deformation and defects of the implant.

4. Risks and warnings

Residual risks and undesirable effects	Below is a list of the potential complications and undesirable side effects associated with the use of the device, surgical procedure, or IOL. Device-related complications associated with the Artificial Iris may include but may not be limited to: <u>Device-related complications:</u> • Elevated intraocular pressure • Decrease in uncorrected distance visual acuity • Decrease in best-corrected distance visual acuity • Worsening of photosensitivity

	• Eye inflammation • Incorrect device positioning, dislocation, and decentration • Secondary (additional) surgical intervention Surgical repositioning, replacement or removal of the device may be necessary to correct device dislocations. Device defects can occur if the device is not handled properly. <u>Surgery-related adverse events:</u> • Cystoid macular edema • Hypopyon • Endophthalmitis • Device migration • Pupillary block • Retinal detachment • Secondary surgical intervention (unplanned) • Corneal edema, persistent at 3 months or later • Chronic iritis/anterior segment inflammation persistent at 3 months or later
Caution	Implantation of the Artificial Iris is not recommended in patients with the following conditions and situations: • Preoperative intraocular pressure (IOP) above 21 mmHg that does not respond to pressure-lowering medication, unless the IOP above 21 mmHg is due to a known underlying condition that is well controlled with glaucoma treatment, such as ocular hypertension or open-angle glaucoma • Patients with severe endothelial corneal dystrophy, because the surgical procedure to implant the Artificial Iris may damage the cornea sufficiently so that the potential benefits of implantation do not outweigh the risks • The Artificial Iris is not designed for placement in the anterior chamber • No useful vision or visual potential in the fellow eye, unless the patient has debilitating visual symptoms, so that the potential benefits of Artificial Iris implantation clearly outweigh the risks • Presence of a condition or finding in the fellow eye that would make it unsafe to implant an Artificial Iris prosthesis in the eye to be treated • Allergy to any of the planned postoperative antibiotic or anti-inflammatory medications, unless a suitable alternative medication can be prescribed • Implantation for cosmetic color changes of the iris • Post-partum women who are nursing or lactating and for whom postoperative medications are contraindicated • Patients with gastric ulcers or diabetes mellitus in whom high doses of orally administered systemic steroids are required postoperatively Any other condition that would interfere with the planned surgical procedure to implant the iris device
Warnings	The Artificial Iris should be used with caution in the following situations:

	• A clear natural crystalline lens – Do not implant the Artificial Iris in phakic eyes • The visual potential of the fellow eye cannot be evaluated preoperatively (e.g. poor visual acuity due to cataract) • Preoperative IOP > 21 mmHg that is known to be stable and well-controlled with glaucoma treatment (e.g. medication, tubes or shunts) • Presence of any other medical condition that might be expected to make the patient an unsuitable candidate for an Artificial Iris implantation • Anticipated complexity of the planned surgical procedure that might increase the potential for complications • Implantation in the fellow eye before stabilization of the first implanted eye (typically 1 month or more) The aperture of the Artificial Iris pupil is fixed at 3.35 mm. In the event, that a larger pupillary opening is required for posterior segment surgery, the Artificial Iris can be explanted and a new Artificial Iris may be implanted after completion of the posterior segment surgery. The technique for a secondary implantation would be determined in the same manner as for a primary implantation.
Other relevant aspects of safety	No Field Safety Corrective Actions (FSCA) have been taken since the products have been launched.

5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)

5.1. Summary of clinical data related to equivalent device

Not applicable. The conformity of the device was not assessed on the basis of equivalence.

5.2. Summary of clinical data from conducted investigations of the device before the CE-marking

Not applicable. All clinical investigations were performed after the first CE marking under the European Medical Device Directive 93/42/EEC (see section 5.3).

5.3. Summary of clinical data from other sources

Clinical data confirming the safety and performance of the Artificial Iris have been derived from the following sources:

- U.S. FDA pre-market approval (PMA) study initiated by the manufacturer
- A comprehensive literature review using PubMed and other databases

Literature search strategies were designed to identify articles relevant to the devices. Both favorable and unfavorable references were identified and summarized.

5.3.1. Clinical studies

The Artificial Iris was evaluated in an U.S. FDA pre-market approval (PMA) study, which started in October 2013. After completion of the study, the FDA approved the Artificial Iris as first prosthetic iris in the United States on the 30th of May 2018. This study represents the first clinical trial of the Artificial Iris and is the largest known case series world-wide of eyes treated with an iris prosthesis. It covers a large number of adult (n=403) and a reasonable number of pediatric (n=44, ≥ 6 years of age) eyes. No methodological flaws or doubts regarding the scientific validity of the contents could be identified by the experts. Overall, the two main common causes of aniridia, traumatic iris defects and congenital aniridia, were represented by more than 70% of the treated eyes, with both variants (“with Fiber” and “Fiber Free”) and implantation site being used approximately at the same frequency. Therefore, this dataset provides a solid basis to demonstrate safety and effectiveness of the Artificial Iris.

Details of the study results are provided in the table below.

U.S. FDA study (AI-001) (Multicentric pre-market approval/post-approval study with the Artificial Iris – United States) [1]

Based on this summary evaluation, the Artificial Iris, when implanted for the treatment of full or partial aniridia, performs as intended as an iris prosthesis. The results exceed all of the AI-001 protocol target safety and efficacy criteria. The incidence of adverse events, complications, and other ocular or vision-related observations was small. Visual symptoms were reduced, and visual function was increased compared to baseline after the Artificial Iris was implanted. There are no emerging safety or efficacy concerns from this cohort of subjects treated with the Artificial Iris. Specifically, at 12 months after Artificial Iris implantation:

- There was an overall improvement in photosensitivity symptoms, quality of life and vision in pediatric and adult subjects.

- The Artificial Iris implant was correctly positioned in at least 95% of the pediatric and adult eyes throughout the 12-month study.
- Self-administered symptom questionnaire: overall, there was a decrease in the proportion of eyes with marked to severe ratings for each of the evaluated visual symptoms: day-time light sensitivity, night-time light sensitivity, difficulty driving at night, reading difficulty, double vision, fluctuation in vision, glare during the day, glare at night, halos during the day, halos at night, starbursts, eye dryness, pain, and foreign body sensation after Artificial Iris implantation compared to before surgery. Statistical significance was achieved for all symptoms except pain at 12 months. For pediatric eyes, there was a decrease in the proportion of eyes with marked to severe ratings for each visual symptom, except eye dryness increased. Statistical significance was achieved in the pediatric eyes for all symptoms except reading difficulty and eye dryness. There were no pediatric eyes with marked to severe pain or foreign body sensation preoperatively or at 12 months.
- There was a statistically significant improvement in health related quality of life, as measured by the NEI VFQ-25, in the pediatric and adult eyes.
- Postoperative cosmesis was aesthetically pleasing, with 94% of the overall cohort and 100% of the pediatric subjects rating the appearance of their eyes as improved, much improved, or very much improved.
- Postoperative inflammation resolved by 1 month in the majority of all eyes treated, as well as the pediatric eyes; 92% or more of the pediatric and adult eyes had no or trace cell or flare, 90% of pediatric and adult eyes had no corneal stromal edema; and, at least 96% of pediatric and adult eyes had no corneal wound edema at 1 month.
- Mean endothelial cell density remained stable with no clinically significant loss over the 12 months of the study in the pediatric and adult eyes.
- Intraocular pressure was well controlled in the majority of eyes; mean postoperative IOP was within 1.5 mmHg of mean preoperative IOP at all visits for the overall cohort and within 2.6 mmHg for the pediatric eyes.
- Spikes of IOP >30 mmHg were the most commonly reported postoperative adverse events. Surgery related IOP increases >30 mmHg occurred in less than 8% of all the treated eyes and in none of the pediatric eyes.
- Although the Artificial Iris is not a refractive device, 65% of all eyes and 55% of pediatric eyes gained 1 or more lines of CDVA; 67% of all eyes and 64% of pediatric eyes gained 1 or more lines of UDVA; and, mean manifest refraction spherical equivalent was -0.36 D for all eyes and -0.83 D for pediatric eyes at 12 months postoperatively.
- Few eyes lost CDVA; 6.5% of pediatric and 7.9% of the overall study eyes lost 2 or more lines of CDVA postoperatively compared to their preoperative vision at 12 months; and, none of the losses were related to the Artificial Iris.
- Visual symptoms generally improved at 1 month postoperatively. Furthermore, these outcome measures were maintained or improved with further follow-up through 12 months postoperatively.
- Pediatric subjects had postoperative safety and efficacy outcomes that were equivalent or better than adults.

All key safety and efficacy endpoints were met and achieved statistical significance for all eyes treated including the pediatric eyes. Results for the key safety and efficacy endpoints and other study outcomes were similar in each of the study cohorts, the pediatric eyes, and the overall study population.

5.3.2. Literature review

A considerable number of articles about the Artificial Iris has been published. Sixteen scientific publications were selected from the literature search and were critically reviewed with a special focus on safety and efficacy. In summary, if the recommendations provided by the manufacturer are followed correctly, the Artificial Iris when implanted for the treatment of full or partial aniridia, achieves excellent visual and cosmetic outcomes and performs as intended as an iris prosthesis. With attentive postoperative management, the Artificial Iris can be implanted into aniridic aphakic or pseudophakic eyes without concerns in terms of risks of severe adverse events. However, postoperative complications are not uncommon due to preoperative comorbidities and/or surgical procedures. In conclusion, the Artificial Iris can be considered as a safe and useful tool to reduce photic symptoms and to improve quality-of-life of aniridic patients.

5.4. Overall summary of the clinical performance and safety

As prosthetic iris devices (PIDs) are intended to replace or complement the defective natural iris (defects may be partial or complete absence of iris tissue), the main criterion for the success of surgery is the restoration of the iris. The restored pseudo-pupil aperture leads to reduced visual symptoms of photosensitivity, in particular glare and photophobia. Together with the improved cosmetic appearance this is expected to improve the quality of life.

PERFORMANCE ANALYSIS

The U.S. FDA PMA study (chapters 5.3.1) demonstrated a significant improvement in photosensitivity symptoms and satisfaction quality of life and vision. Pediatric subjects had postoperative performance (efficacy) outcomes that were equivalent or better than adults. Satisfaction with cosmesis was high with 94% of the overall cohort and 100% of the pediatric subjects rating the appearance of their eyes as improved, much improved, or very much improved. Statistical significance was not tested since the GAIS scale is not designed to test the cosmesis comparison as a statistical hypothesis, and significant improvement does not need to be demonstrated in this CER since the Artificial Iris is not a cosmetic device. Key performance outcomes are presented in Table 1.

.

Table 1: Key performance outcomes

Based on the 12-month postoperative results from the U.S. FDA PMA study

Endpoint	Threshold rate (acceptance criterion)	Clinical outcomes % for n = 339 eyes ("All eyes combined")
Glare and photophobia		p-value
Severity of day-time light sensitivity symptoms	Statistically significant decrease (p<0.05)	<0.0001
Severity of night-time light sensitivity symptoms	Statistically significant decrease (p<0.05)	<0.0001
Severity of glare during the day	Statistically significant decrease (p<0.05)	<0.0001
Severity of glare at night	Statistically significant decrease (p<0.05)	<0.0001
Patient's self-reported satisfaction with quality of life		p-value
NEI VFQ-25 total score	Statistically significant improvement (p<0.05)	<0.0001
Patient's self-reported satisfaction with cosmetic appearance		%
Satisfaction with cosmesis (GAIS rated as improved, much improved, or very much improved)	None since the Artificial Iris is not a cosmetic device	94.0%

The Artificial Iris completely blocks electromagnetic irradiation in the UV, visible and NIR range [2] and its light-blocking function is the same for any kind of aniridia. Therefore, one can assume that it can reduce photosensitivity in any aniridic patient, including those with transillumination defects like bilateral acute iris transillumination (BAIT), ocular albinism or oculocutaneous albinism.

Taken together, these data show that the Artificial Iris fulfils the pre-defined acceptance criteria and achieves the performance intended by the manufacturer as required by the MDR Annex I [3], when used according to the intended use.

SAFETY ANALYSIS

Since the Artificial Iris is a permanent intraocular prosthesis that can be implanted only in eyes which have been or will be concomitantly implanted with an IOL (pseudophakic, aphakic or requiring cataract removal), the risk assessment should not only consider the risks associated with the Artificial Iris but also with the IOL and associated surgeries (IOL-related risks are appropriately reflected in the Instructions for Use [4]).

Patients who are eligible for the implantation of the Artificial Iris often have concomitant ocular comorbidities. There is often a progressive decline in vision, progressive elevation in IOP requiring glaucoma treatments, surgical interventions such as corneal transplants, and/or corneal-limbal stem cell transplants, especially in cases of congenital aniridia. The comorbidities associated with aniridia, and especially congenital aniridia, need to be considered in the safety analysis and profile of the Artificial Iris. The rates of complications and adverse events (AEs) need to be interpreted with respect to the presenting pathologies leading to the surgical

implantation of the Artificial Iris. Nonetheless, the incidence of AEs (including IOL-related AEs), complications, and other ocular or vision-related observations found in the U.S FDA study was small. (**Key safety outcomes are presented in Table 2.** All pre-defined acceptance criteria were fulfilled.) The most relevant among these were increased intraocular pressure (IOP) and device malpositioning. Overall, the experience and surgical skills of the surgeon are key factors for the success of the Artificial Iris implantation.

Table 2: Key safety outcomes for endpoints “stable visual acuity” and “acceptable rate of complications”

Based on the 12-month postoperative results from the U.S. FDA PMA study

Events at 12 months	Threshold rate (acceptance criterion)	Clinical outcomes % for n = 447 eyes (“All eyes combined”)
Device-related adverse events	-	
>2 Line (>10 Letters) Loss of CDVA that is device related	<5%	0.0%
CDVA of worse than 20/40 if the pre-op CDVA was 20/40 or better (device related)	<5%	0.0%
Cumulative Lens Related Adverse Events	<5%	0.0%
Anisometropia	<1%	0.0%
Glare/halos	<1%	0.0%
Diplopia	<1%	0.0%
IOL removal or replacement due to lens power calculation error	<1%	0.0%
Cumulative Surgery Related Adverse Events	-	
Cystoid Macular Edema	18.8%	2.9%
Hypopyon	3.0%	0.0%
Endophthalmitis	3.0%	0.0%
Device migration	5.4%	0.0%
Pupillary block	7.8%	0.0%
Retinal detachment	5.4%	0.7%
Secondary surgical intervention (unplanned)	8.5%	1.3%
Corneal edema, persistent at 3 months or later	4.2%	1.6%
Chronic anterior segment inflammation, persistent at 3 months or later (chronic iritis)	5.0%	1.8%

The clinical data evaluated provide strong evidence that the Artificial Iris is a safe treatment for congenital or acquired aniridia. Although postoperative complications are not uncommon after device implantation, their rates are acceptable since aniridia from either a traumatic or genetic origin is usually accompanied by other ocular defects and malfunctions. The most common complications encountered were prolonged postoperative inflammation, raised IOP/secondary glaucoma, and malpositioning of the implant.

As already mentioned earlier, it should be noted that the rate of complications clearly depends on the experience and surgical skills of the surgeon. Therefore, surgeons have to be certified according to certification process from HumanOptics and must have experience in the surgical management and treatment of aniridia and iris defects. The publication of Mayer et al. from 2018 [5] reported a significant reduction of complications after twelve implantations, this implicates that the procedure is not to be recommended in low-volume settings or for surgeons with little experience.

The endothelial cell count decreased after Artificial Iris implantation but were comparable to the loss in normal cataract surgery.

If an Artificial Iris With Fiber is cut to size, it appears to be essential that a sufficiently sharp trephine is used (e.g. a single-use device). Scanning electron microscopy revealed a variety regarding the smoothness and roughness of the edges, if less sharp trephines were used [5]. Chronic movement of the Artificial Iris on the ciliary body may in such cases lead to chronic irritations of the ciliary body and bleedings [5]. It can be speculated that this was also the root cause in another early publication [6] where an increased rate of certain complications (glaucoma; over-pigmentation of the remaining iris tissue; need for secondary surgery) with the variant with mesh compared to the Fiber Free variant was found.

In one study, a dilation of the residual iris aperture was described [7] in a significant proportion (16.7%) of eyes which had been implanted with the Artificial Iris. Such eyes exhibited a relatively high rate of other complications. However, this “residual iris retraction syndrome” has not been described in any other publications.

Further evidence for the safety of the device under consideration comes from the low rate of complaints and the fact that FSCAs have not been necessary since the start of marketing.

Since 2002 HumanOptics has placed more than 5,000 devices on the market. Over this period of more than 20 years no feedback has been received relating to material degradation. The manufacturer is aware of rare cases of device explantation unrelated to device failure or malfunction. The available information confirms that HumanOptics Artificial Iris is stable over a period of at least 20 years (based on simulated aging tests according to ISO 11979-5). It can be expected that the device is stable indefinitely from the implantation date over the lifetime of the patient.

In summary, one can conclude that the device under consideration Artificial Iris represents the current state of the art and fulfils the pre-defined acceptance criteria and the requirements regarding safety as set out by the MDR Annex I [3]. The Artificial Iris is safe when used according to the intended use.

BENEFIT-RISK ANALYSIS

The **benefits** of a PID implantation are significant. The main reported beneficial effects of prosthetic iris implants include enhanced visual function and improved cosmetic appearance both of which lead to an improvement of the quality of life of aniridic patients [8]. The improvement of vision is attributable to various factors, including reduced spherical and chromatic aberration, increased depth of focus, and reduced glare and photophobia. It is understood that the extent of benefits depends on the potential for improvement, i.e., on the degree of patients’ preoperative discomfort and on ocular comorbidities which may limit the maximum achievable visual results.

With the HumanOptics Artificial Iris the following benefits have been documented:

IMPROVED VISUAL FUNCTION:

- Photophobia, glare and other visual symptoms: The U.S. FDA PMA study demonstrated reduced photosensitivity to daytime and night-time lighting and associated glare. There was also a decrease in the proportion of eyes with marked to severe ratings for each of the other evaluated visual symptoms after Artificial Iris implantation compared to before surgery (difficulty driving at night, reading difficulty, double vision, fluctuation in vision, glare during the day, glare at night, halos during the day, halos at night, starbursts, eye dryness, pain, and foreign body sensation). Statistical significance was achieved for all symptoms except pain at 12 months. Visual symptoms generally improved at 1 month postoperatively and these outcome measures were maintained or improved with further follow-up.

A reduction of subjective symptoms of glare and/or photophobia has been confirmed in large number publications [9-17]. Also a significant improvement of CDVA with glare has been reported [9]. Although that group had reported a significant reduction in subjective glare at daytime and nighttime in 2020 [9], scores did however not significantly improve in a later study from the same group [18] (yet, individual glare scores for daytime and nighttime were improved or stable in the majority of patients). As the latter study was performed with patients who underwent a triple procedure (Artificial Iris + IOL+ penetrating keratoplasty), the authors of the publication suggested that corneal irregularity from the transplant or early corneal failure might have replaced the light and glare sensitivity caused by the iris defect.

- Visual acuity: Although the Artificial Iris is not a refractive device, the procedure improves visual acuity due to presence of an artificial pupil combined with IOL implantation for correction of refractive errors. In the U.S. FDA PMA study, 65% of all eyes and 55% of pediatric eyes gained 1 or more lines of CDVA; 67% of all eyes and 64% of pediatric eyes gained 1 or more lines of UDVA at 12 months postoperatively. Only 7.9% of the overall study eyes and 6.5% of pediatric eyes lost 2 or more lines of CDVA postoperatively and none of the losses were related to the Artificial Iris.

Also other publications (see [6, 9, 11, 13, 16, 19]) did show an improvement of the visual acuity. Of note, one primary investigator reported no improvement in visual acuity in an early study from 2016 [14] but observed a significant positive effect in a later study [13]. This might be related to his increasing surgical experience over the years. In the case series by Crawford et al., 2022 [10], CDVA was improved or stable in all but one eyes (18/19) and the levels of satisfaction with postoperative visual acuity outcomes were high (mean 4.3/5.0) but the improvement in median CDVA value was not significant.

- Contrast sensitivity: Contrast sensitivity was not measured in the U.S. FDA study. But it has been assessed in two publications by Mayer et al. from 2016 and 2020 [13, 14] and was found to be significantly improved under photopic conditions with glare.

IMPROVED COSMETIC APPEARANCE:

In the U.S. FDA PMA study, patients' satisfaction with the postoperative cosmesis was high. 94% of the overall cohort and 100% of the pediatric subjects rating the appearance of their eyes as improved, much improved, or very much improved.

Many publications have confirmed that the patients' cosmesis satisfaction is significantly improved [9, 10, 14-16, 18] and that cosmetic results are very good (e.g. refs. [12, 15, 20]). For example, in the study by Roman et al., 2021 [15], the mean cosmesis score improved by 4.7 points ($P = 0.001$) on a 0 to 10 scale.

IMPROVED QUALITY OF LIFE:

As discussed at the beginning of this chapter, PID implantation can improve the quality of life for aniridic patients through improved vision and cosmetic appearance. Visual health-related quality of life, as measured by the NEI VFQ-25 total score, was found to be significantly improved after implantation of the Artificial Iris in the U.S. FDA PMA study and in the article by Bonnet et al. 2021 [18]. Similarly, Yildirim et al., 2020 [20] reported that the majority of patients (92.4%) would repeat the treatment with the Artificial Iris if they could decide again and that their quality of life did improve by 7.99 ± 1.85 points ($P = 0.001$) on a 1 to 10 scale. A high patient satisfaction with the overall result (8.97 ± 1.42 out of 10 points) was also reported by Mayer et al., 2016 [14].

Regarding **risks**, the observed complications are in line those reported for PID surgery in general. The assessment of risks included the risks associated with the Artificial Iris and surgery but also IOL-related risks (since an IOL has to be implanted before or concurrently with the Artificial Iris surgery).

The rate of postoperative complications is rarely device-related but rather depends on the preoperative ocular comorbidity, which means that PIDs should be implanted into risky patients only after careful assessment of the benefit/risk ratio. Furthermore, the success of Artificial Iris implantation is influenced by the surgical technique and equipment used and involves a learning curve for this surgical intervention.

In summary, the available clinical data do show that the Artificial Iris is a reliable product and that the remaining risks are well controlled.

In summary, a positive benefits-versus-risks profile can be assigned to the medical device Artificial Iris under consideration.

Intraoperative and postoperative complications do exist but occur with a low frequency and are acceptable, when weighed against the benefits. Compared to other PIDs, the HumanOptics Artificial Iris represents a significant advancement in the treatment options for patients who suffer from glare and other visual disturbances secondary to aniridia.

Both the type and frequency of risks and benefits were comparable for adults and children. However, a limitation of this dataset is the fact that the youngest child in the pediatric cohort was 6 years old and the number of pediatric eyes ($n=44$ in the U.S. FDA study) was relatively small. Therefore, it is recommended to collect more data for the use in children.

The manufacturer's own clinical investigation (U.S. FDA study) as well as other published studies with the Artificial Iris provide strong evidence that the device can be successfully used for the treatment of full or partial aniridia. There is sufficient evidence that the intended clinical benefits are achieved. Based upon the fact that this treatment has a low risk profile, it is concluded that a positive benefits versus risks profile can be stated for the device.

5.5. Ongoing or planned post-market clinical follow-up

In order to confirm the long-term safety of the Artificial Iris, the follow up of the U.S. FDA study (AI-001) has been extended into a post-approval study, and long-term real-world data on complications are being collected in a European country. In addition, a clinical investigation with a limited number of eyes is currently running which focusses on the outcomes with the device in an Asian population.

6. Possible diagnostic or therapeutic alternatives

There are different options to treat iris defects, which are iris suturing, the use of colored "iris" contact lenses, and corneal tattooing and implantable prosthetic iris devices.

There are several PID models that can be categorized into 2 groups: iris diaphragms (with or without lens) and endocapsular capsular tension ring (CTR)-based PID.

7. Suggested profile and training for users

The Artificial Iris must be handled by health professionals and implanted by ophthalmic surgeons. A high level of surgical skill and experience in the field of anterior segment surgery is required for Artificial Iris implantation and experience in the surgical management and treatment of aniridia and iris defects is a prerequisite. Before placing the first order, the surgeon must have successfully completed a Certification Course from HumanOptics, which provides current basic knowledge of device applications and surgical techniques. In addition, it is advisable to take part in wet labs with the Artificial Iris and to make use of the service for non-sterile training samples. Application guides may be available on request. Product application videos can be found on www.humanoptics.com in the “application” section for the Artificial Iris.

If you have further questions on the usage of the device, please contact your local sales representative or the manufacturer (customerservice@humanoptics.com).

8. Reference to any harmonised standards and CS applied

Common specification(s) to comply with, if applicable:
N/A
Harmonised standards to apply, if applicable:
DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
DIN EN ISO 14971:2020-07 Medical devices – Application of risk management to medical devices

9. References

- [1] Ayres BD, Fant BS, Landis ZC, Miller KM, Stulting RD, Cionni RJ, et al. Results of the United States Food and Drug Administration Clinical Trial of the CustomFlex Artificial Iris. *Ophthalmology*. 2022;129(6):614-625.
- [2] Yusuf ICH, Fung THM, Perison SN, Patel CK. Artificial Iris Implants for Aniridia: Does Near Infra-Red Light Transmission Risk Retinal Phototoxicity? - Presentation at the Congress of the Association for Research in Vision and Ophthalmology (ARVO). 2014.
- [3] Regulation EU 2017/745 on medical devices (MDR)
- [4] Instructions for Use for the Artificial Iris (V07/2023-10)
- [5] Mayer CS, Laubichler AE, Khoramnia R, Tandogan T, Prahs P, Zapp D, et al. Challenges and complication management in novel artificial iris implantation. *Journal of ophthalmology*. 2018;2018.
- [6] Rickmann A, Szurman P, Januschowski K, Waizel M, Spitzer MS, Boden KT, et al. Long-term results after artificial iris implantation in patients with aniridia. *Graefes Archive for Clinical and Experimental Ophthalmology*. 2016;254(7):1419-1424.
- [7] Mayer CS, Laubichler AE, Masyk M, Prahs P, Zapp D, Khoramnia R. Residual iris retraction syndrome after artificial iris implantation. *American Journal of Ophthalmology*. 2019;199:159-166.
- [8] Srinivasan S, Ting DS, Snyder ME, Prasad S, Koch HR. Prosthetic iris devices. *Can J Ophthalmol*. 2014;49(1):6-17.
- [9] Bonnet C, Miller KM. Safety and efficacy of custom foldable silicone artificial iris implantation: prospective compassionate-use case series. *Journal of Cataract & Refractive Surgery*. 2020;46(6):893-901.
- [10] Crawford AZ, Freundlich SEN, Lim J, McGhee CNJ. Endocapsular artificial iris implantation for iris defects: Reducing symptoms, restoring visual function and improving cosmesis. *Clin Exp Ophthalmol*. 2022;50(5):490-499.
- [11] Figueiredo GB, Snyder ME. Long-term follow-up of a custom-made prosthetic iris device in patients with congenital aniridia. *Journal of Cataract & Refractive Surgery*. 2020;46(6):879-887.
- [12] Khan HM, Eslami M, Khan HN. Long-term outcomes of patients with artificial iris implants. *Can J Ophthalmol*. 2023.
- [13] Mayer CS, Hoffmann AM, Prahs P, Reznicek L, Khoramnia R. Functional outcomes after combined iris and intraocular lens implantation in various iris and lens defects. *BMC Ophthalmol*. 2020;20(1):370.
- [14] Mayer CS, Reznicek L, Hoffmann AE. Pupillary reconstruction and outcome after artificial iris implantation. *Ophthalmology*. 2016;123(5):1011-1018.
- [15] Roman S, Baudouin C. Flexible silicone artificial iris in cases of aniridia and iris deficiencies. *Journal Français d'Ophtalmologie*. 2021;44(9):1387-1395.
- [16] Spitzer MS, Nessmann A, Wagner J, Yoeruek E, Bartz - Schmidt KU, Szurman P, et al. Customized humanoptics silicone iris prosthesis in eyes with posttraumatic iris loss: outcomes and complications. *Acta Ophthalmologica*. 2016;94(3):301-306.
- [17] Wang S, Liu M, Li X, Ge C, Chen W, Zhao Z. Preliminary assessment of the safety and effectiveness of Artificial Iris implantation in China. 19 August 2022, PREPRINT (Version 1) available at Research Square [<https://doi.org/10.21203/rs.3.rs-1894256/v1>]. 2022.
- [18] Bonnet C, Vazirnia P, Deng SX, Aldave AJ, Miller KM. Sutured custom foldable silicone artificial Iris implantation combined with intraocular lens implantation and penetrating keratoplasty: safety and efficacy outcomes. *Cornea*. 2021;40(10):1236-1247.

- [19] Mayer C, Khoramnia R. "Double prosthesis implantation": biometry and refractive outcomes in combined intraocular lens and artificial iris surgery. *Clinical Ophthalmology (Auckland, NZ)*. 2021;15:799-805.
- [20] Yildirim TM, Khoramnia R, Masyk M, Son H-S, Auffarth GU, Mayer CS. Aesthetics of iris reconstruction with a custom-made artificial iris prosthesis. *PLoS One*. 2020;15(8):1-12.

Part 2 – SSCP intended for patients

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your ophthalmologist in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions For Use to provide information on the safe use of the device.

This SSCP is reviewed at least annually and will be updated if needed.

1. Device identification and general information

Device trade name(s)	Artificial Iris • Variant Artificial Iris with Fiber (wF) • Variant Artificial Iris Fiber Free (FF) Different color categories ("CustomFlex [®] ", "Select" and "Standard") are available for both the wF and FF variants.
Basic UDI-DI(s)	04049154_AI_M1_H1_O1_X7
Year when the first certificate (CE) was issued covering the device	2011

2. Intended use of the device

Intended purpose	The Artificial Iris is a small piece of biocompatible silicone intended for use as an iris prosthesis to replace the natural iris which is partly or fully gone, e.g. due to congenital aniridia or due to a trauma on the eye after an accident. It can be implanted either in the capsular bag² or in the ciliary sulcus³ of the eye. The iris prosthesis cannot be implanted if the natural crystalline lens is still present in the eye. See section 3 "mode of action/operation principle".
Indication(s) and intended patient groups	Adults or children (≥6 years) for whom the natural iris is partly or fully gone might be good candidates for the Artificial Iris. In medicine, the lack of natural iris is called "aniridia". The graph on the right below represents a partial aniridia. normal  aniridia 


² **Capsular bag:** The capsule which the natural lens sits in.

³ **Ciliary sulcus:** The space between the iris and the ciliary body (a muscular tissue right behind the iris)


	Note: there are no clinical data available for breastfeeding women and immunocompromised persons in connection with the implantation of the device.
Contraindications	The Artificial Iris is not to be used for purely cosmetic reasons. The device may not be implanted in eyes with any of the following conditions: • Children who are less than 6 years of age because their eyes are still in a stage of major growth development that would be disrupted by ocular surgery • Uncontrolled ocular inflammation or infections • Extremely small eyes • Certain retinal diseases • Untreated chronic glaucoma • Rubeosis of the iris • Problems with your eyes due to diabetic disease • Pregnancy The different contraindications should be discussed with your ophthalmologist.

3. Device description

Description of the device	The Artificial Iris is a small piece of biocompatible silicone with a hole in the middle that represents the pupil (3.35 mm). On the front side the Artificial Iris has a structured and colored surface imitating the structure and appearance of a natural iris. Below, a picture showing a natural iris and an Artificial Iris which has been painted to imitate the iris on the photo. 
Mode of action / operating principle	In everyday life and work, there are many ways to injure your eyes, such as sharp objects and corrosive chemicals. If the natural iris is damaged partially or completely, it is called aniridia. Apart from that, aniridia can also be congenital or have other reasons. Such iris defects severely affect your vision and can have a negative impact on your appearance, reducing your quality of life. For aniridic patients, the treatment of the eye with a medical implant is often the only adequate option. There are various artificial iris implants available on the market for treating aniridia, but traditional implants are primarily functional and often neglect aesthetics. Often, they are only available in a limited selection of standard colors and the tissue structure is

	poorly imitated. Some of these implants are made of rigid materials and require relatively large incisions during implantation, which increases the risk of infection and astigmatism. The Artificial Iris devices from HumanOptics largely overcome these issues. They are made of high-quality medical silicone, foldable, and can be implanted through minimally invasive surgery. Small incision techniques significantly reduce the risk of infection, minimize astigmatism⁴, and speed up the healing process. Additionally, based on the patient's healthy eye's photo template, Artificial Iris experts can precisely transfer the color pattern onto the implant, thus restoring the appearance of the affected eye. Patients can fully engage again in a social life. To reconstruct a defective iris with the Artificial Iris, a surgical procedure is required to restore the eye's aperture function (pupil function), aiming to reduce sensitivity to glare, to improve depth of vision, and to optimize contrast vision. With the defined pupil size (3.35 mm), the iris implant offers a good compromise between day and night vision (the natural pupil is smaller during the day and larger at night). Importantly, the Artificial Iris cannot be implanted if your natural crystalline lens is still present. So, the first step will be to remove the crystalline lens, which is located in the capsular bag of the eye, and to replace it with an intraocular lens (IOL). It is the same procedure as for a cataract surgery. If an IOL is already present in your eye, the Artificial Iris can then be implanted. Below the different steps and surgical options for the implantation of the Artificial Iris. Step 1 – The patient is diagnosed with aniridia, and their natural lens has either been removed or is planned to be removed during the surgery.  The diagram shows two cross-sections of an eye. The left eye, labeled 'normal', has a blue iris and a black pupil. The right eye, labeled 'aniridia', has a blue iris but a large, dark, irregularly shaped area where the pupil would be, indicating the absence of a natural iris.

⁴ **Corneal astigmatism:** This is present if the cornea is not smooth and round in shape but irregularly shaped


Step 2 – The insertion of the Artificial Iris is performed under local or general anesthesia. The flexible Artificial Iris is rolled up and inserted through a few mm incision at the corneal edge of the eye using an injector or forceps.



Step 3 – The Artificial Iris is then unfolded and fixed in the eye. The implant's placement is chosen by the surgeon based on several factors, including the status of your eye. The main types of placement include:

Artificial Iris, IOL and capsular tension ring placed in capsular bag

Artificial Iris placed in sulcus, IOL placed in the capsular bag

	IOL fixed to the Artificial Iris, implantation in the sulcus  As soon as the protective dressing is removed from the eye after the operation, the treatment result is visible and the new iris can assume its light-limiting function.
Materials or substances in contact with the patient's tissues	The Artificial Iris is made of a medical grade hydrophobic silicone elastomer consisting of diphenyl co dimethyl polysiloxane and silicone reinforcing resin. For colorization inorganic pigments are used. The model With Fiber is additionally reinforced by an embedded polymer fiber meshwork consisting of polyethylene terephthalate (PET).
Accessories	None which are relevant to patients.

4. Risks and warnings

Contact your ophthalmologist if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your ophthalmologist if needed.

How potential risks have been controlled or managed	Potential risks have been identified and controlled according to international standards that have been developed by experts from many countries. The corresponding document(s) comprising rules, guidelines and processes are approved and published by a globally recognized organisation.
Remaining risks and undesirable effects	Overall, the Artificial Iris implantation is a safe procedure. However, postoperative complications are not uncommon due to preexisting ocular problems and/or surgical procedures. Please discuss your individual risks with your ophthalmologist. For every ocular surgery, there are some common side effects. The following ones may occur following iris surgery and normally disappear after a short time: • Foreign body sensation caused by the incision through which the Artificial Iris was inserted • Mild eye redness The following symptoms may indicate complications associated with an artificial iris implantation. If these occur you should consult your doctor: • Sudden or gradual deterioration in vision/blurred vision, flashes of light or dark running dots (floaters) • Increased sensitivity to light • Excessive tearing of the eye • Severe redness and/or itching of the eye • Eye pain

	Nausea, vomiting, severe headache
Warnings and precautions	Patients with certain eye diseases and/or malformations have higher-than-normal risks to develop complications during and after Artificial Iris implantation. Before the surgery, your ophthalmologist will perform an assessment of your eye and will discuss with you anything to do with your operation. This includes the risks and benefits of surgery and how long it will take to fully recover after the surgery. Replacement or removal of the iris prosthesis may also be required in case of serious and uncontrolled complications after the surgery.
Summary of any field safety corrective action	No field safety corrective actions ⁵ (FSCA) have been taken since the Artificial Iris has been sold on the market. In other words, there is no safety concern necessitating a FSCA.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1. Clinical background of the device

Aniridia, i.e. the natural iris is damaged partially or completely, can be congenital or acquired. Affected patients suffer from severe glare and their external appearance can be a source of distress. The treatment of the eye with a medical implant is often the only adequate option. To reconstruct a defective iris with Artificial Iris, a short surgical procedure is required to restore the eye's aperture function, aiming to reduce sensitivity to glare, improve depth of vision, and optimize contrast vision.

After you and your ophthalmologist decide to get an Artificial Iris implanted, more examinations will be conducted. The results will help the surgeon pick the best Artificial Iris for you. Additionally, you'll also talk about the different models of Artificial Iris that you can choose from. Further information about the implantation of Artificial Iris can be found under:

<https://www.humanoptics.com/en/patienten/help-with-aniridia/>

The artificial iris will normally remain in your eye for the rest of your life, provided there is no medical contraindication. The product lifetime of the HumanOptics Artificial Iris is proven for twenty years. Regular ophthalmological check-ups are recommended, in consultation with the treating physician.

5.2. The clinical evidence for the CE-marking

The letters 'CE' which appear on products or packages signify that products sold in the European Economic Area (EEA) have been assessed to meet high safety, health, and environmental protection requirements.

Clinical evidence for CE-marking of ocular implants should meet the requirements defined by international standards and is based on clinical data with the device, scientific literature and feedback from the market.

The Artificial Iris has a long history on the market (20 years) and has been well studied in Europe and the U.S.A. The implant has shown a good safety and performance profile as reported in numerous clinical studies including the U.S. study, in which more than 300 patients (adults and children) were evaluated up to 1 year after surgery. Overall, the results indicate that the Artificial

⁵ **Field safety corrective action:** Corrective action taken by a manufacturer for technical or medical reasons to prevent or reduce the risk of a serious incident in relation to a device made available on the market.

Iris performs as intended for the treatment of full or partial aniridia with improvement in photosensitivity symptoms, improvement in quality of life and vision and with excellent cosmetic results.

5.3. Safety

The implantation of the CE-marked Artificial Iris considerably improves photosensitivity symptoms, quality of life and vision in patients with full or partial aniridia. The CE marking for ocular implants is mandatory and means that the manufacturer has verified that the products meet European safety, health or environmental requirements.

No surgery is without risk. Complications, although rare, can occur during the surgery and/or after the surgery. Most of the time, complications are related to the condition of your eye before the surgery. For instance, if your eye is in a very bad condition because of a severe trauma, the surgeon will have more difficulties to perform the Artificial Iris implantation. Patients with congenital aniridia are also more likely to suffer from complications because the eye has several preexisting ocular diseases. Your ophthalmologist can help you understand all the risks associated with the surgical procedure based on your eye condition.

HumanOptics, as a manufacturer, continuously analyses complaints, market feedback and new scientific publications. The results of these activities are currently favorable for the Artificial Iris and support its safety and performance for the treatment of full or partial aniridia. In addition, the U.S. Artificial Iris study has been extended and two other studies/data collections are ongoing to confirm the long-term safety of the device and its usefulness for Asian patients.

6. Possible diagnostic or therapeutic alternatives

When considering alternative treatments, it is recommended to contact your ophthalmologist who can take into account your individual situation. Currently, most aniridia cannot be treated with drugs. Although sunglasses or colored contact lenses can help to shade the light, iris implants are still the most ideal option. So aniridia should be treated by surgically implanting an artificial iris into the eye. The Artificial Iris can be placed in the capsular bag or in the ciliary sulcus depending on the individual patient. If the iris damage is minimal, the pupil can be reconstructed by suturing the iris tissue rather than implanting an iris prosthesis.

There are different types of iris implants on the market, each with their own features. But based on current data, Artificial Iris is quite a valuable option. Of course it depends on several factors including your eye condition and your lifestyle. The different options should be discussed with your ophthalmologist.

7. Suggested training for users

Not applicable to patients since the device is to be used only by healthcare professionals.

SSCP Artificial Iris Rev02 2025-08

Final Audit Report

2025-08-04

Created:	2025-08-01
By:	Dr. Simon Leierseder (leierseder@humanoptics.com)
Status:	Signed
Transaction ID:	CBJCHBCAABAAqT_Mt-9Q6lUDyly8lCYMd7JieKbKSGNG

"SSCP Artificial Iris Rev02 2025-08" History

-  Document created by Dr. Simon Leierseder (leierseder@humanoptics.com)
2025-08-01 - 8:31:07 AM GMT
-  Document e-signed by Dr. Simon Leierseder (leierseder@humanoptics.com)
Signature Date: 2025-08-01 - 8:31:38 AM GMT - Time Source: server
-  Document emailed to Ulrike Bormann (bormann@humanoptics.com) for signature
2025-08-01 - 8:31:40 AM GMT
-  Email viewed by Ulrike Bormann (bormann@humanoptics.com)
2025-08-01 - 8:32:17 AM GMT
-  Document e-signed by Ulrike Bormann (bormann@humanoptics.com)
Signature Date: 2025-08-01 - 8:34:10 AM GMT - Time Source: server
-  Document emailed to Dr. Timo Abbe (abbe@humanoptics.com) for signature
2025-08-01 - 8:34:12 AM GMT
-  Email viewed by Dr. Timo Abbe (abbe@humanoptics.com)
2025-08-04 - 11:06:15 AM GMT
-  Document e-signed by Dr. Timo Abbe (abbe@humanoptics.com)
Signature Date: 2025-08-04 - 11:06:25 AM GMT - Time Source: server
-  Agreement completed.
2025-08-04 - 11:06:25 AM GMT

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9. Referenzen

- [1] Ayres BD, Fant BS, Landis ZC, Miller KM, Stulting RD, Cionni RJ, et al. Results of the United States Food and Drug Administration Clinical Trial of the CustomFlex Artificial Iris. *Ophthalmology*. 2022;129(6):614-625.
- [2] Yusuf ICH, Fung THM, Perison SN, Patel CK. Artificial Iris Implants for Aniridia: Does Near Infra-Red Light Transmission Risk Retinal Phototoxicity? - Presentation at the Congress of the Association for Research in Vision and Ophthalmology (ARVO). 2014.
- [3] Regulation EU 2017/745 on medical devices (MDR)
- [4] Instructions for Use for the Artificial Iris (V07/2023-10)
- [5] Mayer CS, Laubichler AE, Khoramnia R, Tandogan T, Prahs P, Zapp D, et al. Challenges and complication management in novel artificial iris implantation. *Journal of ophthalmology*. 2018;2018.
- [6] Rickmann A, Szurman P, Januschowski K, Waizel M, Spitzer MS, Boden KT, et al. Long-term results after artificial iris implantation in patients with aniridia. *Graefes Archive for Clinical and Experimental Ophthalmology*. 2016;254(7):1419-1424.
- [7] Mayer CS, Laubichler AE, Masyk M, Prahs P, Zapp D, Khoramnia R. Residual iris retraction syndrome after artificial iris implantation. *American Journal of Ophthalmology*. 2019;199:159-166.
- [8] Srinivasan S, Ting DS, Snyder ME, Prasad S, Koch HR. Prosthetic iris devices. *Can J Ophthalmol*. 2014;49(1):6-17.
- [9] Bonnet C, Miller KM. Safety and efficacy of custom foldable silicone artificial iris implantation: prospective compassionate-use case series. *Journal of Cataract & Refractive Surgery*. 2020;46(6):893-901.
- [10] Crawford AZ, Freundlich SEN, Lim J, McGhee CNJ. Endocapsular artificial iris implantation for iris defects: Reducing symptoms, restoring visual function and improving cosmesis. *Clin Exp Ophthalmol*. 2022;50(5):490-499.
- [11] Figueiredo GB, Snyder ME. Long-term follow-up of a custom-made prosthetic iris device in patients with congenital aniridia. *Journal of Cataract & Refractive Surgery*. 2020;46(6):879-887.
- [12] Khan HM, Eslami M, Khan HN. Long-term outcomes of patients with artificial iris implants. *Can J Ophthalmol*. 2023.
- [13] Mayer CS, Hoffmann AM, Prahs P, Reznicek L, Khoramnia R. Functional outcomes after combined iris and intraocular lens implantation in various iris and lens defects. *BMC Ophthalmol*. 2020;20(1):370.
- [14] Mayer CS, Reznicek L, Hoffmann AE. Pupillary reconstruction and outcome after artificial iris implantation. *Ophthalmology*. 2016;123(5):1011-1018.
- [15] Roman S, Baudouin C. Flexible silicone artificial iris in cases of aniridia and iris deficiencies. *Journal Français d'Ophthalmologie*. 2021;44(9):1387-1395.
- [16] Spitzer MS, Nessmann A, Wagner J, Yoeruek E, Bartz - Schmidt KU, Szurman P, et al. Customized humanoptics silicone iris prosthesis in eyes with posttraumatic iris loss: outcomes and complications. *Acta Ophthalmologica*. 2016;94(3):301-306.
- [17] Wang S, Liu M, Li X, Ge C, Chen W, Zhao Z. Preliminary assessment of the safety and effectiveness of Artificial Iris implantation in China. 19 August 2022, PREPRINT (Version 1) available at Research Square [<https://doi.org/10.21203/rs.3.rs-1894256/v1>]. 2022.
- [18] Bonnet C, Vazirnia P, Deng SX, Aldave AJ, Miller KM. Sutured custom foldable silicone artificial Iris implantation combined with intraocular lens implantation and penetrating keratoplasty: safety and efficacy outcomes. *Cornea*. 2021;40(10):1236-1247.

- [19] Mayer C, Khoramnia R. "Double prosthesis implantation": biometry and refractive outcomes in combined intraocular lens and artificial iris surgery. *Clinical Ophthalmology (Auckland, NZ)*. 2021;15:799-805.
- [20] Yildirim TM, Khoramnia R, Masyk M, Son H-S, Auffarth GU, Mayer CS. Aesthetics of iris reconstruction with a custom-made artificial iris prosthesis. *PLoS One*. 2020;15(8):1-12.

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