Summary of safety and clinical performance (SSCP) for the Aspira-aXA/-aXAY

SSCP reference number: P031

Revision history

Rev number	Revision date	Description of change	Revision validated by the Notified Body
Rev.01	2025-03-11	Initial creation	☒ Yes Validation language: English ☐ No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB)

Table of contents

	Page
Part 1 – SSCP intended for healthcare professionals	3
1. Device identification and general information.....	3
2. Intended use of the device.....	4
3. Device description	4
4. Risks and warnings	6
5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)	8
5.1. Summary of clinical data related to equivalent device	8
5.2. Summary of clinical data from conducted investigations of the device before the CE-marking	8
5.3. Summary of clinical data from other sources.....	8
5.4. Overall summary of the clinical performance and safety.....	15
5.5. Ongoing or planned post-market clinical follow-up	19
6. Possible diagnostic or therapeutic alternatives	20
7. Suggested profile and training for users	20
8. Reference to any harmonised standards and CS applied	20
9. References	21
Part 2 – SSCP intended for patients	22
1. Device identification and general information.....	22
2. Intended use of the device.....	22
3. Device description	23
4. Risks and warnings	25
5. Summary of clinical evaluation and post-market clinical follow-up	26
5.1. Clinical background of the device.....	26
5.2. The clinical evidence for the CE-marking	26
5.3. Safety.....	27
6. Possible diagnostic or therapeutic alternatives	27
7. Suggested training for users.....	27

Part 1 – SSCP intended for healthcare professionals

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The SSCP is not intended to replace the Instructions For Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients. This SSCP is reviewed at least annually and will be updated if needed to ensure that any clinical and/or safety information in the SSCP remains correct and complete.

The following information is intended for users/healthcare professionals.

Following this information there is a summary intended for patients.

1. Device identification and general information


Manufacturer's name and address	HumanOptics Holding AG Spardorfer Straße 150 91054 Erlangen Germany
Manufacturer's single registration number (SRN)	DE-MF-000017892
Authorised representative if applicable; name and SRN	N/A
Notified body name and organization number	TÜV NORD CERT GmbH NB 0044

Device trade name(s)	- Aspira-aXA - Aspira-aXAY
Basic UDI-DI(s)	Aspira-aXA/-aXAY: 04049154_PC_M2_H4_O2_CJ
CND/EMDN code(s):	Aspira-aXA/-aXAY: P030102090202 (IOLs, APHAKIC, MONOFOCAL, ASPHERICAL, HYDROPHILIC ACRYLIC)
Class of device	☐ I ☐ IIa ☒ IIb ☐ III
Year when the first certificate (CE) was issued covering the device	- Aspira-aXA: 2016 - Aspira-aXAY: 2020

2. Intended use of the device

Intended purpose	The devices are intended for implantation in the capsular bag to replace the natural lens.
Intended users	The devices must be handled by health professionals and implanted by ophthalmic surgeons.
Use environment	Ophthalmic operating room
Indications	The devices are indicated for the correction of aphakia after surgical removal of the natural lens.
Medical condition(s) to be treated	Aphakia after surgical removal of the natural lens
Intended patient population	Aphakic adult patients Note: children are excluded due to the absence of clinical data and additional risks associated with paediatric cataract surgery. The manufacturer does not hold clinical data relating to pregnant/breastfeeding women or immunocompromised patient populations.
Intended anatomical location	Posterior chamber of the eye
Contraindications	Outside of general contraindications of ocular surgery, there are no specific contraindications for this device.

3. Device description

Short description of the device	One-piece, aspheric aberration-free, monofocal posterior-chamber intraocular lens (IOL) with cut-out haptics and a 360° lens epithelial cell barrier on the posterior optic surface.
Description of the device	The Aspira-aXA/-aXAY lenses are foldable, single-piece posterior-chamber IOLs with non-angulated cut-out (plate/frame) haptics and an aspheric (aberration-free/neutral), monofocal optic. The devices have an overall diameter of 11.0 mm and a body diameter of 7.0 mm. The IOLs are made of hydrophilic acrylic copolymer with UV-absorber, a water content of 26% at 35°C and a refractive index of 1.46. The yellow variant (with "Y") additionally contains a blue light filter that absorbs the high-energy portion of the light between 400 and 500 nm. The posterior optic surface has a 360° lens epithelial cell (LEC) barrier, the aim of which is to counteract the migration of the LEC onto the optic, thereby reducing the occurrence of posterior capsule opacification (PCO). The dioptr range covers -10.0 D to +30.0 D.  Schematic lateral and frontal view of the Aspira-aXA. The dimensions are given in mm. Orange arrows indicate markings for correct orientation

Mode of action / operating principle	The devices are designed to focus light passing through the cornea and pupil onto the retina, like the eye's natural lens. The aspheric optic design aims to provide higher-contrast imaging compared to a spherical design. Monofocal IOLs provide clear vision at one single focal point (usually far focus). Therefore, most patients will need postoperatively to rely on glasses to improve near and intermediate vision.
Materials or substances in contact with the patient's tissues	The hydrophilic acrylic copolymer consisting of 2-hydroxyethyl methacrylate (HEMA) and methyl methacrylate (MMA) with a UV absorber and with a yellow coloring agent (only for models followed by "Y") is permanently in contact with the ocular tissues.
Method of sterilisation and allowable number of uses	The device is sterile (steam sterilised) and designed for single use. Reprocessing of the implant is strictly prohibited.
Previous generations of the device(s)	The design of Aspira-aXA is derived from the long-standing and well-proven MC X11 ASP (HumanOptics). The MC X11 is not sold on the European market any longer.
Variants and/or configurations covered by this document	• Colour: • Colourless: device names ending with "-aXA" • Blue-filtering: device names ending with "-aXAY" • Packaging types: • Compact line (CL) • Safeloader (SL) Note: The Safeloader system is a packaging type and not an accessory.
Accessories which are intended to be used in combination with the device	Not applicable Note: Injector systems are independent medical devices which are intended to be used in combination with the device in question (i.e. these devices are not accessories).
Other devices and products which are intended to be used in combination with the device	Implantation of the device can be performed using either an injector system (for Safeloader and Compact line) or forceps (for Compact line only). Injector system tested for Safeloader and Compact line is the Medice AccuJect™ 2.0-1P (LP 604510 Rev. B (2v8), Medice AG, Switzerland). Note: Injectors designed for use with the HumanOptics Safeloader products feature a yellow plunger. When using a non-tested injector system for implantation of HumanOptics IOLs packaged in the Compact line, the user must refer to the specific instructions for use provided with the injector system to ensure that it is suitable to use with the IOL to be implanted.

4. Risks and warnings

Residual risks and undesirable effects	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Caution	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

Warnings	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Other relevant aspects of safety including a summary of any field safety corrective action (FSCA including FSN), if applicable	There have been no recalls since the devices have been launched. In July 2024, HumanOptics issued a Field Safety Notice affecting all devices in a Safeloader packaging manufactured in the period January 1, 2020 to June 28, 2024. During the standardized incoming goods inspection of returns, HumanOptics had observed very rare cases of damaged outer blister packaging after transportation. The FSN informed the users about the observation and reminded them to check the integrity of the sterile barrier system, as requested in the Instructions for Use, along with a visual instruction how to identify potential cracks. This damage is visible in the form of cracking and can be identified by users before use. It could theoretically cause a contamination of the surgical equipment or the intraocular lens. However, the hazard potential due to crack is considered to be very low. Statistical valuation showed that probability of occurrence of the additional risk of endophthalmitis due to the cracks is by several orders of magnitude lower than the general rate of endophthalmitis after cataract surgery. No patient injury has happened to the manufacturer's best knowledge. Meanwhile, the product packaging has been improved to eliminate this issue completely.

5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)

5.1. Summary of clinical data related to equivalent device

Not applicable. The conformity of the device was not assessed on the basis of equivalence.

5.2. Summary of clinical data from conducted investigations of the device before the CE-marking

Not applicable. All clinical investigations were performed after the first CE marking under the European Medical Device Directive 93/42/EEC (see section 5.3).

5.3. Summary of clinical data from other sources

Clinical data confirming the safety and performance of the Aspira-aXA/-aXAY have been derived from the following sources:

- Several post-market studies
- A comprehensive literature review using PubMed and other databases

Literature search strategies were designed to identify articles relevant to the device. Both favorable and unfavorable references were identified and summarized.

5.3.1. Post-market studies

The manufacturer HumanOptics has initiated several post-market studies with the Aspira-aXA(Y) as shown in Table 1. The studies/activities were labeled with acronyms considering: [device in focus]_[city]_[primary end point].

In the 5 post-market studies, a total of 273 eyes of 172 patients with senile cataract were implanted with the Aspira-aXA and followed for a time period of up to 3 years (see table below). Visual performance outcome was good and similar between study centers, resulting in a high level of patient satisfaction. The low degrees of postoperative lens rotation, tilt and decentration provide further evidence on the overall good safety of the IOL.

Data collected at the 1-year (n=100) and 1.5-year (n=113) and 3-year visits (n=53) also confirm the long-term safety and efficacy of the Aspira-aXA with stable visual acuity and no unanticipated side-effects.

Details of the study results are provided in the subsequent paragraphs.

Table 1: Overview of clinical studies/activities with the Aspira-aXA

Name of activity / Type of activity / Design	Aim	Status / Timelines
Asp-aXA_Ora_ND Post-market study RCT	Comparative study to confirm the long-term (3-year) performance and safety of the Aspira-aXA vs. the Aspira-aA after implantation in cataract patients. To assess the patient-related outcome based on patient satisfaction (including photopic disturbances).	Completed FU (3 years) completed: 03/2021
Asp-aXA_Lin_Rota Post-market study PSAS	A- All study centers: to assess the rotational stability of the Aspira-aXA. according to the requirements of ISO 11979-7:2018 B- Study centers Linz (Lin) and Glauchau (Gla): to evaluate the long-term (1.5-year) lens stability and capsule biocompatibility.	Completed FU (1.5 years) completed: 02/2021
Asp-aXA_Gla_Rota Post-market study PSAS		Completed FU (1.5 years) completed: 12/2020
Asp-aXA_Ber_Rota Post-market study PSAS		Completed FU (4 months) completed: 01/2020 (no 1.5 year visit planned)
Asp-aXA_Hal_VA Post-market study PSAS	To confirm the long-term (1-year) performance and safety of the lens and to provide further data on its rotational stability after implantation in cataract patients. To assess the patient-related outcome based on patient satisfaction (including photopic disturbances).	Completed FU (1 year) completed: 08/2021

FU=Follow-up; RCT=randomized controlled trial; PSAS=prospective single-arm study

Post-market study Asp-aXA_Ora_ND (Aspira-aXA/Aspira-aA – Oranienburg, Germany)

In the course of this study, the Aspira-aXA was implanted in a relatively large number of eyes (n=66). For comparison, a second group of patients was implanted with the well-proven HumanOptics Aspira-aA (n=65), which features a 6.0 mm optic diameter as opposed to 7.0 mm for the Aspira-aXA. Preoperatively, there was no statistically significant difference between the 2 groups in terms of age, axial length, anterior chamber depth, corneal astigmatism, white to white diameter, pupil diameter, subjective sphere and cylinder, spherical equivalent, IOL power and intraocular pressure. The follow-up covered the exceptionally long time span of 3 years.

At the 1-year and 3-year visits, median absolute deviation from target refraction (i.e. difference between expected refraction and actual refraction) was statistically significantly smaller in Aspira-aXA group [1-year; median 0.36 D (0.01 to 1.48); 3-year; median 0.36 D (0.01 to 1.60)] than in the Aspira-aA group [1-year; median 0.52 D (0.04 to 1.95); 3-year; median 0.60 D (0.03 to 1.81), (P<0.01)]. In terms of monocular visual acuity, both IOLs performed very well. At the 1-year follow-up, CDVA was 0.3 logMAR (20/40) or better in 100% of eyes, with 97% of the eyes achieving 0.1 logMAR (20/25) in the Aspira-aXA group and the Aspira-aA group. At the 3-year follow-up, 100% of the Aspira-aXA eyes had CDVA of 0.2 logMAR or better, with 98% achieving 0.1 logMAR or better. In the Aspira-aA group there were 98% achieving 0.2 logMAR or better and 98% 0.1 logMAR or better. Overall, there were no significant differences in the visual acuity outcomes between the Aspira-aXA and Aspira-aA IOLs at any postoperative time point and the values remained fairly constant over time. These data are also in line with those reported in previous studies showing similar visual acuity outcomes with the Aspira-aA and its toric version as well as other modern aspheric hydrophilic and hydrophobic monofocal IOLs.

Monocular examination of contrast sensitivity showed that the Aspira-aA performed overall slightly better than the Aspira-aXA. Nonetheless, and importantly, the scores of both IOLs were on average above age-related normal ranges measured with the same device.

There was no significant difference between the two groups in terms of occurrence of halos and glare sensitivity at the 1-year and 3-year follow-ups. Negative dysphotopsia was completely absent in the Aspira-aXA group at the 1-year visit, although the difference compared to the Aspira-aA group (4 eyes, i.e. 6.5% at least sometimes experienced the phenomena) did not reach significance (p=0.131). After 3 years, the occurrence of negative dysphotopsia decreased in the Aspira-aA group (only 1 patient “always” affected) while it was reported to occur “sometimes” by 2 patients in the Aspira-aXA group. Since these 2 patients developed slight capsule opacification over time and since negative dysphotopsia is acknowledged to occur just after cataract surgery, it is likely that the “unwanted image” perceived by these patients was a consequence of the development of a secondary cataract. At the 3-year visit, the Aspira-aXA also tended to outperform the Aspira-aA in terms of glare sensitivity at night (4 eyes/7.8% vs. 11 eyes/19.3%; p= 0.178).

In terms of patients’ satisfaction with their overall vision, the Aspira-aXA performed significantly better (p=0.024) than the Aspira-aA at the 1-year follow-up, with 90.5% of the patients being “very satisfied” (73.8%) or “satisfied” (16.7%) as opposed to 55.0%, and 42.5% in the Aspira-aA group, respectively. After 3 years of follow-up, the levels of satisfaction increased in both groups and the difference didn’t reach significance anymore (p=0.250). Nonetheless, the outcomes tended to be better with the Aspira-aXA with most patients being “very satisfied” (74.3%) or “satisfied” (20.0%) compared to 57.5% or 37.5% with the Aspira-aA, respectively.

No serious adverse events occurred during the course of the study. All IOLs remained well centered except two Aspira-aA where the optic was in front of the rhexis at the 1-year follow up. These two patients did not receive any additional treatment because they were “very satisfied” and “satisfied” with their overall vision. Capsulotomy with Nd:YAG laser was not necessary for any patient within the first postoperative year. After 3 years, one eye (1.5%) in the Aspira-aXA

group and two eyes (3.1%) in the Aspira-aA group required this treatment due to the development of PCO and associated reduction of visual acuity.

In conclusion, the Aspira-aXA performed well in terms of visual acuity, contrast sensitivity and patients' satisfaction over a long-term follow-up of 3 years and showed to be a safe and effective surgical treatment for cataracts. Overall, there were no clinically relevant differences in the visual function of patients implanted either with the Aspira-aXA or Aspira-aA IOLs, with however a trend for a higher level of patient satisfaction in the Aspira-aXA group.

Post-market study Asp-aXA_Lin_Rota (Aspira-aXA – Linz, Austria)

Wendelstein et al., 2021 [1]

In this study, 74 eyes of 42 cataract patients were implanted with the Aspira-aXA with the purpose to evaluate the lens rotational stability, tilt, decentration and overall visual performance up to 1.5 years postoperatively.

At the 1.5-year follow-up, median absolute deviation from target refraction was 0.69 D (range: 0.05 to 1.58) with no significant difference to the previous visits. Monocular CDVA was good with a mean value of -0.06 ± 0.10 logMAR. UDVA was not evaluated in this report since (16%) had a myopic target refraction.

IOL rotation was assessed at 1 day (1 D), 1 week (1 W), 1 month (1 M), 4 months (4 M) and 1.5 years (1.5 Y) postoperatively. The Aspira-aXA showed good rotational stability with a mean absolute rotation of $2.18^\circ \pm 2.23$ (median 1.42; 0.00 – 10.00) within one day postoperatively [End of Surgery (EoS) vs. 1 D] and decrease to less than 1.5° thereafter (1 D vs 1 W vs 1 M vs 4 M vs 1.5 Y). The only eye with an IOL rotation of 10° (EoS to 1 D) had a long AL (26.45 mm) and an immediate IOP spike after surgery, two conditions which are known to impact negatively postoperative IOL stability. Considering the performance criteria (EoS vs 4 M), the results show that 85% of the lenses had less than 5° rotation and 98% less than 10° rotation.

If rotational stability is particularly relevant for toric IOLs, decentration greater than 0.4 mm and tilt greater than 7.0 degrees can severely lower the visual performance of aberration-correcting aspheric IOLs [2], particularly with MIOs [3]. In this study (analysis based on Scheimpflug images), all IOLs had less than 6.0° of tilt (max = 5.84°) and the majority of them (93.4%) less than 0.4 mm of decentration (max = 0.51 mm). The values did not substantially change over time showing that the lens was very stable in the capsular bag. There was no significant correlation between axial length and tilt or decentration. While interpreting the results, one must consider that there are no baseline preoperative measurements to compare these values to. This is important because there is a strong correlation between the average tilt and decentration values of the crystalline lens and the IOL [4].

Patients' satisfaction with their overall vision was high. Between 92% and 95% were highly satisfied or satisfied when asked 1 month, 4 months and 1.5 years after surgery. Causes for dissatisfaction (where known) were postoperative ocular issues which were not linked to the IOL (Irvine-Gass syndrome and Sicca syndrome).

No serious adverse events occurred during the course of the study. Only PCO treatment was required and was performed on six eyes (8%). Two of them required the Nd:YAG laser treatment between the 4-month and 1.5-year visit, four interventions were performed either the day of the 1.5-year visit or scheduled for a later date. Visual acuity values of these four eyes before treatment were included in the final 1.5-year analysis.

In conclusion, this study shows that the Aspira-aXA performs well in terms of visual acuity, patient satisfaction, safety and IOL capsular stability.

Post-market study Asp-aXA_Gla_Rota (Aspira-aXA – Glauchau, Germany)

Pilger et al., 2021 [5]

Schrecker et al., 2021 [6]

This study follows the same protocol as the study Asp-aXA_Lin_Rota. The difference is that the 1-day follow-up visits were not possible in this center. Fifty-five eyes of 39 cataract patients were implanted with the Aspira-aXA.

At the 1.5-year visit, 96% of the eyes were within ± 1.0 D and 90% within ± 0.50 D deviation from target refraction. Monocular UDVA and CDVA were good with mean values of 0.05 ± 0.13 logMAR and -0.01 ± 0.05 logMAR, respectively. UDVA significantly improved over the complete follow-up time ($p=0.02$), while differences in CDVA ($p=0.40$) were not statistically significant.

Mean IOL rotation was $2.30^\circ \pm 2.83$ (median 1.93; 0.03 – 13.37) within one week postoperatively (EoS to 1 W) and decreased thereafter. Considering the performance criteria (EoS vs 4 M), the results show that 91% of the cases had rotation less than 5° , 97% of the cases rotated less than 10° and only one eye (3%) rotated more than 10° (11.60°).

OCT-based analysis of IOL tilt and decentration showed mean values of $\sim 5^\circ$ and ~ 0.15 mm, respectively, which is in the range previously described for other IOLs (using the same method). There were no significant differences between the measurements made 1 and 4 months after surgery. Notably, the average tilt determined in this study was considerably higher than in the study Asp-aXA_Lin_Rota. Since the protocols and other results of both studies were very similar, the difference can very likely be explained by the fact that two different methods (OCT vs. Scheimpflug imaging) and two different reference axes were used. As mentioned in the summary of the study Asp-aXA_Lin_Rota, while interpreting the results one must consider that the crystalline lens itself is usually tilted and slightly decentered, which affects the postoperative position of the IOL.

At the 1.5-year visit, 97% of patients were highly satisfied or satisfied with their overall vision. One patient (3%) was moderately satisfied which was most likely due to slight fibrosis after 1.5 years.

No serious adverse events occurred during the course of the study. Intraoperatively, a capsular tension ring had to be implanted in one eye due to a capsular rupture. All IOLs remained well centred. Nd:YAG laser treatment due to PCO was required for 3 eyes during the complete follow-up period.

In conclusion, the Aspira-aXA performed well in terms of visual acuity, patient satisfaction, safety and IOL capsular stability.

Post-market study Asp-aXA_Ber_Rota (Aspira-aXA – Charité Berlin, Germany)

Pilger et al., 2021 [5]

The main objective of this study was to evaluate the postoperative IOL rotational stability and visual acuity up to 4-month. The study followed the same protocol as the study Asp-aXA_Lin_Rota. Thirty-six eyes of 19 cataract patients were implanted with the Aspira-aXA.

At the 4-month visit, median absolute deviation from target refraction was 0.21 D with no significant changes over time. Results show good visual acuity performance with mean monocular UDVA and CDVA of -0.03 ± 0.13 logMAR and -0.15 ± 0.06 logMAR, respectively (4-months postoperatively).

The rotational stability was high. Mean IOL rotation was $1.84^\circ \pm 1.42$ (median 1.84; 0.07 – 4.26) within one day postoperatively (EoS vs 1 D). Between all post-operative visits mean IOL rotation was less than 1.0° . None of the IOLs rotated more than 10° .

Patients' satisfaction with their overall vision: at the 4-month visit, all patients were satisfied or very satisfied with their overall vision.

No serious adverse events occurred during the course of the study. All IOLs remained well centred. Nd:YAG laser treatment was not required for any patient during the course of the study.

In conclusion, the Aspira-aXA performed well in terms of visual acuity, patient satisfaction, safety and IOL rotational stability.

Post-market study Asp-aXA_Hal_VA (Aspira-aXA – Halle, Germany)

Pilger et al., 2021 [5]

The main purpose of this study was to assess the clinical performance of the Aspira-aXA up to 1 year, with a focus on the occurrence of positive and negative dysphotopsia. In this study, 42 eyes of 25 cataract patients were implanted with the Aspira-aXA.

At the 1-year follow-up, median absolute deviation from target refraction (Haigis formula) was 0.21 D with no significant changes over time. The Aspira-aXA performed well in terms of monocular UDVA and CDVA with mean values of 0.05 ± 0.10 logMAR and -0.02 ± 0.05 logMAR, respectively. CDVA significantly improved between the 3-month and 1-year ($p=0.02$), whereas differences in UDVA ($p=0.15$) were not statistically significant.

The rotational stability was also good with a mean value of $2.79^\circ \pm 2.98$ (median 1.88; 0.06 to 10.77°) within one week postoperatively (EoS to 1 W) and remains stable thereafter, with a mean IOL rotation of less than 2.15° between all post-operative visits. After 1 year (vs EoS), IOL rotation was less than 5° in 93% of the available eyes, and less than 10° in 100% of the cases.

However, there were three lenses with rotation of $\sim 10^\circ$ in the early postoperative period (EoS to 1 W), for which rotation could not be analysed at the 1-year follow-up due to insufficient image quality. It is possible that a poor friction between the haptic and the capsular bag due to an incomplete ophthalmic viscoelastic clearance influenced lens rotation within the early postoperative period.

Halos were absent for most patients according to the results of the questionnaire. Glare sensitivity was experienced by a 66.7% of patients during the day but was less pronounced during night-time (25.0%). Only one patient experienced negative dysphotopsia sometimes. These data are consistent with those reported in the study Asp-aXA_Ora_ND showing minor or no postoperative unwanted photic phenomena postoperatively.

Most patients (95.8%) were satisfied or very satisfied with their overall vision. The only patient who was unsatisfied could be explained by a post-operative hyperopic refractive error ($SE=+1D$) in one eye, resulting in a poor uncorrected visual acuity. 83.3% reported to be totally free from spectacles for far vision. Interestingly, although the lens has a monofocal design, 58.3% of the patients reported to “never” or “sometimes” use spectacles for intermediate vision.

No serious adverse events occurred during the course of the study. None of the eyes required Nd:YAG laser treatment during the complete follow-up period.

In conclusion, the Aspira-aXA performed well in terms of visual acuity, patient satisfaction and IOL rotational stability.

5.3.2. Literature review

The 4 articles identified through systematic literature review are of various usefulness to assess the safety and performance of the device under consideration.

The most relevant findings are the observations by Bonsemeyer et al. [7], which indicate that the aXA IOL platform may be able to reduce the frequency and extent of dysphotopsias compared to a conventional IOL design with 6 mm optic size, confirming the intended benefit. These findings are consistent with the results of Rozendal et al. [8]. They exchanged the IOLs of patients who suffered from positive and/or negative dysphotopsia with the Aspira-aXA, and all but one patient reported a complete or partial resolution of their complaints. Taken together, the results from these two reports strongly indicate that the design of the Aspira-aXA with a 7 mm optic and plate haptics is associated with a lower rate and severity of positive and negative dysphotopsia compared with conventional monofocal lenses with a 6 mm optic and C-loop haptics.

The two publications by Borkenstein et al. [9, 10] suggest, based on a small number of eyes, that the Aspira-aXA/-aXAY could be a particularly useful option for patients with retinal disease (potentially combined with high myopia), i.e. for eyes which need monitoring of disease progression and may require posterior-segment surgery later on. Although intra- and postoperative complications are frequent in such patients, there were none in their case report/series and the vision significantly improved after IOL implantation, which was associated with a high patient satisfaction.

5.4. Overall summary of the clinical performance and safety

IOLs are intended to replace the focusing power of the natural lens thereby restoring the patient's vision. Accordingly, the main criterion for the success of cataract surgery is the long-term visual result. Since monofocal, non-accommodative lenses like the Aspira-aXA lenses focus the light at a single, fixed distance it is expected that the patients will still require glasses for activities at certain distances, e.g. reading or screen handling. Therefore, the main performance indicators for the devices under consideration are the obtained visual acuity and the reduction of pre-existing astigmatism for the toric models (when available). Key safety indicators are the rate of complications and, in the case of the toric devices, the rate of IOL rotations.

PERFORMANCE ANALYSIS

Based on the results of the post-market studies, in which 172 cataract patients (273 eyes) were enrolled and implanted with the Aspira-aXA, a performance analysis was performed and included all eyes available during a period ranging from 3 months to 3 years post-operatively. The primary performance endpoint was met with all eyes, except one at the 4-month follow-up (251/252; 99.6%), achieving CDVA of 0.3 logMAR or better, exceeding the primary safety and performance endpoints (SPE) rates as defined in ISO 11979-7:2018-08 for posterior-chamber lenses (92.5% overall). **Key performance CDVA outcomes are presented in Table 2.** The visual performance remained stable over the investigated study period of 3 years, supporting the good performance of the Aspira-aXA, with even 98% of the eyes achieving a CDVA of 0.1 logMAR or better at the 3-years follow-up.

Table 2: Evaluation of performance outcomes according to DIN EN ISO 11979-7: 2018-08 [11]

Endpoint	Acceptance criterion	Clinical outcome	Data source
Monocular CDVA (% of eyes reaching 0.3 logMAR or better)	≥ 92.5% of the eyes should achieve CDVA 0.3 logMAR or better [SPE rate]		Post-market studies:
		100% [-0.20 to 0.15] (n=59) at 3-month 100% [-0.18 to 0.15] (n=59) at 1-year 100% [-0.18 to 0.15] (n=53) at 3-year	Asp-aXA_Ora_ND
		98.5% [-0.30 to 0.60] (n=66) at 4-month 100% [-0.20 to 0.20] (n=63) at 1.5-years	Asp-aXA_Lin_Rota
		100% [-0.10 to 0.22] (n=54) at 4-month 100% [-0.10 to 0.10] (n=54) at 1.5-years	Asp-aXA_Gla_Rota
		100% [-0.20 to 0.00] (n=34) at 4-month	Asp-aXA_Ber_Rota
		100% [-0.10 to 0.22] (n=39) at 3-month 100% [-0.10 to 0.10] (n=40) at 1-year	Asp-aXA_Hal_VA

The scientific publications regarding Aspira-aXA/-aXAY add the following information. Even in eyes with comorbidities such as retinal pathologies, the postoperative visual acuity was improved substantially compared to preoperative levels [9, 10].

Two datasets allow a comparison with other monofocal devices: In the above-mentioned RCT Asp-aXA_Ora_ND, the Aspira-aXA provided a UDVA and CDVA which was very similar to the comparator Aspira-aA, a monofocal IOL from HumanOptics with a 6 mm optic diameter. The case series presented by Rozendal et al., 2021 [8] in which a standard IOL was explanted due to negative and/or positive dysphotopia and replaced by the Aspira-aXA, allows to compare visual acuities within the same eye i.e. Aspira-aXA *versus* competitor standard monofocal IOLs. Preoperatively, decimal CDVA was 1.22 ± 0.2 (-0.09 logMAR) compared to 1.21 ± 0.2 (-0.08 logMAR) postoperatively. Although the names of the explanted IOLs were not provided, this finding indicates that the Aspira-aXA does not substantially differ from standard competitor monofocal IOLs in terms of visual acuity.

Taken together, these data show that the Aspira-aXA/-aXAY fulfils the pre-defined acceptance criteria and achieves the performance intended by the manufacturer as required by the MDR Annex I, when used according to the intended use.

SAFETY ANALYSIS

Several complications and undesirable effects have been reported in the literature for foldable posterior-chamber IOLs but the incidence is overall very small, particularly, if one takes into account that eligible patients often have several ocular comorbidities. In most cases these side-effects are related to the surgery and not device-related. Regarding the Aspira-aXA/-aXAY, the reported outcomes show only a small percentage of complications and undesirable effects which have been already described in the literature.

Based on the results of the post-market studies for the Aspira-aXA/-aXAY (Table 1) a safety analysis was performed. The primary safety endpoint was met as the cumulative and persistent adverse events (complication) rates were within limits compared to the historical SPE rates reported in ISO 11979-7:2018. **Key safety outcomes are presented below in Table 3.**

Table 3: Calculation of complication rates according to DIN EN ISO 11979-7: 2018-08 [11]

Endpoint	Acceptance criterion (SPE rate) %	Threshold rate % for n = 100	Threshold rate % for n = 300	Clinical outcomes % for n = 273	Data source
Cumulative					Post-market studies (n = 273): Asp-aXA_Ora_ND Asp-aXA_Lin_Rota Asp-aXA_Gla_Rota Asp-aXA_Ber_Rota Asp-aXA_Hal_VA
Cystoid macular oedema	3.0	8.9	6.0	2.93	
Hypopyon	0.3	3.0	1.8	0.00	
Endophthalmitis	0.1	3.0	1.0	0.00	
Lens dislocated from posterior chamber	0.1	3.0	1.0	0.00	
Pupillary block	0.1	3.0	1.0	0.00	
Retinal detachment	0.3	3.0	1.8	0.00	
Secondary surgical intervention	0.8	4.2	2.6	0.00	

Persistent				
Corneal stroma oedema	0.3	3.0	1.8	0.00
Cystoid macular oedema	0.5	4.2	2.2	1.10 *
Iritis	0.3	3.0	1.8	0.00
Raised IOP requiring treatment	0.4	4.2	1.8	0.00

* value is not statistically significantly lower than SPE rate (threshold for a study with n=300 is not exceeded)

Regarding complaints related to intra- and postoperative complications, few have been reported for the Aspira-aXA/-aXAY and neither serious incidents nor reportable serious adverse events have been recorded.

Clinical studies with HumanOptics IOLs have been conducted using both the Compact Line as well as the Safeloader-packaged IOLs and confirmed the performance and safety of the tested injectors. However, common injector-associated complications, such as trapped haptics or misconfiguration of the IOL haptics after injector release, have also been observed in a few cases. These were typically associated with handling errors.

Since 2003 HumanOptics has placed several million hydrophilic acrylic IOLs on the market. Over this period of more than 20 years no feedback has been received relating to material degradation other than very rare cases (< 5 ppm/year) of suspected IOL opacification (calcification). The latter is caused by crystallisation of calcium on the IOL surface, usually resulting from environmental circumstances, i.e. it does not constitute a degradation due to aging processes of IOLs manufactured by HumanOptics. The available information confirms that HumanOptics hydrophilic acrylic IOLs are stable over a period of at least 20 years (based on simulated aging tests according to ISO 11979-5). It can be expected that the device is stable indefinitely from the implantation date over the lifetime of the patient.

Based on the above-mentioned results and the outcomes in eyes with comorbidities such as retinal pathologies [9, 10], it can be concluded that the Aspira-aXA/-aXAY is safe for implantation in the capsular bag for the visual correction of aphakia, and no unanticipated risks associated with Aspira-aXA/aXAY IOLs were identified. In summary, one can conclude that the device under consideration represents the current state of the art and fulfils the pre-defined acceptance criteria and the requirements regarding safety as set out by the MDR Annex I.

BENEFIT-RISK ANALYSIS

The main **benefits** of cataract surgery and IOL implantation in general are improved visual function and quality of life. Important aspects involved are:

- Improved vision (visual acuity, contrast sensitivity, depth perception, etc.)
- Reduction in falls, fractures and vehicle crashes
- Reduced eyeglass dependence (through correction of refractive errors and, in case of TIOs, of astigmatism)
- Improvement of quality of life (activity, confidence, mood)
- Reduction of anxiety and depression

In terms of improved vision, the different IOL types currently available on the European market can address single or a combination of vision problems. For instance, **monofocal IOLs** are designed to improve vision at just one set distance (generally far vision); **multifocal IOLs** are designed to additionally improve vision at near and/or intermediate distances and **toric IOLs** specifically correct pre-existing corneal astigmatism (they are available as monofocal toric or multifocal toric IOLs).

Three types of intended clinical benefits can be attributed to the Aspira-aXA/-aXAY.

The patient benefit “**restored vision**” is based on achieving good CDVA after IOL implantation, which has been demonstrated in the previous chapters. Furthermore, the lens showed good contrast sensitivity outcomes along with high levels of patients’ satisfaction with more than 90% of them being “very satisfied” or “satisfied” with their overall vision. A previous study has also shown that the blue light filter added in HumanOptics yellow IOLs has no negative impact on visual acuity, contrast sensitivity under photopic and mesopic conditions with and without glare, color perception, patient satisfaction, and rate of complications [12].

There is clinical evidence that the large optic diameter (7 mm) of the Aspira-aXA/-aXAY may provide additional benefits. First, it is expected to **reduce the occurrence of visually impairing optical phenomena such as dysphotopsia**. In fact, the article by Bonsemeyer et al. [7] (based on an extension of the study Asp-aXA_Ora_ND) reported a significantly lower prevalence of negative dysphotopsia 1 month after surgery for the Aspira-aXA compared to a lens with 6 mm optic size (Aspira-aA). Similar results were obtained for the perceived extent and frequency of the photopic phenomena, although significance was mostly not reached. A limitation of this study is the fact that the two IOL models in the comparison not only differed in optic size but also in the haptics design, which may have influenced dysphotopsias. Additional support for the positive effect of a large optic on negative dysphotopsia comes from several instances of customer feedback referring to patients who suffered from dysphotopsia after cataract surgery. Those disturbances disappeared, when the original IOL (from varying competitors) was explanted and replaced with an Aspira-aXA (data not shown). This positive effect on the resolution of dysphotopsia was also reported by Rozendal et al. [8] at the ESCRS congress 2021. In 12 patients (15 eyes) complaining about negative and/or positive dysphotopsia, they exchanged the original IOL bag-to-bag with the Aspira-aXA. Complete (n=11/15; 73.3%) or partial (n=3/15, 20.0%) resolution of positive and negative dysphotopsia after IOL exchange with the Aspira-aXA was achieved in all but one patient [this patient reported complete resolution after more than 1 year]. Taken together, the results from these two reports strongly indicate that the design of the Aspira-aXA with a 7 mm optic and plate haptics is associated with a lower rate and severity of positive and negative dysphotopsia compared with conventional monofocal lenses with a 6 mm optic and C-loop haptics.

A second benefit of the large optic diameter (and concomitantly a larger capsulorhexis) is an **extended view of the fundus**. With the aXA IOL platform, the view of the retina is considered to be better than that with standard-sized 6.0-mm IOLs [9, 10]. Although this might not be of very much interest to the surgeon focused on the anterior segment, it can be very relevant in the presence of retinal diseases. The findings from a small case series [10] and a case report [9] confirmed that the Aspira-aXA can be very useful in this context. Many patients suffering from central visual loss rely on their peripheral vision and even adapt to using a retinal location that is not affected by their disease. So it is important to keep this region as healthy as possible. The enhanced view of the fundus provided by the aXA platform helps the treating physician to monitor the progression of the retinal disease and might also facilitate further treatments (e.g. anti-VEGF injections or laser photocoagulation) to keep the peripheral region as healthy as possible [10]. Furthermore, the wide fundus view might be beneficial for all types of posterior-segment surgeries [10]. Thus, the device may also be particularly suitable for patients (e.g. with high myopia) who have an elevated risk to develop retinal or vitreous diseases at some point in their life.

It is understood that the extent of benefits depends on the degree of patients' preoperative discomfort and ocular comorbidity which may limit the performance of the device as, for instance, in patients with retinal diseases. In the case of patients suffering from severe cataract, who are considerably visually impaired, the implantation of IOLs can mean a huge relief and improvement in quality of life.

Regarding **risks**, few complications have been reported for the Aspira-aXA/-aXAY, which is in line with the low frequency of complications reported for cataract surgery in general. Nonetheless, complications can always occur but only in extremely rare cases complications may lead to serious and permanent impairment of the patient. Some of the most relevant risks are ocular inflammation and cystoid macular edema.

As mentioned earlier, postoperative complications are rarely device-related but rather depend on the preoperative ocular comorbidity, which means that IOLs should be implanted into risky patients only after careful assessment of the benefit/risk ratio. Furthermore, the success of cataract surgery is influenced by the surgical technique, equipment used, and the surgeon's experience. In addition, there is a learning curve when performing implantations with a IOLs with novel designs. Nevertheless, the implantation of IOLs has to be considered as a routine procedure for the target user group, i.e. ophthalmic surgeons. Since the essential aspects of the handling and implantation process are adequately described in the IFU, a special training for the Aspira-aXA/-aXAY is not required. However, device handling errors cannot be completely excluded but appear to occur very rarely.

In summary, a large body of evidence proves that the Aspira-aXA/-aXAY is a reliable product and that the remaining risks are well-controlled.

A positive benefits-versus-risks profile can be assigned to the Aspira-aXA/-aXAY, particularly considering the current knowledge and state of the art in the medical fields concerned. The intended clinical benefits are fulfilled. The probable benefits outweigh the probable risks for the intended use of the device.

5.5. Ongoing or planned post-market clinical follow-up

The manufacturer has initiated an additional post-market study with the goal of assessing the clinical outcomes with the Aspira-aXA after cataract surgery combined with pars plana vitrectomy. The background for this is that cataract and vitreoretinal diseases often occur together at an advanced age. Here, the Aspira-aXA with its large optic could provide a better view of the retinal periphery and could also reduce the risk of postoperative "iris capture".

6. Possible diagnostic or therapeutic alternatives

At present, the only well-tolerated treatment available to replace the focusing power of the natural lens after surgical extraction is the implantation of an IOL. The preferred placement of an IOL is in the capsular bag. Mild/early cataracts may also be managed by prescribing new eyeglasses to correct a lens-induced change in refractive errors. High-power eye glasses (still in use in developing countries) and contact lenses (only for high myopic eyes) are alternative treatments for aphakia but do not achieve visual results as good as IOLs.

Due to the large optic diameter (7.0 mm) of the Aspira-aXA/-aXAY, its implantation in the capsular bag might be more challenging compared to standard IOLs, which typically have a 6.0 mm optic. Otherwise, there are currently no known disadvantages in comparison to other foldable monofocal posterior-chamber IOLs. In contrast, a large optic diameter might provide certain advantages such as prevention of post-operative dysphotopsia and better view of the retinal periphery.

7. Suggested profile and training for users

IOLs must be handled by health professionals and implanted by ophthalmic surgeons. A high level of surgical skill is required for IOL implantation. The surgeon should have observed and/or assisted in a large number of implantations and successfully completed one or more courses on IOL implantation before attempting a surgery. It is also strongly recommended that the surgeon and his/her team receive a training on the specific product and accessories by the respective sales representative before the first implantation of the device. Application guides (PDF) may be available on request. Product application videos can be found on www.humanoptics.com in the “application” section for the respective device.

If you have further questions on the usage of the device, please contact your local sales representative or the manufacturer (application@humanoptics.com).

8. Reference to any harmonised standards and CS applied

Common specification(s) to comply with, if applicable:
N/A
Harmonised standards to apply, if applicable:
DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
DIN EN ISO 14971:2020-07 Medical devices – Application of risk management to medical devices
Relevant adopted monographs of the European Pharmacopoeia
Not applicable

9. References

- [1] Wendelstein J, Laubichler P, Fischinger I, Mariacher S, Beka S, Mursch-Edlmayr A, et al. Rotational Stability, Tilt and Decentration of a New IOL with a 7.0 mm Optic. *Curr Eye Res.* 2021;46(11):1673-1680.
- [2] Holladay JT, Piers PA, Koranyi G, van der Mooren M, Norrby NE. A new intraocular lens design to reduce spherical aberration of pseudophakic eyes. *J Refract Surg.* 2002;18(6):683-691.
- [3] Sieburth R, Chen M. Intraocular lens correction of presbyopia. *Taiwan J Ophthalmol.* 2019;9(1):4-17.
- [4] Kimura S, Morizane Y, Shiode Y, Hirano M, Doi S, Toshima S, et al. Assessment of tilt and decentration of crystalline lens and intraocular lens relative to the corneal topographic axis using anterior segment optical coherence tomography. *PLoS One.* 2017;12(9):e0184066.
- [5] Pilger D, Bertelmann E, Brockmann T, Duncker T, Duncker GIW, Schrecker J. Postoperative Lens Rotation of a 7.0 mm Optic IOL with Plate Haptics. *J Ophthalmic Vis Res.* 2021;5(1):014-020.
- [6] Schrecker J, Seitz B, Langenbucher A. Performance einer neuen 7-mm-Intraokularlinse mit Nachbeobachtung über 1,5 Jahre [Performance of a new 7 mm intraocular lens with follow-up over 1.5 years]. *Ophthalmologe.* 2022;119(4):367-373.
- [7] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. *J Cataract Refract Surg.* 2022;48(1):75-82.
- [8] Rozendal LRW, Van Vught L, Beenakker JWM, Luyten GPM. Bag-to-bag intraocular lens exchange with a wide optic intraocular lens for treatment of positive and negative dysphotopsia. 39th Congress of the ESCRS 2021.
- [9] Borkenstein AF, Borkenstein EM. Creating Hybrid Monovision with 7.0 mm XL Optic and High-Add AMD Intraocular Lenses (XL-MAGS) in a Patient with Retinitis Pigmentosa. *Case Rep Ophthalmol.* 2019;10(3):304-311.
- [10] Borkenstein AF, Borkenstein EM. Efficacy of Large Optic Intraocular Lenses in Myopic Eyes with Posterior Segment Pathology. *Ophthalmol Ther.* 2022;11(1):443-452.
- [11] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia
- [12] Kühle M. Comparison of visual function with aspheric yellow, aspheric clear and spherical clear intraocular lenses. *J Emmetropia* 2013;4:123-130.

Part 2 – SSCP intended for patients

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.


The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your ophthalmologist in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions For Use to provide information on the safe use of the device.

[This SSCP is reviewed at least annually and will be updated if needed.](#)

1. Device identification and general information

Device trade name(s)	- Aspira-aXA - Aspira-aXAY
Manufacturer's name and address	HumanOptics Holding AG Spardorfer Str. 150 91054 Erlangen Germany
Basic UDI-DI(s)	Aspira-aXA/-aXAY: 04049154_PC_M2_H4_O2_CJ
Year when the first certificate (CE) was issued covering the device	- Aspira-aXA: 2016 - Aspira-aXAY: 2020



2. Intended use of the device

Intended purpose	The device is a tiny intraocular lens (IOL) intended to be implanted in the capsular bag ¹ of your eye to replace the natural lens (crystalline lens).
Indication(s) and intended patient groups	The IOL is used to replace your natural lens after removal by an ophthalmic surgeon. Usually, the removal of your natural lens is because your lens has become cloudy with the age. This is called "cataract". Cataract affects your vision.  Implantation of an IOL is indicated for adult patients only. Note: there are no clinical data available for children, pregnant or breastfeeding women, and immunocompromised persons in connection with the implantation of intraocular lenses from HumanOptics.




¹ **Capsular bag:** The capsule which the natural lens sits in.

	Depending on your eye anatomy, life-style and preferences, different types of IOLs can be offered to you by your ophthalmologist (see section 3 “mode of action/operation principle”).
Contraindications	There are general contraindications for ocular surgery but no additional contraindications which are specific to the implantation of HumanOptics IOLs.

3. Device description

Description of the device	The Aspira-aXA/-aXAY is a tiny artificial lens made of plastic, called an intraocular lens (IOL), which is implanted into the eye to replace the natural lens. In most cases, the IOL will be used to replace the clouded lens (cataract) of a person. The IOLs have a diameter of just a few millimetres and fit perfectly into the eye due to their shape and the flexibility of their material. IOLs are composed of an “optic” (the round part of the IOL) and the “haptics”. The optic has a dioptric power that permits the focusing of light onto the retina and the haptics which are attached to the edge of the optic will help maintain the IOL in place in the eye. The Aspira-aXA/-aXAY is available in colourless or yellow versions. The yellow variant (with “-Y”) absorbs a portion of the blue light, which may help protect the retina² over the long term.  Orange arrows indicate markings for correct orientation Schematic frontal and lateral view of the Aspira-aXA. The IOLs exist also in yellow (Aspira-aXAY). The dimensions are given in mm.
Mode of action / operating principle	If you have been diagnosed with a cataract, it is because your natural lens has become cloudy and this affects your vision. Most commonly, cataract occurs in older people, and usually develops gradually. Cataract removal is one of the most common operations worldwide. During the procedure, usually an outpatient operation, the natural eye lens is removed and replaced by an IOL. The short procedure is usually painless and causes very little discomfort for the patient (see also chapter 5.1). Step 1 – You have been diagnosed with a cataract 

² **Retina:** The innermost, light-sensitive layer of tissue of the eye that transmits that information to the brain via the optic nerve

	Step 2 – The cataract is break it into tiny pieces using ultrasonic waves or a laser, which are then suctioned out of the eye with a vacuum. The procedure takes usually only a few minutes and is performed under local anesthesia  Step 3 – After removal of the cataract, an IOL is inserted with an injector into the capsular bag of your eye to restore your vision  Step 4 – The IOL is now in place in the capsular bag of your eye  There is a variety of IOL types to choose from — all with different benefits. Choosing the right IOL can be difficult and it depends on several factors including your eye condition and your lifestyle. The different options should be discussed with your ophthalmologist. Among the different IOL types, the monofocal IOLs are currently, the most widely used lenses. Monofocal IOLs like the Aspira-aXA are designed to provide clear vision at a fixed distance (usually far distance). Therefore, it is expected that your far vision will significantly improve postoperatively but glasses or contact lenses might be needed for activities at near (like reading) and intermediate distances (like computer work). The Aspira-aXA has a relatively large optic (7 mm) compared to the majority of other IOLs (6 mm), meaning that patients with large pupils for instance might perceive less disturbing optical phenomena such as glare and halos around point sources of light (called also positive dysphotopsia). The Aspira-aXA might also help preventing shadows in the peripheral field of vision (called also negative dysphotopsia).
Materials or substances in contact with the patient's tissues	All IOLs manufactured by HumanOptics are made of biocompatible plastic (hydrophilic acrylic material) with an integrated UV absorber. Models with names ending with "Y" additionally contain a yellow coloring agent intended to filter the blue light.
Information about medicinal substances in the device, if any	Not applicable.
Accessories	None which are relevant for patients

4. Risks and warnings

Contact your ophthalmologist if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your ophthalmologist if needed.

How potential risks have been controlled or managed	Potential risks have been identified and controlled according to international standards that have been developed by experts from many countries. The corresponding document(s) comprising rules, guidelines and processes are approved and published by a globally recognized organisation.
Remaining risks and undesirable effects	As with any surgical intervention, routine cataract surgery and associated IOL implantation are not entirely risk-free. However, this procedure is among the most common and most successful ones performed nowadays in the world. There are no known risks which are specific to the Aspira-aXA/-aXAY and which have not been previously reported for similar IOLs. The following side effects may occur following the IOL implantation but they normally disappear shortly after the surgery: • Foreign body sensation caused by the small incision through which the lens was inserted • Halos, which can form around sources of light, particularly at dusk or in darkness • You may experience increased sensitivity to light in the first few weeks after surgery The following symptoms may indicate complications associated with an IOL implantation. If these occur, you should consult your ophthalmologist: • Sudden or gradual deterioration in vision / blurred vision • Severe redness of the eye • Eye pain Months or years after IOL implantation, your vision might deteriorate and become blurry or hazy. This might be due to posterior capsule opacification (PCO), which occurs when a cloudy layer of scar tissue forms behind the IOL. This condition is not dangerous and can be treated by performing an opening in the capsule (or bag) that contains your IOL. Opening the capsule allows light to enter the eye and improve the clarity of the vision, as long as the rest of the eye is healthy. This procedure is called “Nd:YAG laser capsulotomy”. It is fast (a few minutes), painless and is an easy way for restoring your vision. There is no other treatment for PCO apart from this procedure, and the condition will not improve by itself.
Warnings and precautions	Patients with certain eye diseases and/or malformations have higher-than-normal risks to develop complications during and after IOL implantation. Before the surgery, your ophthalmologist will perform an assessment of your eye and will discuss with you anything to do with your operation. This includes the risks and benefits of surgery and how long it will take to fully recover after the surgery. Difficulties during surgery may also require choosing another lens than planned.
Summary of any field safety corrective action (FSCA including FSN), if applicable	No devices have been recalled since the Aspira-aXA/-aXAY have been sold on the market.

	In July 2024, HumanOptics issued a Field Safety Notice (FSN)³ for healthcare professionals. The notice affected all devices in a certain packaging which had been manufactured in the period January 1, 2020 to June 28, 2024. The manufacturer had observed very rare cases of damaged outer blister packaging after transportation. The associated hazard potential is considered to be very low and no patient injury has happened to the manufacturer's knowledge. The FSN served to inform the users of this possible risk and to tell them how a damaged packaging can be recognized. No specific action is needed in the case of patients who have been treated with potentially affected products. Meanwhile, the product packaging has been improved.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1. Clinical background of the device

Cataract, i.e. the clouding of the natural lens of the eye, is a very common condition and most commonly affect adults as a result of ageing. Cataract causes deterioration of your vision and usually gets slowly worse over time. Surgery to replace the cloudy lens is the only way to improve your vision. Implantation of an IOL is the only well-tolerated treatment available to replace the cloudy lens. This procedure that usually takes 30 minutes is typically associated with rapid visual recovery, good visual outcomes and minimal complications.

After you and your ophthalmologist have decided that you should have your cataract removed, the length of your eye and the shape of your cornea will be measured. These measurements will help your surgeon choose the best IOL for you. You will also discuss the various lens options available to you. Further information about cataract surgery can be found under:

<https://www.humanoptics.com/en/patienten/help-with-cataracts/>

After implantation, the IOL will normally remain in your eye for the rest of your life, provided there is no medical contraindication. The product lifetime of HumanOptics IOLs is proven for twenty years. Regular ophthalmological check-ups are recommended, as for patients with natural lenses.

5.2. The clinical evidence for the CE-marking

The letters 'CE' which appear on products signify that products sold in the European Economic Area [EEA] have been assessed to meet high safety, health, and environmental protection requirements.

Clinical evidence for CE-marking IOLs should meet the requirements defined by international standards and is based on clinical data with the device, scientific literature and feedback from the market.

Clinical evidence with the Aspira-aXA/-aXAY was based on the data from 5 post-market studies. Patients were examined up to 3 years after IOL implantation. All studies showed that the IOLs achieved the desired improvement in vision after surgery. In each of the studies, more than 90% of the patients were satisfied or very satisfied with their overall vision. Those who were not satisfied usually experienced complications of cataract surgery or suffered from health

³ A field safety notice (FSN) is an important means of communicating safety-related information to users.

conditions which interfered with their vision. Overall, there was no safety concern regarding the implantation of the lens in patients undergoing cataract surgery.

In summary the Aspira-aXA/-aXAY performs well in terms of visual acuity, patient satisfaction and safety. Its stability in the eye is good and comparable to other IOLs on the market. It shows to be safe and effective option for the surgical treatment of cataract.

5.3. Safety

The implantation of the CE-marked Aspira-aXA/-aXAY IOL considerably improves vision in cataract patients. The CE marking for IOLs is mandatory and means that the manufacturer has verified that the products meet European safety, health or environmental requirements.

The side effects of IOLs are well-known and are similar for all IOLs on the market. They generally occur rarely and are not very harmful (see chapter 4 “Remaining risks and undesirable effects”). Therefore, the benefit-risk ratio of the device is very positive. The manufacturer continuously analyses complaints, other market feedback and new scientific publications. The results of these activities have been favourable, which further supports the safety of the device, besides the other data which were evaluated for the CE marking.

The manufacturer is currently conducting an additional clinical study which focuses on the clinical outcomes of the Aspira-aXA IOL after implantation in cataract patients who additionally have vitreoretinal diseases; conditions which affect the structures in the eye called the retina and the vitreous⁴.

6. Possible diagnostic or therapeutic alternatives

When considering alternative treatments, it is recommended to contact your ophthalmologist who can take into account your individual situation.

Currently, there are no drugs which can stop or reverse the clouding of the natural lens of the eye. In the case of a pronounced cataract, the clouded lens must therefore be removed and replaced by an IOL. The preferred placement of an IOL is in the capsular bag. If you have a mild/early cataract, your ophthalmologist may prescribe only eyeglasses until your cataract becomes more mature.

There are different types of IOLs on the market, each having specific advantages and disadvantages, which your eye doctor can inform you about. The Aspira-aXA IOLs (aXA IOL platform) have a larger optic diameter than most other IOLs on the market. This means that patients who are implanted with these devices may perceive less often or less strongly disturbing optical phenomena which can occur after IOL implantation. Moreover, the IOLs offer your ophthalmologist a wide view of your retina after surgery, which can help find and watch potential other ocular diseases.

7. Suggested training for users

Not applicable to patients since the device is to be used only by healthcare professionals.

⁴ **Vitreous:** A clear gel that fills the space between the lens (in the front of the eye) and the retina

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Restrisiken und unerwünschte Wirkungen	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Vorsicht	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to

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Common specification(s) to comply with, if applicable:
N/A
Harmonised standards to apply, if applicable:
DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
DIN EN ISO 14971:2022-04 Medical devices – Application of risk management to medical devices
Relevant adopted monographs of the European Pharmacopoeia
Not applicable

9. Referenzen

- [1] Wendelstein J, Laubichler P, Fischinger I, Mariacher S, Beka S, Mursch-Edlmayr A, et al. Rotational Stability, Tilt and Decentration of a New IOL with a 7.0 mm Optic. *Curr Eye Res.* 2021;46(11):1673-1680.
- [2] Holladay JT, Piers PA, Koranyi G, van der Mooren M, Norrby NE. A new intraocular lens design to reduce spherical aberration of pseudophakic eyes. *J Refract Surg.* 2002;18(6):683-691.
- [3] Sieburth R, Chen M. Intraocular lens correction of presbyopia. *Taiwan J Ophthalmol.* 2019;9(1):4-17.
- [4] Kimura S, Morizane Y, Shiode Y, Hirano M, Doi S, Toshima S, et al. Assessment of tilt and decentration of crystalline lens and intraocular lens relative to the corneal topographic axis using anterior segment optical coherence tomography. *PLoS One.* 2017;12(9):e0184066.
- [5] Pilger D, Bertelmann E, Brockmann T, Duncker T, Duncker GIW, Schrecker J. Postoperative Lens Rotation of a 7.0 mm Optic IOL with Plate Haptics. *J Ophthalmic Vis Res.* 2021;5(1):014-020.
- [6] Schrecker J, Seitz B, Langenbucher A. Performance einer neuen 7-mm-Intraokularlinse mit Nachbeobachtung über 1,5 Jahre [Performance of a new 7 mm intraocular lens with follow-up over 1.5 years]. *Ophthalmologe.* 2022;119(4):367-373.
- [7] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. *J Cataract Refract Surg.* 2022;48(1):75-82.
- [8] Rozendal LRW, Van Vught L, Beenakker JWM, Luyten GPM. Bag-to-bag intraocular lens exchange with a wide optic intraocular lens for treatment of positive and negative dysphotopsia. 39th Congress of the ESCRS 2021.
- [9] Borkenstein AF, Borkenstein EM. Creating Hybrid Monovision with 7.0 mm XL Optic and High-Add AMD Intraocular Lenses (XL-MAGS) in a Patient with Retinitis Pigmentosa. *Case Rep Ophthalmol.* 2019;10(3):304-311.
- [10] Borkenstein AF, Borkenstein EM. Efficacy of Large Optic Intraocular Lenses in Myopic Eyes with Posterior Segment Pathology. *Ophthalmol Ther.* 2022;11(1):443-452.
- [11] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia
- [12] Küchle M. Comparison of visual function with aspheric yellow, aspheric clear and spherical clear intraocular lenses. *J Emmetropia* 2013;4:123-130.

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Остатъчни рискове и нежелани ефекти	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Внимание	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to

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Zbytková rizika a nežádoucí účinky	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Upozornění	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Υπολειπόμενοι κίνδυνοι και ανεπιθύμητες ενέργειες	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Προσοχή	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Riesgos residuales y efectos indeseables	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Precaución	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Risques résiduels et effets indésirables	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric <u>IOLs</u>: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Mise en garde	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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4. Rizici i upozorenja

Preostali rizici i neželjeni učinci	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Oprez	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss)

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Maradék kockázatok és nemkívánatos hatások	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to <u>note</u> in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Figyelem	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Rischi residui ed effetti indesiderati	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Attenzione	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Liekančios rizikos ir nepageidaujami reiškiniai	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection <u>with</u> toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Atsargiai	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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4. Risico's en waarschuwingen

Resterende risico's en ongewenste effecten	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Waarschuwing	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to

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Ryzyko resztkowe i działania niepożądane	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Ostrzeżenie	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to

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4. Riscos e avisos

Riscos residuais e efeitos indesejáveis	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A_postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Precaução	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to

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Riscuri reziduale și efecte nedorite	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Atenție	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to

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Kvarstående risker och biverkningar	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Varning	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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