Summary of safety and clinical performance (SSCP) for the Monofocal Aspira-aA Family

SSCP reference number: P034

Revision history

Rev number	Revision date	Description of change	Revision validated by the Notified Body
Rev.01	2025-03-11	Initial creation	☒ Yes Validation language: English ☐ No (only applicable for class IIa or some IIb implantable devices (MDR Article 52(4) 2nd paragraph) for which the SSCP is not yet validated by the NB)
Rev.02	2025-02-18	Addition of new models Aspira-aA+ and Aspira-aAY+	☒ Yes Validation language: English ☐ No

Table of contents

	Page
Part 1 – SSCP intended for healthcare professionals	3
1. Device identification and general information	3
2. Intended use of the device	5
3. Device description	5
4. Risks and warnings	8
5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)	10
5.1. Summary of clinical data related to equivalent device	10
5.2. Summary of clinical data from conducted investigations of the device before the CE-marking	10
5.3. Summary of clinical data from other sources	11
5.4. Overall summary of the clinical performance and safety	10
5.5. Ongoing or planned post-market clinical follow-up	16
6. Possible diagnostic or therapeutic alternatives	17
7. Suggested profile and training for users	17
8. Reference to any harmonised standards and CS applied	17
9. References	18
Part 2 – SSCP intended for patients	20
1. Device identification and general information	20
2. Intended use of the device	21
3. Device description	21
4. Risks and warnings	23
5. Summary of clinical evaluation and post-market clinical follow-up	25
5.1. Clinical background of the device	25
5.2. The clinical evidence for the CE-marking	25
5.3. Safety	25
6. Possible diagnostic or therapeutic alternatives	26
7. Suggested training for users	26

Part 1 – SSCP intended for healthcare professionals

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The SSCP is not intended to replace the Instructions For Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients. This SSCP is reviewed at least annually and will be updated if needed to ensure that any clinical and/or safety information in the SSCP remains correct and complete.

The following information is intended for users/healthcare professionals.

Following this information there is a summary intended for patients.

1. Device identification and general information

Manufacturer's name and address	HumanOptics Holding AG Spardorfer Straße 150 91054 Erlangen Germany
Manufacturer's single registration number (SRN)	DE-MF-000017892
Authorised representative if applicable; name and SRN	N/A
Notified body name and organization number	TÜV NORD CERT GmbH NB 0044

Device trade name(s)	• Aspira-aA • Aspira-aAY • Aspira-aA+ • Aspira-aAY+ • Torica-aA • Torica-aAY (collectively referred to as "Monofocal Aspira-aA Family")
Basic UDI-DI(s)	• Aspira-aA/-aAY and Aspira-aA+/-aAY+: 04049154_PC_M2_H2_O2_BU • Torica-aA/-aAY: 04049154_PC_M2_H2_O3_BX
CND/EMDN code(s):	• Aspira-aA/-aAY and Aspira-aA+/-aAY+: P030102090202 (IOLs, APHAKIC, MONOFOCAL, ASPHERICAL, HYDROPHILIC ACRYLIC) • Torica-aA/-aAY: P030102090302 (IOLs, APHAKIC, MONOFOCAL, TORIC, HYDROPHILIC ACRYLIC)
Class of device	☐ I ☐ IIa ☒ IIb ☐ III

Year when the first certificate (CE) was issued covering the device

- Aspira-aA/-aAY: 2007
- Torica-aA/-aAY: 2012
- Aspira-aA+/-aAY+: Not applicable at this stage, as the product has been submitted for CE-marking.


2. Intended use of the device

Intended purpose	The devices are intended for implantation in the capsular bag to replace the natural lens.
Intended users	The devices must be handled by health professionals and implanted by ophthalmic surgeons.
Use environment	Ophthalmic operating room
Indications	The devices are indicated for the correction of aphakia after surgical removal of the natural lens. In addition, toric IOLs (Torica-aA/-aAY) are indicated to correct pre-existing regular, stable corneal astigmatism.
Medical condition(s) to be treated	Aphakia after surgical removal of the natural lens
Intended patient population	Aphakic adult patients Note: children are excluded due to the absence of clinical data and additional risks associated with paediatric cataract surgery. The manufacturer does not hold clinical data relating to pregnant/breastfeeding women or immunocompromised patient populations.
Intended anatomical location	Posterior chamber of the eye
Contraindications	Outside of general contraindications of ocular surgery, there are no specific contraindications for this device.

3. Device description

Short description of the devices	One-piece, monofocal, aspheric, posterior-chamber intraocular lens (IOL) with C-loop haptics and 360° lens epithelial cell barrier on the posterior optic surface - Aspira-aA/-aAY: with aberration-free¹ optic - Torica-aA/-aAY: with toric, aberration-free¹ optic - Aspira-aA+/-aAY+: with adjusted aspheric optic
Description of the devices	Aspira-aA/-aAY, Aspira-aA+/-aAY+ and Torica-aA/-aAY are foldable, single-piece posterior chamber IOLs with non-angulated C-loop haptics and a monofocal, aspheric optic. The toric Torica-aA/-aAY is additionally fitted with a toric anterior optic surface (astigmatism correction) and has additionally markings in the periphery on the anterior surface, indicating the axis of the plus cylinder. The model Aspira-aA+/-aAY+ has an adjusted aspheric optic that distributes a portion of the light from the far focus toward the intermediate region to increase the subjective depth of field compared to the Aspira-aA/-aAY. The IOLs have an overall diameter of 12.5 mm and an optical diameter of 6.0 mm. The IOLs are made of hydrophilic acrylic copolymer with UV-absorber, a water content of 26% at 35°C and a refractive index of 1.46. The yellow variants (with "Y") additionally contain a blue light filter that absorbs the high-energy portion of the blue light between 400 and 500 nm.

¹ i.e. neutral with respect to spherical aberration

	The posterior surface of the optic features a 360° lens epithelial cell (LEC) barrier that is aimed to counteract the migration of the LEC, thereby reducing the frequency of posterior capsule opacification (PCO). For the Aspira-aA/-aAY and the Aspira-aA+/-aAY+, the dioptric range covers -20.0 D to +60.0 D and 0.0 D to +30.0 D, respectively. For the Torica-aA/-aAY, the dioptric range covers spherical equivalents between -20.0 D to +60.0 D and cylinders between +1.0 D and +30.0 D (max. total power per meridian: ≥ -20.5 D and $\leq +60.5$ D).  <i>Schematic lateral and frontal view of the non-toric (left) and toric models (right). The dimensions are given in mm.</i>
Mode of action / operating principle	The devices are designed to focus light passing through the cornea and pupil onto the retina, like the eye's natural lens. Monofocal IOLs provide clear vision at one single focal point (usually far focus). The aberration-free aspheric optic design of the models Aspira-aA and Aspira-aAY aims to provide high-contrast imaging compared to a spherical design. The models with ending aA+ and aAY+ have an adjusted aspheric optic that distributes a portion of the light towards the intermediate region to increase the subjective depth of field. Toric IOLs correct corneal astigmatism. Their function can be described as follows: corneal astigmatism means that the anterior corneal surface and/or the posterior corneal surface are not rotationally symmetrical. Therefore, such a cornea has a meridian of strongest refractive power and one of weakest refractive power, both of which are usually orientated perpendicular to each other (regular astigmatism). As a result of its toric anterior optic surface, the Torica-aA/-aAY also has meridians of strongest and weakest refractive powers. Those are arranged opposite to the corneal effect in order to correct the astigmatism.
Materials or substances in contact with the patient's tissues	The hydrophilic acrylic copolymer consisting of 2-hydroxyethyl methacrylate (HEMA) and methyl methacrylate (MMA) with a UV absorber and with a yellow coloring agent (only for models followed by "Y") is permanently in contact with the ocular tissues.
Method of sterilisation and allowable number of uses	The devices are sterile (steam sterilised) and designed for single use. Reprocessing of an implant is strictly prohibited.
Previous generations of the device(s)	The precursor model of the Aspira-aA/-aAY is the AS IOL, a lens with an optical diameter of 5.75 mm, a total diameter of 12.0 mm and with a spherical anterior and posterior surfaces. The Aspira-aA+/-aAY+ is based on the parent model Aspira-aA/-aAY and equivalent to the latter according to Annex XIV, Section 3 of Regulation (EU) 2017/745 (MDR) and MDCG 2020-5.

	The Torica-aA is the result of a combination of the Aspira-aA platform and the toric optic of the former Torica-s, a 3-piece posterior-chamber IOL made of a silicone optic and PMMA haptics.
Variants and/or configurations covered by this document	• Colour: • Colourless: device names ending with “-aA” • Blue-filtering: device names ending with “-aAY” • Packaging types: • Compact line (CL) • Safeloader (SL) Note: The Safeloader system is a packaging type and not an accessory.
Accessories which are intended to be used in combination with the device	Not applicable Note: Injector systems are independent medical devices which are intended to be used in combination with the device in question (i.e. these devices are not accessories).
Other devices and products which are intended to be used in combination with the device	Implantation of the device can be performed using either an injector system (for Safeloader and Compact line) or forceps (for Compact line only). A list of suitable injection systems can be found under www.humanoptics.com. When using a non-tested injector system for implantation of HumanOptics IOLs packaged in the Compact line, the user must refer to the specific instructions for use provided with the injector system to ensure that it is suitable to use with the IOL to be implanted.

4. Risks and warnings

Residual risks and undesirable effects	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Caution	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

Warnings	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Special considerations prior to toric IOL implantation	• The dioptric power of toric IOLs is indicated as “spherical equivalent (SE) and cylinder (cyl)”. Please check the label carefully. • Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism. • Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation • If possible, posterior corneal astigmatism (tomography) should be considered, especially with astigmatism < 2.0 D.
Other relevant aspects of safety including a summary of any field safety corrective action (FSCA including FSN), if applicable	There have been no recalls since the devices have been launched. In July 2024, HumanOptics issued a Field Safety Notice affecting all devices in a Safeloader packaging manufactured in the period January 1, 2020 to June 28, 2024. During the standardized incoming goods inspection of returns, HumanOptics had observed very rare cases of damaged outer blister packaging after transportation. The FSN informed the users about the observation and reminded them to check the integrity of the sterile barrier system, as requested in the Instructions for Use, along with a visual instruction how to identify potential cracks. This damage is visible in the form of cracking and can be identified by users before use. It could theoretically cause a contamination of the surgical equipment or the intraocular lens. However, the hazard potential due to crack is considered to be very low. Statistical valuation showed that probability of occurrence of the additional risk of endophthalmitis due to the cracks is by several orders of magnitude lower than the general rate of endophthalmitis after cataract surgery. No patient injury has happened to the manufacturer’s best knowledge. Meanwhile, the product packaging has been improved to eliminate this issue completely.

5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)

5.1. Summary of clinical data related to equivalent device

The Aspira-aA+/aAY+ is based on the parent model Aspira-aA/-aAY and equivalent to the latter according to Annex XIV, Section 3 of the Regulation (EU) 2017/745 (MDR) and MDCG 2020-5. Clinical data for the Aspira-aA/-aAY are presented in Chapter 5.3.

5.2. Summary of clinical data from conducted investigations of the device before the CE-marking

Not applicable. All clinical investigations were performed after the first CE marking under the European Medical Device Directive 93/42/EEC (see section 5.3).

5.3. Summary of clinical data from other sources

Clinical data confirming the safety and performance of the Monofocal Aspira-aA Family have been derived from the following sources:

- Several post-market studies
- A comprehensive literature review using PubMed and other databases

Literature search strategies were designed to identify articles relevant to the devices. Both favorable and unfavorable references were identified and summarized.

5.3.1. Post-market studies

The manufacturer HumanOptics has initiated several post-market studies with the Aspira-aA(Y) and Torica-aA/-aAY as shown in Table 1. The studies/activities were labeled with acronyms considering: [device in focus]_[city]_[primary end point].

In 7 post-market studies performed by the manufacturer 533 eyes belonging to 409 patients were implanted with the Aspira-aA or the Aspira-aAY. The results collected over follow-up periods of up to 3 years demonstrate a very low rate of (mostly mild) complications and a good visual performance as well as a high level of patient satisfaction. Moreover, the data indicate that the Aspira-aA and the Aspira-aAY exhibit a comparable visual performance and safety, confirming that the safety of the blue light filter. Another clinical study demonstrated that the Aspira-aA can be safely used in eyes with pseudoexfoliation syndrome (PEX), which is known to entail an increased risk of complications. With regards to the manufacturing finishing process, the Aspira-aA showed a similarly good visual performance in both the polished and non-polished version. However, the non-polished IOLs with the sharper edges slightly outperformed the polished IOLs in terms of LEC outgrowth, capsular opacification and contrast sensitivity. Due to the overall low rates of clinically relevant PCO and the absence of other complications in both groups, those results nevertheless support a high level of safety and a good performance for both versions. For the Aspira-aA+ and Aspira-aAY+ currently available clinical data are based on product equivalence with the Aspira-aA and Aspira-aAY, only. After CE-marking, product-specific clinical data will be systematically collected.

In terms of IOL platform stability, the results of 3 completed post-market studies involving 240 eyes of 143 patients, one conducted with a Aspira-aA version with toric markings and two with the Torica-aA, reveal a good rotational stability of the lens at the 1-month interim analysis and the 6-month follow-up. Additionally, one of the Torica-aA studies validates the good long-term safety and performance of the lens, including IOL rotational stability, visual acuity, astigmatism correction, and patient satisfaction over a period of 1.5 years.

Details of the study results are provided in the subsequent paragraphs.

Table 1: Overview of clinical studies/activities with the Monofocal Aspira-aA Family

Name of activity / Type of activity / Design	Aim	Study status / Postoperative follow-up	Indication / Number of <u>eyes / patients</u> (<u>subjects</u>) implanted with the <u>device under</u> <u>consideration</u>	Patient population (vulnerable groups)
Aspira-aA/-aAY				
Asp-aA_Erl_VA Post-market study RCT	Comparative study to confirm the long-term (1-year) performance and safety of the Aspira-aA and Aspira-aAY vs. AS after implantation in cataract patients. To assess the patient-related outcome based on patient satisfaction.	Completed 1 year (completed: 2011)	Cataract 60 / 60	Age: 53 to 90 years Pregnant or lactating subjects: N/A Children: N/A
Asp-aA_Wie_Rota Post-market study PSAS	Single-arm study to evaluate the mid-term (6-month) rotational stability of the Aspira-aA (former name: MC 6125 AS) after implantation in cataract patients. Besides IOL orientation, no other data were collected.	Completed 6 months (completed: 2011)	Cataract 150 / 84	Age: 50 to 90 years Pregnant or lactating subjects: N/A Children: N/A
Asp-aA_Hal_Rota Post-market study PNCS	Comparative study to evaluate the mid-term (8-10 weeks) in-the-bag stability and visual performance of two monofocal single-piece IOLs with different haptic designs (Aspira-aA and Aspira-aQA).	Completed 8-10 weeks (completed: 2014)	Cataract 20 / 17	Age: 40 to 82 years Pregnant or lactating subjects: N/A Children: N/A
Inv-aA_Gla_SA Post-market study RCT	Comparative study to evaluate the visual performance and safety of a new custom-designed aspheric IOL (Invidua-aA) vs. the standard Aspira-aA.	Completed 3 months (completed: 2017)	Cataract 29 / 27	Age: 62 to 80 years Pregnant or lactating subjects: N/A Children: N/A
Asp-aA_Gla_Poli Post-market study RCT	Comparative study to assess the long-term (2-year) impact of tumble-polishing on the performance and safety of the Aspira-aA.	Completed 2 years (completed: 2018)	Cataract 29 / 24 (tumble-polished device) 31/25 (polish-free device)	Age: 52 to 90 years Pregnant or lactating subjects: N/A Children: N/A
Asp-aA_Lub_PEX Post-market study PNCS	Comparative study to evaluate the capsular and uveal biocompatibility (up to 1 year) of the Aspira-aA in a group of eyes with pseudoexfoliation syndrome (PEX) and to compare the results with those of a group of healthy eyes.	Completed 1 year (completed: 2019)	Cataract 75 / 65 (PEX) 74 / 64 (healthy)	Age: 44 to 90 years Pregnant or lactating subjects: N/A Children: N/A

Asp-aXA_Ora_ND Post-market study RCT	Comparative study to confirm the long-term (3-year) performance and safety of the Aspira-aXA vs. the Aspira-aA after implantation in cataract patients. To assess the patient-related outcome based on patient satisfaction (including photopic disturbances).	Completed 3 years (completed: 2021)	Cataract 65 / 43	Age: 55 to 79 years Pregnant or lactating subjects: N/A Children: N/A
Torica-aA/-aAY				
Tor-aA_Lin_Rota Post-market study PSAS	To confirm the long-term (1.5-year) performance (incl. astigmatism correction and visual outcomes) and safety (including rotational stability) of the Torica-aA after implantation in cataract patients with pre-existing astigmatism. To assess the patient-related outcome based on patient satisfaction.	Completed 1.5 years (completed: 2015)	Cataract 40 / 26	Age: 54 to 85 years Pregnant or lactating subjects: N/A Children: N/A
Tor-aA_Ros_Rota Post-market study PSAS	To compare the rotational stability of the Torica-aA between eyes with axial length (AL) below and above 25 mm. Other parameters related to visual performance and safety will be evaluated as well. To assess the patient-related outcome based on patient satisfaction.	Completed 4-6 months (completed: 2023)	Cataract 25 / 16 (AL <25.0 mm) 25 / 17 (AL ≥25.0 mm)	Age: 50 to 82 years Pregnant or lactating subjects: N/A Children: N/A

FU=Follow-up; RCT=randomized controlled trial; PSAS=prospective single-arm study; PNCS=prospective non-randomized comparative study; M=month

Post-market study Asp-aA_Erl_VA (Aspira-aA vs. Aspira-aAY vs. AS – Erlangen, Germany) [1]

In this randomized controlled clinical study, the visual performance of the Aspira-aA and Aspira-aAY (and also the AS) was evaluated up to 1 year postoperatively. These two product variants are identical to each other in every way except for the blue-light filter contained in the Aspira-aAY. Both groups performed well and without significant differences in terms of monocular distance visual acuity and contrast sensitivity. The data are also in line with those reported for other modern aspheric hydrophilic and hydrophobic monofocal IOLs. As expected for monofocal IOLs at plano target refraction the intermediate and near visual acuity was lower. Notably, the blue-light filter also did not appear to affect color perception under mesopic and photopic conditions. The level of patient satisfaction was high in both groups (all patients reported to be “satisfied” or “very satisfied”).

Out of 30 patients/eyes operated per group, 28 and 19 could be followed over the whole study period for Aspira-aAY and Aspira-aA, respectively. Except for the intraoperative exchange of an IOL due to a broken haptic, there were no intraoperative complications. In each group one patient required Nd:YAG laser capsulotomy for PCO treatment. No other postoperative complications were reported. Taken together, the data demonstrate a good visual performance and safety of both Aspira-aA and Aspira-aAY, with no group-specific differences. These findings are in line with other clinical studies, which, taken together, could not conclusively demonstrate a clinically relevant difference between blue-filtering and colourless IOLs.

In conclusion, the presence of the yellow dye does not affect the clinical performance and safety of the medical devices under consideration. Clinical results obtained with one model can be regarded as valid for both variants (colourless or yellow).

Post-market study Asp-aA_Wie_Rota (Aspira-aA – Wien, Austria)

The purpose of this prospective single-arm study was to evaluate the rotational stability and anterior chamber depth of the Aspira-aA as IOL platform after implantation in the capsular bag (n=150 eyes). A strength of this clinical study was the introduction of an objective photograph-assisted method to measure IOL rotation. The follow-up covered a time span of 6 months.

The mean IOL rotation was $3.3^{\circ} \pm 4.0^{\circ}$ at 6-months postoperatively compared to the 1-hour baseline state. IOL rotation was $\leq 5^{\circ}$ in 75.9%, $< 10^{\circ}$ in 94.8% and $< 20^{\circ}$ in 98.3% (with one eye exceeding a rotation of 20°). The study also revealed that IOL rotation was independent of axial length or initial axis alignment.

In conclusion, the Aspira-aA performed very well in terms of rotational stability. As the obtained results overachieve the acceptance criteria of the DIN EN ISO 11979-7:2018-08, the Aspira-aA platform (C-loop haptics) was determined to be suitable for the development of a toric variant. Since the study was focused exclusively on rotational stability, no data were reported concerning complications and visual performance.

Post-market study Asp-aA_Hal_Rota (Aspira-aA vs. Aspira-aQA – Halle, Germany)

This prospective non-randomized comparative study was primarily designed to assess the capsular stability of the Aspira-aA and another monofocal IOL (Aspira-aQA; not on the market any longer) in order to determine, whether they were suitable as platforms for the development of a customized IOL. Although this primary endpoint could not be addressed (due to insufficient mydriasis and poor image quality), the remaining data are useful for assessing the safety and performance of the Aspira-aA. 18 out of 20 eyes implanted with the device were available at the last follow-up 8-10 weeks postoperatively.

The results confirm that the Aspira-aA performed well in terms of visual acuity with a mean monocular UDVA and CDVA (logMAR) of 0.20 ± 0.14 and 0.10 ± 0.11 , respectively (8-10 weeks postoperatively).

The eyes implanted with Aspira-aA displayed a relatively high but still acceptable absolute deviation from target SE of 0.61 ± 0.43 and a mild hyperopic shift. The latter finding emphasizes the need for each surgeon to optimize the manufacturer's A-constants for themselves. Despite the above-mentioned refraction errors, the monocular UDVA after 8-10 weeks was good overall (mean 0.20 ± 0.14 logMAR). Moreover, the eye with the worst UDVA (1.0 logMAR at visit 1 and 0.4 at visit 2) had an extreme axial length of 31.3 mm, which was probably causative for the large deviation from target refraction (1.51 D).

The vast majority of patients (77.8%) were satisfied or very satisfied with their vision. For three eyes however only "slight" satisfaction was expressed. This can be explained by complications or ocular diseases (macular edema and posterior vitreous detachments), which occurred in these cases, and is very likely unrelated to the IOL.

In terms of complications mild PCO was observed in two eyes implanted with the Aspira-aA (11%). However, it has to be stressed out that these findings were not clinically relevant as none of the eyes required Nd:YAG laser capsulotomy. Postoperative not IOL-related ocular complications (macular edema and posterior vitreous detachments) which were noticed in 3 out of 20 eyes, did not affect the visual performance.

In conclusion, the results support a good visual performance of Aspira-aA, although the number of eyes (20) is limited. The number of adverse events observed was also low but the small number of cases and the short follow-up period (up to 10 weeks) do not allow a definite conclusion.

Post-market study Inv-aA_Gla_SA (Aspira-aA vs. Invidua-aA – Glauchau, Germany) [2]

The goal of this randomized controlled trial was to compare the visual performance of a custom-made Aspira-aA variant designed with individual aberration correction (Invidua-aA) with that of the standard aspheric aberration-free Aspira-aA of otherwise identical design. Since analysis of complete study is still ongoing, interim results for a follow-up period of 3 months (28 eyes) are discussed below.

The Aspira-aA achieved a very good monocular UDVA and CDVA 3 months postoperatively with logMAR values of 0.03 ± 0.09 and -0.01 ± 0.03 , respectively. There was no significant difference between the 2 groups. As expected for monofocal IOLs with plano target refraction, the intermediate and near visual acuities were lower compared to distance visual acuity. Refractive accuracy was overall high with both lenses. However, the Invidua-aA outperformed the Aspira-aA in terms of contrast sensitivity under photopic conditions and for some spatial frequencies under mesopic conditions. This reflects the successful individual compensation for corneal spherical aberration in eyes treated with the Invidua-aA. Nevertheless, the aberration-free Aspira-aA also delivered acceptable results in terms of contrast sensitivity, which confirms the observations made in the post-market studies Asp-aA_Erl_VA, Asp-aA_Gla_Poli and Asp-

aXA_Ora_ND. The monocular defocus curves were similar and not significantly different in both groups, with a plateau ranging from +0.50 D to -0.50 D of defocusing, and a median visual acuity of 0.00 logMAR.

Three eyes (5%) in the Invidua-aA group and one eye (4%) in the Aspira-aA developed PCO and required Nd:YAG laser capsulotomy 3 months postoperatively. No other adverse events were reported.

In conclusion, the Aspira-aA performed well in terms of visual acuity and contrast sensitivity and exhibited a very low complication rate over the currently available follow-up period of 3 months.

Post-market study Asp-aA_Gla_Poli (Aspira-aA – Glauchau, Germany)

The purpose of this randomized comparative study was to evaluate whether the difference in the finishing process (tumble-polished vs. non-polished) of an otherwise identical hydrophilic acrylic IOL (Aspira-aA) could influence the clinical performance and the optical quality of the lenses and consequently quality of vision. A strength of this clinical study is the relatively long follow-up period (2 years) and the low number of lost-to-follow-up subjects (26 out of 29 eyes implanted with the polished and 29 out of 31 eyes implanted with the non-polished version were available until the last follow-up, respectively).

The data demonstrate that the Aspira-aA performed very well in both groups in terms of visual acuity 2 years postoperatively with mean monocular UDVA and CDVA of 0.09 ± 0.12 logMAR and 0.00 ± 0.08 logMAR for polished IOLs and 0.06 ± 0.11 logMAR and -0.01 ± 0.03 logMAR for non-polished IOLs, respectively. Intermediate visual acuity (UIVA/CIVA) was slightly lower, as expected for a monofocal IOL at plano target refraction. Overall, the results show stable distance visual acuity over time with no significant differences between the polished and non-polished IOLs. Moreover, the absolute deviation from target refraction was low in both groups (median 0.24 D (0.01 to 1.03) for polished and 0.25 D (0.01 to 1.74) for non-polished lenses, respectively).

LEC outgrowth at the capsulorhexis edge was almost non-existent in this study despite other published findings on hydrophilic lens materials. Even after 2 years, there were only occasional occurrences of LEC leading to a mean score of 0.50 and 0.30 in the polished and non-polished group, respectively. There was however a slight tendency to less LEC outgrowth in the non-polished IOLs (IOLs with a sharper edge) than in the polished IOLs. Similarly, the non-polished IOL group exhibited significantly overall lower PCO scores at both follow-ups. These findings are in line with the available literature reporting that a sharp posterior IOL edge prevent the development of PCO. However, the incidence of clinically relevant PCO (requiring Nd:YAG laser treatment) was low in both groups (two eyes in the polished and one eye in the non-polished group).

In agreement with the PCO scores, monocular examination of contrast sensitivity showed a trend towards better contrast sensitivity in the non-polished IOLs compared to the polished IOLs under photopic and mesopic conditions. However, in both groups, values were mainly within or above the age-related normal ranges measured with the same device. In terms of scattered light, there were no statistically significant differences between the two groups, i.e. the omission of the polishing step did not result in an increase of stray light. Finally, mean ocular spherical and coma aberrations were similar in both groups with a mean internal aberration close to zero, confirming the aberration-free design of the Aspira-aA.

In terms of complications, no ocular adverse events or other complications were reported except the three PCO cases mentioned above.

In conclusion, the Aspira-aA showed a similarly good visual performance in both the polished and non-polished version over a period of two years. However, the non-polished IOLs with the sharper edges slightly outperformed the polished IOLs in terms of LEC outgrowth, capsular

opacification and contrast sensitivity. Due to the overall low rates of clinically relevant PCO and the absence of other complications in both groups, these clinical study results nevertheless support a high level of safety and a good performance for both versions.

Post-market study Asp-aA_Lub_PEX (Aspira-aA – Lublin, Poland) [3]

The main purpose of this comparative prospective clinical study was to evaluate the capsular and uveal biocompatibility of the Aspira-aA in a group of PEX (pseudoexfoliation syndrome) eyes up to 1-year post-operatively and to compare the results with those of a group of healthy eyes. Seventy-five and 74 eyes per group were implanted with the device, respectively.

The Aspira-aA performed very well in both groups in terms of monocular visual acuity with mean UDVA and CDVA 1-year postoperatively of 0.30 ± 0.24 logMAR and 0.02 ± 0.18 logMAR for the PEX eyes and 0.16 ± 0.18 logMAR and -0.03 ± 0.16 logMAR for the healthy eyes, respectively. As expected, the outcomes were significantly better in the healthy group compared to the PEX group due to the ocular pathology of the PEX eyes. The CDVA in the PEX group was also similar to what has been reported in the literature with the same IOL type. At the 1-year visit median absolute deviation from target refraction was around 0.5 D in both groups, which is reflected in general by better CDVA vs UDVA values.

As PEX eyes have a stronger immune reaction to IOLs, they have been proposed as a test model for IOL biocompatibility. In this study, flare values (inflammatory parameter indicating the amount of proteins in the aqueous humour) were overall very low in both groups, indicating a high biocompatibility of the Aspira-aA. Flare values were higher in PEX eyes than in healthy eyes at 1-week and 3-months postoperative visit but smaller at 1-year. The decrease of flare values for PEX eyes at 1-year might be caused by the anti-glaucomatous medication, which was received by a higher portion of eyes in the PEX group compared to the healthy group (see discussion in the final clinical study report).

Both ACO and PCO scores increased over time. Nonetheless, the PCO scores were very low with all mean score values below 1.0. Most eyes in both groups showed no LEC outgrowth at the capsulorhexis edge. PCO scores tended to be higher in healthy eyes, whereas ACO scores were higher in the PEX group at some follow-up visits. It can be doubted, if these differences are clinically meaningful.

Cataract surgery in PEX eyes is usually challenging because of poor pupillary dilation and progressive zonular weakness and is known to be associated with a greater number of complications (in particular postoperative inflammation, capsule opacification, and IOL luxation). Nevertheless, no intraoperative complications occurred during the study and all IOLs remained well centred (subjective assessment). Postoperatively, there were no ocular adverse events except for some expected cases of PCO and ACO which required Nd:YAG laser capsulotomy. One year after surgery, 2 eyes in the PEX group and 1 eye in the healthy group received this treatment due to ACO and PCO, respectively.

In conclusion, the Aspira-aA performed similarly in both, the PEX and healthy groups. The results showed no concerns in terms of uveal and capsular biocompatibility; therefore, the IOL can be considered safe even in eyes with ocular comorbidities which are at higher risk to develop complications. Moreover, it is worth mentioning in this context that the hydrophilic acrylic material used for manufacturing the present lenses was ActiCryl H26, whereas in the other clinical studies the ContaFlex materials (CI26 and CI26Y) were used.

Post-market study Asp-aXA_Ora_ND (Aspira-aXA/Aspira-aA – Oranienburg, Germany)

In the course of this randomized controlled trial, the visual performance of the Aspira-aXA, a new IOL with a 7.0 mm optic (n=66) was compared with that of the well-proven Aspira-aA (n=65) with features a 6.0 mm optic. The follow-up covered the exceptionally long time span of 3 years.

The results show that both IOLs performed very well in terms of monocular visual acuity and that there were no significant differences between the IOL groups at any postoperative time point and the values remained fairly constant over time. The data are also in line with those reported in previous studies showing similar visual acuity outcomes with the Aspira-aA and its toric version as well as other modern aspheric hydrophilic and hydrophobic monofocal IOLs. The absolute deviation from target refraction was slightly smaller in the Aspira-aXA group than in the Aspira-aA group ($p=0.006$) with median values of 0.36 D (0.01 to 1.60) and 0.60 D (0.03 to 1.81), respectively.

In this study, monocular examination of contrast sensitivity showed that the Aspira-aA performed overall slightly better than the Aspira-aXA. Nonetheless, and importantly, the scores of both IOLs were on average above age-related normal ranges measured with the same device.

Regarding photic phenomena, there was no significant difference between the two groups in terms of occurrence of halos and glare sensitivity at the 1-year and 3-year follow-ups. Negative dysphotopsia was completely absent in the Aspira-aXA group at the 1-year visit, although the difference compared to the Aspira-aA group (4 eyes, i.e. 6.5% at least sometimes experienced the phenomena) did not reach significance ($p=0.131$). After 3 years, the occurrence of negative dysphotopsia decreased in the Aspira-aA group (only 1 patient “always” affected) while it was reported to occur “sometimes” by 2 patients in the Aspira-aXA group. Since these 2 patients developed slight capsule opacification over time and since negative dysphotopsia is acknowledged to occur just after cataract surgery, it is likely that the “unwanted image” perceived by these patients was a consequence of the development of a secondary cataract. At the 3-year visit, the Aspira-aXA also tended to outperform the Aspira-aA in terms of glare sensitivity in the night (4 eyes/7.8% vs. 11 eyes/19.3%; $p=0.178$).

In terms of patients’ satisfaction with their overall vision, the Aspira-aXA performed significantly better ($p=0.024$) than the Aspira-aA at the 1-year follow-up, with 90.5% of the patients being “very satisfied” (73.8%) or “satisfied” (16.7%) as opposed to 55.0%, and 42.5% in the Aspira-aA group, respectively. After 3 years of follow-up, the levels of satisfaction increased in both groups and the difference didn’t reach significance anymore ($p=0.250$). Nonetheless, the outcomes tended to be better with the Aspira-aXA with most patients being “very satisfied” (74.3%) or “satisfied” (20.0%) compared to 57.5% or 37.5% with the Aspira-aA, respectively.

No IOL-related adverse events occurred during the course of the study and all IOLs remained well centred except two Aspira-aA where the optic was in front of the rhexis at the 1-year follow up. These two patients did not receive any additional treatment because they were “very satisfied” and “satisfied” with their overall vision. Capsulotomy with Nd:YAG laser was not necessary for any patient within the first postoperative year. After 3 years, only one eye (1.5%) in the Aspira-aXA group and two eyes (3.1%) in the Aspira-aA group required this treatment due to the development of PCO and associated reduction of visual acuity.

In conclusion, both IOLs performed well in terms of visual acuity, contrast sensitivity and patients’ satisfaction over a long-term follow-up of 3 years and showed to be a safe and effective surgical treatment for cataracts. Overall, there were no clinically relevant differences in the visual function of patients implanted either with the Aspira-aXA or Aspira-aA IOLs, with however a trend for a higher level of patient satisfaction in the Aspira-aXA group.

Post-market study Tor-aA_Lin_Rota (Torica-aA – Linz, Austria) [4, 5]

The purpose of this prospective single-arm study was to evaluate rotational stability, astigmatism correction, visual outcomes, and patient satisfaction after implantation of the Torica-aA IOL following cataract extraction. A strength of this bipartite clinical study is the relatively long follow-up period (1.5 years) and the fact that all 26 patients (40 eyes) attended all planned visits.

The results show that the IOL performs very well in terms of monocular UDVA/UNVA (depending on target refraction) and CDVA 6 months and 1.5 years postoperatively. The eyes targeted for emmetropia (n=36) achieved a mean UDVA of 0.09 ± 0.12 logMAR with 94% of eyes having ≤ 0.3 and 74% of eyes having ≤ 0.1 . The three eyes with myopic target refraction achieved a UNVA of ≤ 0.0 logMAR. Regarding CDVA, 97% of eyes (n=38) achieved ≤ 0.2 logMAR (including 79% of eyes having ≤ 0.0). Median absolute deviation from target refraction was 0.29 D (range 0.00 to 1.24) and thus acceptable.

Mean residual refractive astigmatism was -0.56 ± 0.34 D and was within ± 0.5 D in 67% of eyes and within 1.00 D in 92% of eyes. Vector analysis showed a significant reduction in magnitude of the J0 vector after surgery ($p < 0.0001$) with a mean absolute change in astigmatism of 0.77 ± 0.35 D. The mean J45 vector did not change significantly from preoperative values, which can be explained by very low preoperative values (only 11 out of 40 eyes had oblique astigmatism).

The mean IOL rotation was $1.81^\circ \pm 1.87^\circ$ at 1.5-years. IOL rotation was $< 3^\circ$ in 85% of eyes, $\leq 5^\circ$ in 94% of eyes, and no IOL rotated more than 7.2° . Between 6-month and 1-year, the mean IOL rotation was minimal ($0.84^\circ \pm 0.82^\circ$) with IOL rotation $< 2^\circ$ in 96% of eyes. The outcome at 1.5-year was not significantly different from the 6-month visit ($1.99^\circ \pm 1.88^\circ$, $P = 0.58$), indicating that the IOL was stable 6 months postoperatively.

The vast majority of patients (96.0%) were satisfied or very satisfied with their vision and all stated that they would choose the same lens model again. All patients bilaterally implanted were totally free from glasses for distance vision (target for emmetropia) or near vision (target for myopia).

In terms of complications, one case presented intraoperatively with a broken IOL haptic and the lens was replaced during surgery. Postoperatively, one eye of a patient with bilateral IOL implantation developed severe retinal detachment between the 6-month and 1.5-year period resulting in significant decreased visual acuity. Occasional halos at night (but not during the day) were reported by two patients, with the discomfort being judged as very mild or moderate. One patient reported glare rarely during the day and rated this discomfort as very mild. None of the eyes required Nd:YAG laser capsulotomy (due to PCO) or secondary IOL repositioning (due to IOL rotation) during the course of the study.

In conclusion the results of this clinical study show that implantation of the Torica-aA is a safe and effective procedure for the treatment of aphakia and the correction of low to moderate pre-existing astigmatism with excellent outcomes over 1.5 years of follow-up. Besides the successful correction of astigmatism its toric optical design enables a good visual performance which is coupled with a high degree of patient satisfaction.

Post-market study Tor-aA_Ros_Rota (Torica-aA – Rostock, Germany) [6]

The objective of this prospective study, initiated by the manufacturer, was to evaluate the rotational stability, astigmatism correction, visual performance, and patient satisfaction subsequent to the implantation of the Torica-aA IOL. Besides obtaining safety and performance data for “normal” eyes with an axial length (AL) of <25.0 mm (Group A), this study also provides data about the applicability in large eyes with an AL ≥25.0 mm (Group B). These data are especially of clinical interest with regards to the implantation of the Torica-aA in Asian eyes, where large ocular dimensions are quite common.

33 patients (50 eyes) with senile cataract and a mean age of ~67 years were recruited and implanted with the Torica-aA. The data of one patient with two eyes that were assigned to group A were excluded from the analysis as this patient died before the end of the study. Therefore, 23 eyes and 25 eyes were available for analysis in Group A and Group B, respectively. Target refraction was emmetropia for 35 eyes (73%), including all eyes from group A and 12 eyes from group B. The remaining 13 eyes (27%) of group B had a myopic target refraction of -2.5 D. Follow-up visits were scheduled at 1 day, 1 week, 1 month and 4 to 6 months after surgery (referred to as 6-month follow-up visit below).

At the final follow up, the mean absolute deviation of the subjective SE from the target refraction was 0.02 ± 0.52 D in Group A and 0.06 ± 0.48 D in Group B. Mean monocular CDVA was 0.06 ± 0.12 and 0.02 ± 0.07 , respectively. All eyes in both groups achieved a CDVA of 0.3 logMAR or better. Mean monocular UDVA was 0.14 ± 0.13 logMAR in group A and 0.08 ± 0.09 logMAR in group B (not considering eyes with a myopic target refraction).

The absolute rotation of the IOL was assessed at the “end of surgery” (EoS) and at all postoperative follow-up visits. In the early postoperative period, i.e. from EoS to 1-day and to 1-week, average rotation was higher than afterwards. In group A, the mean absolute rotation from EoS to 6-months postoperatively was $2.3 \pm 2.2^\circ$ with a maximum of 8° rotation in one eye. In group B, the value was $3.2 \pm 6.5^\circ$ with three maximum outliers of 30° , 14° and 10° rotation. Nevertheless, the outliers achieved excellent visual outcomes (uncorrected ≤ 0.22 logMAR, corrected 0.0 logMAR). There was no statistical difference between group A and group B regarding the mean absolute rotation from EoS to 6 months postoperatively. The Torica-aA fulfilled the ISO 11979-7 acceptance criteria for IOL rotation with total absolute IOL rotation of $\leq 5^\circ$ in 92% of eyes ($n=44/48$), $<10^\circ$ in 94% of eyes ($n=45/48$) and $<20^\circ$ in 98% of eyes ($n=47/48$) combined for group A and group B. The ISO 11979-7 acceptance criteria for IOL rotation are defined as a maximum of 90% for a rotation $<10^\circ$ and a maximum of 95% for a rotation $<20^\circ$.

Preoperative astigmatism was significantly reduced 6 months postoperatively, in group A from 2.26 ± 1.32 D to -0.72 ± 0.49 D ($p < 0.000$) and in group B from -1.98 ± 1.42 D to -0.65 ± 0.37 D ($p < 0.000$). The evaluation of the astigmatic power vectors (J0, J45) showed a significant reduction of the J0 vector (group A: $p=0.038$; group B: $p=0.037$) between the preoperative visit (corneal keratometric astigmatism) and the 6-month follow-up (refractive astigmatism). The preoperative magnitude of the J45 vector was comparably low and its postoperative reduction did not reach statistical significance in both groups.

No serious adverse events occurred during the study. A minor intraoperative complication was observed when the haptics of one IOL broke off during release from the injector, likely due to a user error. Four cases of transient cystoid macular edema were observed and successfully treated. Nd:YAG laser treatment due to posterior capsule opacification (PCO) was required for only 2 eyes across both groups during the complete follow-up period.

Patient-reported subjective satisfaction and experience of dysphotopsia (halos or glare) were assessed using a questionnaire. All patients were “highly satisfied” or “satisfied” with their vision and would choose the same IOL again. Most patients never experienced dysphotopsia (halos or glare) during the day or at night.

In conclusion, the results from the Tor-aA_Ros_Rota study confirmed that Torica-aA implantation is a safe and effective procedure for the treatment of aphakia and the correction of moderate astigmatism with good visual results and a good safety profile in terms of IOL rotational stability and complication rate. Furthermore, the data show that the implantation is also effective for long eyes with an AL ≥ 25.0 mm. This eye condition is often seen in Asian eyes and is thus of clinical relevance. The size of the data set was too small to analyze correlations, but the highest outliers in terms of rotation occurred in the AL ≥ 25 mm group, which could indicate a possible influence of AL in certain cases.

5.3.2. Literature review

The 18 articles identified through a systematic literature review are of various levels of usefulness to assess the safety and performance of the devices of the Monofocal Aspira-aA Family [7-24].

In 16 articles, a total of 1,716 eyes were implanted with the Aspira-aA or Aspira-aAY, of which 1089 eyes were evaluated in a retrospective collection of practice data [20]. Patients were followed over periods of >2 months and up to 1 year postoperatively. Overall, the data support a very good clinical performance of both IOLs and indicate a high degree of safety, although the most follow-up periods reported here are too short to exclude adverse effects that may occur over a long term (>1 year). When compared to an aberration-neutral or aberration-correcting monofocal IOL (Akreos Adapt Advanced or Tecnis ZCB00) [8, 11] or to two multifocal IOLs (Oculentis Mplus and LS-313 MF30) [12, 18], the aberration-free Aspira-aA displays a distance visual performance which is equal to, or even better to other IOLs. In addition, total and internal high-order aberrations were reported to be significantly lower in patients implanted with the Aspira-aA compared to the ones from the multifocal group, resulting in a better objective visual quality in the Aspira-aA patients. In comparison to another aberration-neutral monofocal IOL, the total, coma, spherical and IOL aberrations were comparable [8]. One publication compared the Aspira-aA to the HumanOptics Aspira-aXA, a hydrophilic acrylic IOL with an enlarged optic [21]. In direct comparison the Aspira-aXA showed a reduced frequency and extent of dysphotopsias. However, the overall frequency and extent of dysphotopsias decreased over time resulting in a general high patient satisfaction for both IOLs. Even when implanted in patients with retinal comorbidities either in a combined cataract and PPV surgery (called: phacovitrectomy) or in consecutive surgeries, the Aspira-aA showed a good operative and visual recovery [17, 22]. Good clinical performance and safety were also achieved when implanted in the ciliary sulcus (off-label use) of patients with PPV [17] as well as when combined with a sulcus-fixated multifocal add-on lens [7, 9]. One large-scale collection of real-world data (1089 eyes) revealed that the distance visual performance of the Aspira-aAY is comparable to the current gold-standard IOLs on the market [20]. Moreover this article reported a favourable clinical safety profile of the Aspira-aAY, with secondary cataract (PCO) occurring in only approx. 2.5% of cases, which was comparable to or even lower than the rates seen with other gold-standard IOLs.

Two publications dealt with the implantation of the standard-sized Aspira-aA instead of a custom-made device into nanophthalmic eyes [10, 15]. Here, a relatively high risk of intraoperative IOL cracks/breakage during implantation was reported [10], which was presumably due to the difficulties to inject an IOL with a bulky optic (large dioptric powers used: +39 to +58 D [10, 15]). As both reports were based on a very small number of subjects (5 eyes/3 patients and 2 eyes/1 patients), these results do not question the safety of the device. Moreover, it should be emphasized that the overall results were satisfactory, considering the special characteristics of nanophthalmic eyes, as manifest refraction improved greatly in all but one eye [10, 15]. In addition, no intra- or postoperative complications and a stable IOP were noted in the case report [15]. In two publications, 46 eyes were implanted with the Aspira-aAY in combination with the HumanOptics Artificial Iris [16, 19]. The results revealed an improvement in contrast sensitivity, CDVA and patient satisfaction, which is likely a combined effect of iris repair and IOL implantation. Taken together, the results indicate that the Aspira-aAY is safe and effective when implanted along with an iris prosthesis.

Overall, the results clearly demonstrate the safety and good clinical performance, even in patients with comorbidities that lead to challenging demands on the device under consideration.

Regarding the Torica-aA/-aAY two articles (one clinical trial, one case report; 36 eyes/patients in total) provide evidence for the high rotational stability over a period of up to 3 months, thus

confirming the findings made in the two post-market studies Asp-aA_Wie-Rota with the Aspira-aA and the Tor-aA_Lin_Rota with the Torica-aA. Furthermore, the side-by-side comparison with the hydrophobic IOL AcrySof IQ Toric [14] showed that the Torica-aAY was on a par with this well-established device and gold standard lens in terms of rotational stability. No intraoperative complications occurred in any of the cases. However, postoperative rotation which required repositioning was reported in two cases, in one likely caused by premature Nd: YAG laser treatment [13]. The ability to correct astigmatism and to provide good visual acuity was confirmed in the case report. In conclusion, the results clearly support the safety of the device under consideration and provide limited support (since just from a single case) for its performance.

5.4. Overall summary of the clinical performance and safety

IOLs are intended to replace the focusing power of the natural lens thereby restoring the patient's vision. Accordingly, the main criterion for the success of cataract surgery is the long-term visual result. Since monofocal, non-accommodative lenses like the Monofocal Aspira-aA Family lenses focus the light at a single, fixed distance it is expected that the patients will still require glasses for activities at certain distances, e.g. reading or screen handling. Therefore, the main performance indicators for the devices under consideration are the obtained distance visual acuity and the reduction of pre-existing astigmatism for the toric models. Key safety indicators are the rate of complications and, in the case of the toric IOLs, the rate of IOL rotations.

PERFORMANCE ANALYSIS

Based on the results of the post-market studies for the Monofocal Aspira-aA Family (Table 1), a performance analysis was performed in accordance with the ISO 11979-7:2018-08. This analysis included all available eyes during a period ranging from 2 months to 3 years post-operatively. **Key performance CDVA outcomes are presented in Table 2.** All healthy eyes, except one at the 3-month follow-up of study Asp-aXA_Ora_ND, achieved a CDVA of 0.3 logMAR or better. Even among subjects with pseudoexfoliation syndrome (PEX; study Asp-aA_Lub_PEX), 94% and 93.4% of eyes achieved 0.3 logMAR or better at 3 months and 1 year postoperatively, respectively, matching the primary safety and performance endpoints (SPE) rates. Thus, the acceptance criterion, which is based on the SPE rates defined in ISO 11979-7:2018-08 for posterior-chamber lenses (92.5% overall), was met in every single study.

Similar findings regarding monocular CDVA in healthy eyes have been reported in the literature [8, 11, 13, 18] (based on mean and standard deviation since the percentage of eyes reaching 0.3 logMAR has not been reported).

The studies also indicate that the visual performance (even colour perception) was not affected by the presence of a blue-light filter (Aspira-aA vs. Aspira-aAY; post-market study Asp-aA_Erl_VA). Moreover, the visual performance of the polished and non-polished variant was similar, although the latter slightly outperformed the former in terms of contrast sensitivity. The level of patient satisfaction was overall high and the contrast sensitivity was good in all studies, in which such data were collected.

Table 2: Evaluation of performance CDVA outcomes according to DIN EN ISO 11979-7:2018-08

Endpoint	Acceptance criterion	Clinical outcome	Data source
Monocular CDVA (% of eyes reaching 0.3 logMAR or better)	≥ 92.5% of the eyes should achieve CDVA 0.3 logMAR or better [SPE rate]		Post-market studies:
		Aspira-aA/-aAY	
		98.4% [-0.18 to 0.40] (n=63) at 3-month 100% [-0.18 to 0.22] (n=62) at 1-year 100% [-0.18 to 0.30] (n=58) at 3-year	Asp-aXA_Ora_ND
		100% [0.00 to 0.30] (n=18) at 2-month	Asp-aA_Hal_Rota
		100% [-0.10 to 0.10] (n=28) at 3-month 100% [-0.10 to 0.10] (n=26) at 2-year	Asp-aA_Gla_Poli Polished Aspira-aA
		100% [-0.10 to 0.10] (n=31) at 3-month 100% [-0.10 to 0.10] (n=29) at 2-year	Non-polished Aspira-aA
		100% [-0.20 to 0.30] (n=70) at 3-month 98.5% [-0.20 to 1.0] (n=67) at 1-year	Asp-aA_Lub_PEX Non-PEX (healthy) patients
		94% [-0.20 to 0.90] (n=67) at 3-month 93.4% [-0.20 to 0.9] (n=61) at 1-year	PEX patients
		100% [0.00 to 0.22] (n=19) at 3-month 100% [-0.08 to 0.20] (n=19) at 1-year	Asp-aA_Erl_VA Aspira-aA
		100% [-0.08 to 0.10] (n=28) at 3-month 100% [-0.08 to 0.22] (n=19) at 1-year	Aspira-aAY
		100% [-0.10 to 0.10] (n=26) at 3-month 100% [-0.10 to 0.10] (n=26) at 1-year	Inv-aA_Gla_SA
		Torica-aA/-aAY	
		100% [-0.20 to 0.22] (n=40) at 3-month 97.4% [-0.28 to 0.36] (n=39) at 1.5-year	Tor-aA_Lin_Rota
		100% [-0.10 to 0.22] (n=50) at 1-month	Tor-aA_Ros_Rota AL <25.0 mm (min 21.6 mm)

		100% [0.00 to 0.20] (n=50) at 1-month	AL $\geq$ 25.0 mm (max 28.1 mm)

Regarding the performance endpoint “reduction of pre-existing astigmatism”, mean (absolute) astigmatism in the post-market study Tor-aA_Lin_Rota was reduced from $1.62 \text{ D} \pm 0.39$ (topographic measurement) preoperatively to $0.56 \pm 0.34 \text{ D}$ (refractive measurement) 6 months after surgery. Results were similar 1.5 years postoperatively. Vector analysis showed a significant reduction in magnitude of the J0 vector after surgery ($p < 0.0001$) with a mean absolute change in astigmatism of $0.77 \pm 0.35 \text{ D}$. The mean J45 vector did not change significantly from preoperative values, which can be explained by very low preoperative values. Similarly, in the study Tor-aA_Ros_Rota [6], the J0 vector was significantly reduced at 6 months postoperatively, while the postoperative reduction of the J45 vector (starting from comparably low preoperative values) did not reach statistical significance. In addition, the results from this study showed that the correction of astigmatism was similar irrespective of whether the IOL was implanted in eyes with an AL $< 25.0 \text{ mm}$ or an AL $\geq 25.0 \text{ mm}$. A reduction of the refractive astigmatism subsequent to the implantation of the Torica-aA has also been demonstrated in a case report by Kaindlstorfer et al. [13] (from 3.75 D preoperatively to 0.25 D at one week after surgery). These data are confirmed by results with the related model ToricaDiff-aA (same toric surface design as Torica-aA but additionally multifocal properties), where residual refractive astigmatism was significantly reduced from preoperative values (post-market study TDI-aA_Hal_VA; data not shown). In conclusion, the above findings show that pre-existing regular, stable corneal astigmatism can be corrected with the Torica-aA/-aAY as intended by the manufacturer.

The scientific publications regarding Aspira-aA/-aAY add the following information. In a randomized controlled trial [11] the aberration-free lens Aspira-aA under consideration was compared to a popular aberration-correcting monofocal IOL (Tecnis ZCB00). The results indicate that the Aspira-aA achieved an at least equal visual performance compared to the Tecnis IOL. Although the measurements of optical quality were slightly in favor of the aberration-correcting IOL at large pupil sizes, the Aspira-aA outperformed the other lens in terms of contrast sensitivity at high spatial frequencies. This was the case even under mesopic conditions, i.e. with dilated pupil, where spherical aberrations should have the strongest effect. When the aberration-free Aspira-aA was compared to another aberration-free monofocal IOL (Akreos Adapt Advanced), both IOLs performed well, with no differences in visual performance in terms of distance visual acuity as well as total, coma, spherical and IOL aberrations [8]. In two studies, the monofocal Aspira-aA was compared to multifocal IOLs (Mplus LS MF30 and LS-313) [12, 18] and achieved an UDVA that was comparable to the multifocal IOLs. One large-scale collection of real-world data confirmed the good performance of the Aspira-aA in a large number of eyes (1,089), with UDVA values comparable to competitors on the market, including the gold standard lenses from Alcon, Rayner and Hoya [20]. Two publications [7, 9] demonstrate furthermore that the Aspira-aA can be successfully combined with a sulcus-fixated multifocal add-on lens (mean/median UDVA $\leq 0.1 \text{ logMAR}$), albeit implanting a single capsular-bag multifocal IOL instead appears as a more straightforward choice today. Several publications report on the outcome of Aspira-aA/-aAY implantation in eyes with certain ocular pathologies (nanophthalmus, PPV or aniridia) [10, 15-17, 19, 22], sometimes implanting the IOL outside of the capsular bag (off-label use). Although the visual capacity of those eyes is often limited, the publications showed that a CDVA of 0.0 logMAR can be reached in principle.

In the case of the Torica-aA, a case report [13] confirmed the highly satisfying results in terms of visual performance. Although focusing on a postoperative complication (IOL rotation caused Nd:YAG capsulotomy two weeks after implantation in the absence of posterior capsule opacification), the early postoperative data demonstrated good visual acuity (UDVA: 0.0

logMAR, CDVA: -0.1 logMAR), good refraction (-2.50 D, as targeted) and a low residual refractive astigmatism (-0.25 D).

Taken together, these data show that the Monofocal Aspira-aA Family fulfil the pre-defined acceptance criteria and achieve the performance intended by the manufacturer as required by the MDR Annex I, when used according to the intended use.

SAFETY ANALYSIS

Several complications and undesirable effects have been reported in the literature for foldable posterior-chamber IOLs (including toric IOLs = TIOLs) but the incidence is overall very small, particularly, if one takes into account that eligible patients often have several ocular comorbidities. In most cases these side-effects are related to the surgery and not device-related. Regarding the Monofocal Aspira-aA Family lenses, the reported outcomes show only a small percentage of complications and undesirable effects which have been already described in the literature.

Based on the results of the post-market studies for the Monofocal Aspira-aA Family (Table 1) a safety analysis was performed. The primary safety endpoint was met as the cumulative and persistent adverse events (complication) rates were within limits compared to the historical SPE rates reported in ISO 11979-7:2018. **Key safety outcomes are presented below in Table 3.**

In terms of rotational stability, data analysis showed that the Monofocal Aspira-aA Family IOL platform meets the safety requirements as described in DIN EN ISO 11979-7:2018-08 (Table 4).

Table 3: Calculation of complication rates according to DIN EN ISO 11979-7:2018-08 [25]

Endpoint	Acceptance criterion (SPE rate) %	Threshold rate % for n = 100	Threshold rate % for n = 300	Clinical outcomes % for n = 458	Data source
Cumulative					<u>Post-market studies Monofocal Aspira-aA Family (n = 458 eyes):</u> <u>Aspira-aA/-aAY (n = 370)</u> Asp-aA_Erl_VA Asp-aA_Hal_Rota Inv-aA_Gla_SA Asp-aA_Gla_Poli Asp-aA_Lub_PEX Asp-aXA_Ora_ND <u>Torica-aA/-aAY (n = 88)</u> Tor-aA_Lin_Rota Tor-aA_Ros_Rota
Cystoid macular oedema (CME)	3.0	8.9	6.0	1.75	
Hypopyon	0.3	3.0	1.8	0.00	
Endophthalmitis	0.1	3.0	1.0	0.00	
Lens dislocated from posterior chamber	0.1	3.0	1.0	0.44*	
Pupillary block	0.1	3.0	1.0	0.00	
Retinal detachment	0.3	3.0	1.8	0.66*	
Secondary surgical intervention	0.8	4.2	2.6	0.00	
Persistent					
Corneal stroma oedema	0.3	3.0	1.8	0.00	

Cystoid macular oedema (CME)	0.5	4.2	2.2	0.22	
Iritis	0.3	3.0	1.8	0.00	
Raised IOP requiring treatment	0.4	4.2	1.8	0.00	

* value is not statistically significantly higher than SPE rate (threshold for a study with n=300 is not exceeded)

Table 4: Calculation of rotational stability according to DIN EN ISO 11979-7:2018-08

Endpoint	Acceptance criterion %	Clinical outcomes % for n = 140 analyzable eyes	Data source
Rotation < 10°	90	95.7	Post-market studies (n = 140): Asp-aA_Wie_Rota Tor-aA_Lin_Rota Tor-aA_Ros_Rota
Rotation < 20°	95	98.6	

Data for mid-term follow-up (1-6 months postoperatively) vs end of surgery (within 24h) are shown.

Regarding complaints related to complications, very few have been reported for the Monofocal Aspira-aA Family lenses and the only type of serious incidents concerned few cases of IOL calcification. In all cases, when an IOL calcification was confirmed, it appeared to be caused by secondary calcification due to environmental circumstances, i.e. they were not related to the IOL itself.

Clinical studies with HumanOptics IOLs have been conducted using both the Compact Line as well as the Safeloader-packaged IOLs and confirmed the performance and safety of the tested injectors. However, common injector-associated complications, such as trapped haptics or misconfiguration of the IOL haptics after injector release, have also been observed in a few cases. These were typically associated with handling errors.

In summary, one can conclude that the devices under consideration Aspira-aA/-aAY and Torica-aA/-aAY represents the current state of the art and fulfil the pre-defined acceptance criteria and the requirements regarding safety as set out by the MDR Annex I. The devices are safe for implantation in the capsular bag for the visual correction of aphakia (and corneal astigmatism), and no unanticipated risks associated with the IOLs were identified.

Since 2003 HumanOptics has placed several million hydrophilic acrylic IOLs on the market. Over this period of more than 20 years no feedback has been received relating to material degradation other than very rare cases (< 5 ppm/year) of suspected IOL opacification (calcification). The latter is caused by crystallisation of calcium on the IOL surface, usually resulting from environmental circumstances, i.e. it does not constitute a degradation due to aging processes of IOLs manufactured by HumanOptics. The available information confirms that HumanOptics hydrophilic acrylic IOLs are stable over a period of at least 20 years (based on simulated aging tests according to ISO 11979-5). It can be expected that the device is stable indefinitely from the implantation date over the lifetime of the patient.

BENEFIT-RISK ANALYSIS

The main **benefits** of cataract surgery and IOL implantation in general are improved visual function and quality of life. Important aspects involved are:

- Improved vision (visual acuity, contrast sensitivity, depth perception, etc.)
- Reduction in falls, fractures and vehicle crashes
- Reduced eyeglass dependence (through correction of refractive errors and, in case of TIOLs, of astigmatism)
- Improvement of quality of life (activity, confidence, mood)
- Reduction of anxiety and depression

In terms of improved vision, the different IOL types currently available on the European market can address single or a combination of vision problems. For instance, **monofocal IOLs** are designed to improve vision at just one set distance (generally far vision); **multifocal IOLs** are designed to additionally improve vision at near and/or intermediate distances and **toric IOLs** specifically correct pre-existing corneal astigmatism (they are available as monofocal toric or multifocal toric IOLs).

Two types of intended clinical benefits can be attributed to the Monofocal Aspira-aA Family IOLs. The patient benefit “**restored vision**” is based on achieving good CDVA after IOL implantation, which has been demonstrated in the previous chapters. Accordingly, in the post-market studies where patient satisfaction was analyzed, most patients reported to be “very satisfied” or “satisfied” with the vision (studies Asp-aA_Erl_VA, Asp-aA_Hal_Rota, Asp-aXA_Ora_ND, Tor-aA_Lin_Rota and Tor-aA_Ros_Rota; range: 77.8% to 100%). In the remaining cases, ocular diseases (e.g. macular edema), not IOL-related, were the likely cause for the lack of high levels of satisfaction.

In the case of TIOLs (Torica-aA/-aAY), the “**effective reduction of pre-existing astigmatism**” can lead to another benefit i.e a reduced dependence on eyeglasses. In line with this, all patients bilaterally implanted with the device in the post-market study Tor-aA_Lin_Rota achieved complete spectacle independence for the targeted vision distance (far or near).

It is understood that the extent of benefits depends on the degree of patients’ preoperative discomfort and ocular comorbidity which may limit the performance of the device as, for instance, in patients with retinal diseases. In the case of patients suffering from severe cataract, who are considerably visually impaired, the implantation of IOLs can mean a huge relief and improvement in quality of life.

Regarding **risks**, very few complications have been reported for the Monofocal Aspira-aA Family, which is in line with the low frequency of complications reported for cataract surgery in general. Nonetheless, complications can always occur but only in extremely rare cases complications may lead to serious and permanent impairment of the patient. Some of the most relevant risks are ocular inflammation and cystoid macular edema.

As mentioned earlier, postoperative complications are rarely device-related but rather depend on the preoperative ocular comorbidity, which means that IOLs should be implanted into risky patients only after careful assessment of the benefit/risk ratio. Furthermore, the success of cataract surgery is influenced by the surgical technique, equipment used, and the surgeon’s experience. In addition, there is a learning curve when performing implantations with a IOLs with novel designs. Nevertheless, the implantation of IOLs has to be considered as a routine procedure for the target user group, i.e. ophthalmic surgeons. Since the essential aspects of the handling and implantation process are adequately described in the IFU, a special training for the Monofocal Aspira-aA Family is not required. However, device handling errors cannot be completely excluded but appear to occur very rarely. A risk that specifically affects TIOLs

compared to Non-TIOLs is postoperative rotation which might sometimes necessitate a secondary procedure to realign the lens. Currently, the respective data for the devices in question compared favourably to the thresholds set by the applicable standards.

In summary, a large body of evidence proves that the Aspira-aA/-aAY and Torica-aA/-aAY are reliable products and that the remaining risks are well-controlled.

A positive benefits-versus-risks profile can be assigned to the devices under consideration Aspira-aA/-aAY and Torica-aA/-aAY, particularly considering the current knowledge and state of the art in the medical fields concerned. There is sufficient evidence that the intended clinical benefits are achieved. The probable benefits related to restoring vision secondary to cataract extraction outweigh the probable risks for the intended use of the device. The benefits-versus-risks profile of the new Aspira-aA+ and Aspira aAY+ is non-inferior to the equivalent Aspira-aA/-aAY since the risks are the same and the benefits are the same or may be greater.

5.5. Ongoing or planned post-market clinical follow-up

Clinical evidence for the safety and performance of the Monofocal Aspira-aA Family has been based on the manufacturer's own clinical data and literature results, showing that implantation of the IOLs in the capsular bag is a safe treatment for the correction of aphakia and pre-existing astigmatism for the toric version, providing good and stable visual performance as well as high levels of patient satisfaction.

Furthermore, the Aspira-aA/-aAY and Torica-aA/-aAY are well established in the medical field of cataract surgery, and have been successfully marketed for more than 10 years without apparent risks which would have arisen in the course of the PMS activities of the manufacturer. Therefore, we could not identify any unanswered question or uncertainties which should be subject to a future specific PMCF activity for the Monofocal Aspira-aA Family lenses. Nonetheless, as required by the MDR in cases where CE marking is based on product equivalence, HumanOptics will perform a post-market study with the Aspira-aA+ (or Aspira-aAY+) to evaluate its visual performance.

In addition, a user survey is currently being conducted to continuously collect user data referring to safety, performance and usability of HumanOptics products.

6. Possible diagnostic or therapeutic alternatives

Implantation of an IOL is the only well-tolerated treatment available to replace the focusing power of the natural lens after surgical extraction. The preferred placement of an IOL is in the capsular bag. Mild/early cataracts may also be managed by prescribing new eyeglasses to correct a lens-induced change in refractive errors. High-power eyeglasses (still in use in developing countries) and contact lenses (only for high myopic eyes) are alternative treatments for aphakia but do not achieve visual results as good as IOLs.

7. Suggested profile and training for users

IOLs must be handled by health professionals and implanted by ophthalmic surgeons. A high level of surgical skill is required for IOL implantation. The surgeon should have observed and/or assisted in a large number of implantations and successfully completed one or more courses on IOL implantation before attempting a surgery. It is also strongly recommended that the surgeon and his/her team receive a training on the specific product and accessories by the respective sales representative before the first implantation of the device. Application guides (PDF) may be available on request. Product application videos can be found on www.humanoptics.com in the “application” section for the respective device.

If you have further questions on the usage of the device, please contact your local sales representative or the manufacturer (application@humanoptics.com).

8. Reference to any harmonised standards and CS applied

Common specification(s) to comply with, if applicable:
N/A
Harmonised standards to apply, if applicable:
DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
DIN EN ISO 14971:2020-07 Medical devices – Application of risk management to medical devices
Relevant adopted monographs of the European Pharmacopoeia
Not applicable

9. References

- [1] Kühle M. Comparison of visual function with aspheric yellow, aspheric clear and spherical clear intraocular lenses. *J Emmetropia* 2013;4:123-130.
- [2] Schrecker J, Langenbucher A, Seitz B, Eppig T. First results with a new intraocular lens design for the individual correction of spherical aberration. *J Cataract Refract Surg.* 2018;44(10):1211-1219.
- [3] Rozegnal-Madej A, Zarnowski T. Biocompatibility of acrylic intraocular lenses in eyes with pseudoexfoliation syndrome - Aspira-aA (HumanOptics AG) vs AcrySof SA60AT (Alcon Lab) – one year follow up. 38th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; virtual Session: Pseudophakic IOLs: Monofocal). 2020.
- [4] Gyöngyössy B, Jirak P, Schönherr U. Rotational stability and patient satisfaction after implantation of a new toric IOL. *Eur J Ophthalmol.* 2016;26(4):321-327.
- [5] Gyöngyössy B, Jirak P, Schönherr U. Long-term rotational stability and visual outcomes of a single-piece hydrophilic acrylic toric IOL: a 1.5-year follow-up. *Int J Ophthalmol.* 2017;10(4):573-578.
- [6] Rusch W, Anwar M, Altas C, Steirat M, Fuchsluger T, Walckling M, et al. Influence Of Axial Length On Rotational Stability And Other Clinical Outcomes After Implantation Of A Hydrophilic Toric Intraocular Lens. 41th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; Session: IOLs and Evaluation Testing). 2023.
- [7] Schrecker J, Blass S, Langenbucher A. Silicone-diffractive versus acrylic-refractive supplementary IOLs: visual performance and manual handling. *J Refract Surg.* 2014;30(1):41-48.
- [8] Rajabi MT, Korouji S, Farjadnia M, Naderan M, Rajabi MB, Khosravi B, et al. Higher order aberration comparison between two aspherical intraocular lenses: MC6125AS and Akreos advanced optics. *Int J Ophthalmol.* 2015;8(3):565-568.
- [9] Liekfeld A, Ehmer A, Schröter U. Visual function and reading speed after bilateral implantation of 2 types of diffractive multifocal intraocular lenses: Add-on versus capsular bag design. *J Cataract Refract Surg.* 2015;41(10):2107-2114.
- [10] Chang JS, Ng JC, Chan VK, Law AK. Cataract Surgery with a New Fluidics Control Phacoemulsification System in Nanophthalmic Eyes. *Case Rep Ophthalmol.* 2016;7(3):218-226.
- [11] Lasta M, Miháltz K, Kovács I, Vécsei-Marlovits PV. Effect of Spherical Aberration on the Optical Quality after Implantation of Two Different Aspherical Intraocular Lenses. *J Ophthalmol.* 2017;2017:8039719.
- [12] Liu S-Y, Yang X. Visual outcome with the Oculentis Mplus intraocular lens. *Int Eye Sci.* 2017;17(11):2134-2136.
- [13] Kaindlstorfer C, Kneifl M, Reinelt P, Schönherr U. Rotation of a toric intraocular lens from neodymium:YAG laser posterior capsulotomy. *J Cataract Refract Surg.* 2018;44(4):510-511.
- [14] Shukhaiev SV, Ulitina AY, Shukhaeva EN. Comparative evaluation of rotational stability of hydrophobic (AcrySof IQ Toric) and hydrophilic (Torica-aAY) toric IOLs (Russian full-text; English abstract and translation available). *Point of View East – West.* 2019;1:16-19.
- [15] Naujokaitis T, Scharf D, Baur I, Khoramnia R, Auffarth GU. Bilateral implantation of +56 and +58 diopter custom-made intraocular lenses in patient with extreme nanophthalmos. *Am J Ophthalmol Case Rep.* 2020;20:100963.
- [16] Mayer CS, Hoffmann AM, Prahs P, Reznicek L, Khoramnia R. Functional outcomes after combined iris and intraocular lens implantation in various iris and lens defects. *BMC Ophthalmol.* 2020;20(1):370.

- [17] Sun Y, Zhang Z, Wei Y, Chi W, Zhang S. Intraocular Lens Fixation Technique Without Corneal Incision in Minimally Invasive Vitrectomized Eyes. *Ophthalmol Ther.* 2022;11(2):729-737.
- [18] Yang J, Hai-Bo C, Ming-Bing Z, Jun Y. Error analysis of visual quality in segmented refractive multifocal intraocular lenses implanted eyes. *International Eye Science.* 2022;447-451.
- [19] Mayer CS, Baur ID, Storr J, Khoramnia R. Traumatische Aniridie: Konservativer oder chirurgischer Therapieansatz? *Ophthalmologe.* 2022;119(4):395-399.
- [20] Chuprov AD, Startseva MI, Original, Чупров А, Старцева М. Comparative analysis of clinical and economic efficiency in implantation of various IOL models; Original Сравнительный анализ клинико-экономической эффективности при имплантации различных моделей ИОЛ. *Journal "Medicine"; Original Журнал «Медицина»* №. 2022;2:1.
- [21] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. *J Cataract Refract Surg.* 2022;48(1):75-82.
- [22] Liu X, Shao Y, Lin H, Liu C, Shen J, Zhang L, et al. Corneal densitometry: an innovative method to quantitatively evaluate corneal changes after phacovitrectomy. *BMC Ophthalmol.* 2023;23.
- [23] Ma Y, Liu D, Lv X, Lin T. Study of wavefront aberration after multi-focal intraocular lens implantation. *International Journal of Frontiers in Medicine.* 2024;6(5).
- [24] Wlaź A, Kustra A, Aung T, Żarnowski T. Evaluation of changes of anterior segment parameters in patients with pseudoexfoliation syndrome after cataract surgery using anterior segment optical coherence tomography. *Sci Rep.* 2024;14(1):8279.
- [25] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia

Part 2 – SSCP intended for patients

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.


The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your ophthalmologist in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions For Use to provide information on the safe use of the device.

This SSCP is reviewed at least annually and will be updated if needed.

1. Device identification and general information

Device trade name(s)	• Aspira-aA • Aspira-aAY • Aspira-aA+ • Aspira-aAY+ • Torica-aA • Torica-aAY (collectively referred to as “Monofocal Aspira-aA Family”)
Manufacturer’s name and address	HumanOptics Holding AG Spardorfer Str. 150 91054 Erlangen Germany
Basic UDI-DI(s)	• Aspira-aA/-aAY and Aspira-aA+/-aAY+: 04049154_PC_M2_H2_O2_BU • Torica-aA/-aAY: 04049154_PC_M2_H2_O3_BX
Year when the first certificate (CE) was issued covering the device	• Aspira-aA/-aAY: 2007 • Torica-aA/-aAY: 2012 • Aspira-aA+/-aAY+: Not applicable at this stage, as the product has been submitted for CE-marking.


2. Intended use of the device

Intended purpose	The device is a tiny intraocular lens (IOL) intended to be implanted in the capsular bag ² of your eye to replace the natural lens (crystalline lens).
Indication(s) and intended patient groups	The IOL is used to replace your natural lens after removal by an ophthalmic surgeon. Usually, the removal of your natural lens is because your lens has become cloudy with the age. This is called “cataract”. Cataract affects your vision.  Implantation of an IOL is indicated for adult patients only. Note: there are no clinical data available for children, pregnant or breastfeeding women, and immunocompromised persons in connection with the implantation of intraocular lenses from HumanOptics. Depending on your eye anatomy, lifestyle and preferences, different types of IOLs can be offered to you by your ophthalmologist (see section 3 “mode of action/operation principle”).
Contraindications	There are general contraindications for ocular surgery but no additional contraindications which are specific to the implantation of HumanOptics IOLs.

3. Device description

Description of the device	The Aspira-aA/-aAY, Aspira-aA+/-aAY+ and Torica-aA/-aAY are tiny artificial lenses made of plastic, called intraocular lenses (IOLs), which are implanted into the eye to replace the natural lens. They are designed to give you a clear view at a fixed distance (usually far distance), but the Aspira-aA+/-aAY+ may also help you see more clearly at intermediate distances (see section 3 “mode of action/operation principle”). The Torica-aA/-aAY is designed like the Aspira-aA/-aAY but can additionally reduce or eliminate corneal astigmatism (see section 3). In most cases, the IOL will be used to replace the clouded lens (cataract) of a person. If lens opacity increases, surgery is the only treatment available to restore good vision. Currently, no alternative or drug therapies are available. The IOLs have a diameter of just a few millimeters and fit perfectly into the eye due to their shape and the flexibility of their material. IOLs are composed of an “optic” (the round part of the IOL) and the “haptics”. The optic has a dioptric power that permits the focusing of light onto the retina and the haptics which are attached to the edge of the optic will help maintain the IOL in place in the eye.

² **Capsular bag:** The capsule which the natural lens sits in.

	The IOLs are available in colourless or yellow versions. The yellow variant (with “-Y”) absorbs a portion of the blue light, which may help protect the retina³ over the long term.  Schematic lateral and frontal view of the Aspira-aA (left) and Torica-aA (right). The dimensions are given in mm.
Mode of action / operating principle	The device is a static refractive lens substituting the natural crystalline lens. It focuses incoming light onto the retina. Your doctor may also having you diagnosed with corneal astigmatism (irregular curvature of the cornea). Astigmatism is not a disease and relatively common: about 20% of cataract patients have moderate to severe astigmatism, which can cause blurred or distorted vision. A good option to reduce or eliminate your corneal astigmatism is a toric IOL such as the Torica-aA/-aAY.
Materials or substances in contact with the patient's tissues	All IOLs manufactured by HumanOptics are made of biocompatible plastic (hydrophilic acrylic material) with an integrated UV absorber. Models with names ending with “Y” additionally contain a yellow coloring agent intended to filter the blue light.
Information about medicinal substances in the device, if any	Not applicable.
Accessories	None which are relevant to patients

³ **Retina:** The innermost, light-sensitive layer of tissue of the eye that transmits that information to the brain via the optic nerve

4. Risks and warnings

Contact your ophthalmologist if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your ophthalmologist if needed.

How potential risks have been controlled or managed	Potential risks have been identified and controlled according to international standards that have been developed by experts from many countries. The corresponding document(s) comprising rules, guidelines and processes are approved and published by a globally recognized organisation.
Remaining risks and undesirable effects	As with any surgical intervention, routine cataract surgery and associated IOL implantation are not entirely risk-free. However, this procedure is among the most common and most successful ones performed nowadays in the world. There are no known risks which are specific to the Aspira-aA/-aAY, Aspira-aA+/-aAY+ and Torica-aA/-aAY and which have not been previously reported for similar IOLs. The following side effects may occur following the IOL implantation but they normally disappear shortly after the surgery: • Foreign body sensation caused by the small incision through which the lens was inserted • Halos, which can form around sources of light, particularly at dusk or in darkness • You may experience increased sensitivity to light in the first few weeks after surgery The following symptoms may indicate complications associated with an IOL implantation. If these occur, you should consult your ophthalmologist: • Sudden or gradual deterioration in vision / blurred vision • Severe redness of the eye • Eye pain Months or years after IOL implantation, your vision might deteriorate and become blurry or hazy. This might be due to posterior capsule opacification (PCO), which occurs when a cloudy layer of scar tissue forms behind the IOL. This condition is not dangerous and can be treated by performing an opening in the capsule (or bag) that contains your IOL. Opening the capsule allows light to enter the eye and improve the clarity of the vision, as long as the rest of the eye is healthy. This procedure is called "Nd:YAG laser capsulotomy". It is fast (a few minutes), painless and is an easy way for restoring your vision. There is no other treatment for PCO and the condition will not improve by itself.
Warnings and precautions	Patients with certain eye diseases and/or malformations have higher-than-normal risks to develop complications during and after IOL implantation. Before the surgery, your ophthalmologist will perform an assessment of your eye and will discuss with you anything to do with your operation. This includes the risks and benefits of surgery and how long it will take to fully recover after the surgery. Difficulties during surgery may also require choosing another lens than planned.
Summary of any field safety corrective action (FSCA including FSN), if applicable	No devices have been recalled since the Aspira-aA/-aAY, Aspira-aA+/-aAY+ and Torica-aA/-aAY have been sold on the market.

	In July 2024, HumanOptics issued a Field Safety Notice (FSN)⁴ for healthcare professionals. The notice affected all devices in a certain packaging which had been manufactured in the period January 1, 2020 to June 28, 2024. The manufacturer had observed very rare cases of damaged outer blister packaging after transportation. The associated hazard potential is considered to be very low and no patient injury has happened to the manufacturer's knowledge. The FSN served to inform the users of this possible risk and to tell them how a damaged packaging can be recognized. No specific action is needed in the case of patients who have been treated with potentially affected products. Meanwhile, the product packaging has been improved.

⁴ A field safety notice (FSN) is an important means of communicating safety-related information to users.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1. Clinical background of the device

Cataract, i.e. the clouding of the natural lens of the eye, is a very common condition and most commonly affect adults as a result of ageing. Cataract causes deterioration of your vision and usually gets slowly worse over time. Surgery to replace the cloudy lens is the only way to improve your vision. Implantation of an IOL is the only well-tolerated treatment available to replace the cloudy lens. This procedure that usually takes 30 minutes is typically associated with rapid visual recovery, good visual outcomes and minimal complications.

After you and your ophthalmologist have decided that you should have your cataract removed, the length of your eye and the shape of your cornea will be measured. These measurements will help your surgeon choose the best IOL for you. You will also discuss the various lens options available to you. Further information about cataract surgery can be found under:

<https://www.humanoptics.com/en/patienten/help-with-cataracts/>

After implantation, the IOL will normally remain in your eye for the rest of your life, provided there is no medical contraindication. The product lifetime of HumanOptics IOLs is proven for twenty years. Regular ophthalmological check-ups are recommended, as for patients with natural lenses.

5.2. The clinical evidence for the CE-marking

The letters 'CE' which appear on products or packages signify that products sold in the European Economic Area (EEA) have been assessed to meet high safety, health, and environmental protection requirements.

Clinical evidence for CE-marking IOLs should meet the requirements defined by international standards and is based on clinical data with the device, scientific literature and feedback from the market.

The safety and good performance of the Monofocal Aspira-aA Family has been shown *via* the results of several clinical studies and scientific articles. The implantation of the IOLs in the capsular bag is a safe treatment to restore vision after cataract surgery resulting in high levels of patient satisfaction. The Aspira-aA/-aAY and Torica-aA/-aAY are well established in the medical field of cataract surgery, and have been successfully marketed for more than 10 years without safety concerns. The newer model Aspira-aA+/-aAY+ is almost the same as the Aspira-aA/-aAY, so no clinical data had to be collected before the CE marking. Although it is thought of as a slightly improved version, it does not cause any extra risks.

In summary, the Aspira-aA/-aAY, Aspira-aA+/-aAY+ and Torica-aA/-aAY perform well in terms of visual acuity, astigmatism correction for the toric model, patient satisfaction and safety. The clinical outcomes are good and comparable to other modern IOLs on the market. The lenses show to be a safe and effective option for the surgical treatment of cataract.

5.3. Safety

The implantation of the CE-marked Aspira-aA/-aAY, Aspira-aA+/-aAY+ and Torica-aA/-aAY IOLs considerably improve vision in cataract patients. The CE marking for IOLs is mandatory and means that the manufacturer has verified that the products meet European safety, health or environmental requirements.

The side effects of IOLs are well-known and are similar for all IOLs on the market. They generally occur rarely and are not very harmful (see chapter 4 "Remaining risks and undesirable

effects”). Therefore, the benefit-risk ratio of the devices is very positive. The current clinical data are sufficient to demonstrate the performance and safety of the Monofocal Aspira-aA Family and therefore no future post-market studies are required.

The manufacturer continuously analyses complaints, other market feedback and new scientific publications. The results of these activities have been favorable. This further supports the safety of the devices, besides the other data which were evaluated for the CE marking.

6. Possible diagnostic or therapeutic alternatives

When considering alternative treatments, it is recommended to contact your ophthalmologist who can take into account your individual situation.

Currently, there are no drugs which can stop or reverse the clouding of the natural lens of the eye. In the case of a pronounced cataract, the clouded lens must therefore be removed and replaced by an IOL. The preferred placement of an IOL is in the capsular bag. If you have a mild/early cataract, your ophthalmologist may prescribe only eyeglasses until your cataract becomes more mature.

There are different types of IOLs on the market, each having specific advantages and disadvantages. As mentioned earlier, choosing the right IOL can be difficult and it depends on several factors including your eye condition and your lifestyle. The different options should be discussed with your ophthalmologist.

7. Suggested training for users

Not applicable to patients since the device is to be used only by healthcare professionals.

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Restrisiken und unerwünschte Wirkungen	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Vorsicht	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Common specification(s) to comply with, if applicable:
N/A
Harmonised standards to apply, if applicable:
DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
DIN EN ISO 14971:2022-04 Medical devices – Application of risk management to medical devices
Relevant adopted monographs of the European Pharmacopoeia
Not applicable

9. Referenzen

- [1] Kühle M. Comparison of visual function with aspheric yellow, aspheric clear and spherical clear intraocular lenses. J Emmetropia 2013;4:123-130.
- [2] Schrecker J, Langenbucher A, Seitz B, Eppig T. First results with a new intraocular lens design for the individual correction of spherical aberration. J Cataract Refract Surg. 2018;44(10):1211-1219.
- [3] Rozegnal-Madej A, Zarnowski T. Biocompatibility of acrylic intraocular lenses in eyes with pseudoexfoliation syndrome - Aspira-aA (HumanOptics AG) vs AcrySof SA60AT (Alcon Lab) – one year follow up. 38th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; virtual Session: Pseudophakic IOLs: Monofocal). 2020.
- [4] Gyöngyössy B, Jirak P, Schönherr U. Rotational stability and patient satisfaction after implantation of a new toric IOL. Eur J Ophthalmol. 2016;26(4):321-327.
- [5] Gyöngyössy B, Jirak P, Schönherr U. Long-term rotational stability and visual outcomes of a single-piece hydrophilic acrylic toric IOL: a 1.5-year follow-up. Int J Ophthalmol. 2017;10(4):573-578.
- [6] Rusch W, Anwar M, Altas C, Steirat M, Fuchsluger T, Walckling M, et al. Influence Of Axial Length On Rotational Stability And Other Clinical Outcomes After Implantation Of A Hydrophilic Toric Intraocular Lens. 41th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; Session: IOLs and Evaluation Testing). 2023.
- [7] Schrecker J, Blass S, Langenbucher A. Silicone-diffractive versus acrylic-refractive supplementary IOLs: visual performance and manual handling. J Refract Surg. 2014;30(1):41-48.
- [8] Rajabi MT, Korouji S, Farjadnia M, Naderan M, Rajabi MB, Khosravi B, et al. Higher order aberration comparison between two aspherical intraocular lenses: MC6125AS and Akreos advanced optics. Int J Ophthalmol. 2015;8(3):565-568.
- [9] Liekfeld A, Ehmer A, Schröter U. Visual function and reading speed after bilateral implantation of 2 types of diffractive multifocal intraocular lenses: Add-on versus capsular bag design. J Cataract Refract Surg. 2015;41(10):2107-2114.
- [10] Chang JS, Ng JC, Chan VK, Law AK. Cataract Surgery with a New Fluidics Control Phacoemulsification System in Nanophthalmic Eyes. Case Rep Ophthalmol. 2016;7(3):218-226.
- [11] Lasta M, Miháltz K, Kovács I, Vécsei-Marlovits PV. Effect of Spherical Aberration on the Optical Quality after Implantation of Two Different Aspherical Intraocular Lenses. J Ophthalmol. 2017;2017:8039719.
- [12] Liu S-Y, Yang X. Visual outcome with the Oculentis Mplus intraocular lens. Int Eye Sci. 2017;17(11):2134-2136.
- [13] Naujokaitis T, Scharf D, Baur I, Khoramnia R, Auffarth GU. Bilateral implantation of +56 and +58 diopter custom-made intraocular lenses in patient with extreme nanophthalmos. Am J Ophthalmol Case Rep. 2020;20:100963.
- [14] Mayer CS, Hoffmann AM, Prahs P, Reznicek L, Khoramnia R. Functional outcomes after combined iris and intraocular lens implantation in various iris and lens defects. BMC Ophthalmol. 2020;20(1):370.
- [15] Sun Y, Zhang Z, Wei Y, Chi W, Zhang S. Intraocular Lens Fixation Technique Without Corneal Incision in Minimally Invasive Vitrectomized Eyes. Ophthalmol Ther. 2022;11(2):729-737.
- [16] Yang J, Hai-Bo C, Ming-Bing Z, Jun Y. Error analysis of visual quality in segmented refractive multifocal intraocular lenses implanted eyes. International Eye Science. 2022:447-451.
- [17] Mayer CS, Baur ID, Storr J, Khoramnia R. Traumatische Aniridie: Konservativer oder chirurgischer Therapieansatz? Ophthalmologe. 2022;119(4):395-399.

- [18] Chuprov AD, Startseva MI, Original, Чупров А, Старцева М. Comparative analysis of clinical and economic efficiency in implantation of various IOL models; Original Сравнительный анализ клинико-экономической эффективности при имплантации различных моделей ИОЛ. Journal "Medicine"; Original Журнал «Медицина» №. 2022;2:1.
- [19] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. J Cataract Refract Surg. 2022;48(1):75-82.
- [20] Liu X, Shao Y, Lin H, Liu C, Shen J, Zhang L, et al. Corneal densitometry: an innovative method to quantitatively evaluate corneal changes after phacovitrectomy. BMC Ophthalmol. 2023;23.
- [21] Kaindlstorfer C, Kneifl M, Reinelt P, Schonherr U. Rotation of a toric intraocular lens from neodymium:YAG laser posterior capsulotomy. J Cataract Refract Surg. 2018;44(4):510-511.
- [22] Shukhaiev SV, Ulitina AY, Shukhaeva EN. Comparative evaluation of rotational stability of hydrophobic (AcrySof IQ Toric) and hydrophilic (Torica-aAY) toric IOLs (Russian full-text; English abstract and translation available). Point of View East – West. 2019;1:16-19.
- [23] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia

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Остатъчни рискове и нежелани ефекти	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Внимание	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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- [1] Kühle M. Comparison of visual function with aspheric yellow, aspheric clear and spherical clear intraocular lenses. *J Emmetropia* 2013;4:123-130.
- [2] Schrecker J, Langenbucher A, Seitz B, Eppig T. First results with a new intraocular lens design for the individual correction of spherical aberration. *J Cataract Refract Surg.* 2018;44(10):1211-1219.
- [3] Rozegnal-Madej A, Zarnowski T. Biocompatibility of acrylic intraocular lenses in eyes with pseudoexfoliation syndrome - Aspira-aA (HumanOptics AG) vs AcrySof SA60AT (Alcon Lab) – one year follow up. 38th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; virtual Session: Pseudophakic IOLs: Monofocal). 2020.
- [4] Gyöngyössy B, Jirak P, Schönherr U. Rotational stability and patient satisfaction after implantation of a new toric IOL. *Eur J Ophthalmol.* 2016;26(4):321-327.
- [5] Gyöngyössy B, Jirak P, Schönherr U. Long-term rotational stability and visual outcomes of a single-piece hydrophilic acrylic toric IOL: a 1.5-year follow-up. *Int J Ophthalmol.* 2017;10(4):573-578.
- [6] Rusch W, Anwar M, Altas C, Steirat M, Fuchsluger T, Walckling M, et al. Influence Of Axial Length On Rotational Stability And Other Clinical Outcomes After Implantation Of A Hydrophilic Toric Intraocular Lens. 41th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; Session: IOLs and Evaluation Testing). 2023.
- [7] Schrecker J, Blass S, Langenbucher A. Silicone-diffractive versus acrylic-refractive supplementary IOLs: visual performance and manual handling. *J Refract Surg.* 2014;30(1):41-48.
- [8] Rajabi MT, Korouji S, Farjadnia M, Naderan M, Rajabi MB, Khosravi B, et al. Higher order aberration comparison between two aspherical intraocular lenses: MC6125AS and Akreos advanced optics. *Int J Ophthalmol.* 2015;8(3):565-568.
- [9] Liekfeld A, Ehmer A, Schröter U. Visual function and reading speed after bilateral implantation of 2 types of diffractive multifocal intraocular lenses: Add-on versus capsular bag design. *J Cataract Refract Surg.* 2015;41(10):2107-2114.
- [10] Chang JS, Ng JC, Chan VK, Law AK. Cataract Surgery with a New Fluidics Control Phacoemulsification System in Nanophthalmic Eyes. *Case Rep Ophthalmol.* 2016;7(3):218-226.
- [11] Lasta M, Miháltz K, Kovács I, Vécsei-Marlovits PV. Effect of Spherical Aberration on the Optical Quality after Implantation of Two Different Aspherical Intraocular Lenses. *J Ophthalmol.* 2017;2017:8039719.
- [12] Liu S-Y, Yang X. Visual outcome with the Oculentis Mplus intraocular lens. *Int Eye Sci.* 2017;17(11):2134-2136.
- [13] Naujokaitis T, Scharf D, Baur I, Khoramnia R, Auffarth GU. Bilateral implantation of +56 and +58 diopter custom-made intraocular lenses in patient with extreme nanophthalmos. *Am J Ophthalmol Case Rep.* 2020;20:100963.
- [14] Mayer CS, Hoffmann AM, Prahs P, Reznicek L, Khoramnia R. Functional outcomes after combined iris and intraocular lens implantation in various iris and lens defects. *BMC Ophthalmol.* 2020;20(1):370.
- [15] Sun Y, Zhang Z, Wei Y, Chi W, Zhang S. Intraocular Lens Fixation Technique Without Corneal Incision in Minimally Invasive Vitrectomized Eyes. *Ophthalmol Ther.* 2022;11(2):729-737.
- [16] Yang J, Hai-Bo C, Ming-Bing Z, Jun Y. Error analysis of visual quality in segmented refractive multifocal intraocular lenses implanted eyes. *International Eye Science.* 2022:447-451.
- [17] Mayer CS, Baur ID, Storr J, Khoramnia R. Traumatische Aniridie: Konservativer oder chirurgischer Therapieansatz? *Ophthalmologe.* 2022;119(4):395-399.

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Zbytková rizika a nežádoucí účinky	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Upozornění	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Υπολειπόμενοι κίνδυνοι και ανεπιθύμητες ενέργειες	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Προσοχή	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Riesgos residuales y efectos indeseables	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Precaución	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Risques résiduels et effets indésirables	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment Cystoid macular edema Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric <u>IOLs</u>: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Mise en garde	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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
4. Rizici i upozorenja

Preostali rizici i neželjeni učinci	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Oprez	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Maradék kockázatok és nemkívánatos hatások	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to <u>note</u> in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Figyelem	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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- [18] Chuprov AD, Startseva MI, Original, Чупров А, Старцева М. Comparative analysis of clinical and economic efficiency in implantation of various IOL models; Original Сравнительный анализ клинко-экономической эффективности при имплантации различных моделей ИОЛ. Journal "Medicine"; Original Журнал «Медицина» №. 2022;2:1.
- [19] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. J Cataract Refract Surg. 2022;48(1):75-82.
- [20] Liu X, Shao Y, Lin H, Liu C, Shen J, Zhang L, et al. Corneal densitometry: an innovative method to quantitatively evaluate corneal changes after phacovitrectomy. BMC Ophthalmol. 2023;23.
- [21] Kaindlstorfer C, Kneifl M, Reinelt P, Schonherr U. Rotation of a toric intraocular lens from neodymium:YAG laser posterior capsulotomy. J Cataract Refract Surg. 2018;44(4):510-511.
- [22] Shukhaiev SV, Ulitina AY, Shukhaeva EN. Comparative evaluation of rotational stability of hydrophobic (AcrySof IQ Toric) and hydrophilic (Torica-aAY) toric IOLs (Russian full-text; English abstract and translation available). Point of View East – West. 2019;1:16-19.
- [23] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia

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Rischi residui ed effetti indesiderati	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Attenzione	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Liekančios rizikos ir nepageidaujami reiškiniai	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection <u>with</u> toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Atsargiai	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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4. Risico's en waarschuwingen

Resterende risico's en ongewenste effecten	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Waarschuwing	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Ryzyko resztkowe i działania niepożądane	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Ostrzeżenie	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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4. Riscos e avisos

Riscos residuais e efeitos indesejáveis	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Precaução	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Riscuri reziduale și efecte nedorite	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Atenție	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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Kvarstående risker och oönskade effekter	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Varning	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation

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