Summary of safety and clinical performance (SSCP) for the Multifocal Triva-aA Family

SSCP reference number: P037

Revision history

Rev number	Revision date	Description of change	Revision validated by the Notified Body
Rev.01	2025-03-11	Initial creation	☒ Yes Validation language: English ☐ No (only applicable for class IIa or some IIb implantable devices (MDR, Article 52 (4) 2nd paragraph) for which the SSCP is not yet validated by the NB)

Table of contents

	Page
Part 1 – SSCP intended for healthcare professionals	3
1. Device identification and general information.....	3
2. Intended use of the device.....	4
3. Device description	4
4. Risks and warnings	7
5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)	9
5.1. Summary of clinical data related to equivalent device	9
5.2. Summary of clinical data from conducted investigations of the device before the CE-marking	9
5.3. Summary of clinical data from other sources.....	9
5.4. Overall summary of the clinical performance and safety.....	17
5.5. Ongoing or planned post-market clinical follow-up	25
6. Possible diagnostic or therapeutic alternatives	25
7. Suggested profile and training for users	26
8. Reference to any harmonised standards and CS applied	26
9. References	27
Part 2 – SSCP intended for patients	29
1. Device identification and general information.....	29
2. Intended use of the device.....	29
3. Device description	30
4. Risks and warnings	33
5. Summary of clinical evaluation and post-market clinical follow-up	35
5.1. Clinical background of the device.....	35
5.2. The clinical evidence for the CE-marking	35
5.3. Safety.....	36
6. Possible diagnostic or therapeutic alternatives	36
7. Suggested training for users.....	36

Part 1 – SSCP intended for healthcare professionals

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The SSCP is not intended to replace the Instructions For Use as the main document to ensure the safe use of the device, nor is it intended to provide diagnostic or therapeutic suggestions to intended users or patients. This SSCP is reviewed at least annually and will be updated if needed to ensure that any clinical and/or safety information in the SSCP remains correct and complete.

The following information is intended for users/healthcare professionals. Following this information there is a summary intended for patients.

1. Device identification and general information

Manufacturer's name and address	HumanOptics Holding AG Spardorfer Straße 150 91054 Erlangen Germany
Manufacturer's single registration number (SRN)	DE-MF-000017892
Authorised representative if applicable; name and SRN	N/A
Notified body name and organization number	TÜV NORD CERT GmbH NB 0044


Device trade name(s)	• Diff-aA (trade name: Diffractiva) • Diff-aAY (trade name: Diffractiva) • Triva-aA • Triva-aAY • TrivaT-aA • TrivaT-aAY (collectively referred to as "Multifocal Triva-aA Family")
Basic UDI-DI(s)	• Diff-aA/-aAY, Triva-aA/-aAY: 04049154_PC_M2_H2_O4_C2 • TrivaT-aA/-aAY: 04049154_PC_M2_H2_O5_C5
CND/EMDN code(s):	• Diff-aA/-aAY, Triva-aA/-aAY: P030102100202 (IOLs, APHAKIC, MULTIFOCAL, ASPHERICAL, HYDROPHILIC ACRYLIC) • TrivaT-aA/-aAY: P030102100302 (IOLs, APHAKIC, MULTIFOCAL, TORIC, HYDROPHILIC ACRYLIC)
Class of device	☐ I ☐ IIa ☒ IIb ☐ III
Year when the first certificate (CE) was issued covering the device	• Diff-aA/-aAY: 2010 • Triva-aA/-aAY: 2020 • TrivaT-aA/-aAY: 2021

2. Intended use of the device

Intended purpose	The devices are intended for implantation in the capsular bag to replace the natural lens.
Intended users	The devices must be handled by health professionals and implanted by ophthalmic surgeons.
Use environment	Ophthalmic operating room
Indications	The devices are indicated for the correction of aphakia after surgical removal of the natural lens. As multifocal IOLs, they are also indicated for patients who desire near, intermediate and distance vision with increased spectacle independence. However, accommodation will not be restored. Toric-multifocal IOLs can further correct preexisting regular, stable corneal astigmatism.
Medical condition(s) to be treated	Aphakia after surgical removal of the natural lens
Intended patient population	Aphakic adult patients Note: children are excluded due to the absence of clinical data and additional risks associated with paediatric cataract surgery. The manufacturer does not hold clinical data relating to pregnant/breastfeeding women or immunocompromised patient populations.
Intended anatomical location	Posterior chamber of the eye
Contraindications	Outside of general contraindications of ocular surgery, there are no specific contraindications for this device.

3. Device description

Short description of the devices	One-piece, aspheric aberration-free, multifocal diffractive posterior-chamber intraocular lenses (IOLs) with non-angulated C-loop haptics and a 360° lens epithelial cell barrier on the posterior optic surface. The toric versions are additionally fitted with a toric posterior optic surface (astigmatism correction).
Description of the devices	Diff-aA/-aAY, Triva-aA/-aAY and TrivaT-aA/-aAY are foldable, single-piece, posterior chamber IOLs with non-angulated C-loop haptics and an aspheric (aberration-free/neutral), multifocal optic. The anterior surface has a 3.5 mm apodised diffractive central zone consisting of 9 diffractive steps for the Diff-aA/-aAY and 7 steps for the Triva-aA/-aAY and the TrivaT-aA/-aAY. The peripheral refractive zone is comparable to a monofocal implant. The near addition power (ADD) for all lenses is +3.5 D, which corresponds to a distance of approximately 40 cm from the eye. Additionally, the Triva-aA/-aAY and TrivaT-aA/-aAY have an intermediate addition power (ADD) of +1.75 D corresponding to a distance of approximately 80 cm from the eye. The toric TrivaT-aA/-aAY has additionally markings in the periphery on the posterior surface, indicating the axis of the plus cylinder. The IOLs have an overall diameter of 12.5 mm and an optical diameter of 6.0 mm. The IOLs are made of hydrophilic acrylic copolymer with UV-absorber, a water content of 26% at 35°C and a refractive index of 1.46. The yellow variants (with "Y") additionally contain a blue light filter that absorbs the high-energy portion of the blue light between 400 and 500 nm.

	The posterior surface of the optic features a 360° lens epithelial cells (LEC) barrier that is aimed to counteract the migration of the LEC, thereby reducing the frequency of posterior capsule opacification (PCO). For the Diff-aA/-aAY and Triva-aA/-aAY the dioptr range covers +10.0 D to +30.0 D. For the TrivaT-aA/-aAY the dioptr range covers spherical equivalents between +10.0 D and +30.0 D and cylinders between +1.0 D and +6.0 D.  Schematic frontal and lateral view of the Diff-aA (left) and TrivaT-aA (right) with its toric markings. The IOLs exist also in yellow (Diff-aAY, Triva-aAY and TrivaT-aAY). The dimensions are given in mm.
Mode of action / operating principle	The devices are designed to focus light passing through the cornea and pupil onto the retina, like the eye's natural lens. The aspheric optic design aims to provide higher-contrast imaging compared to a spherical design. Multifocal IOLs are designed to provide improved vision at more than one distance, giving the patient a chance to be more independent from spectacles. While bifocal IOLs such as the Diff-aA/-aAY have two foci, a far and near focus, trifocal IOLs such as the Triva-aA/-aAY and TrivaT-aA/-aAY have an additional intermediate focus. Therefore, successful implantation of a trifocal IOL may result in full spectacle independence. The astigmatism-correcting function of the toric version TrivaT-aA/-aAY can be described as follows: corneal astigmatism means that the anterior corneal surface and/or the posterior corneal surface are not rotationally symmetrical. Therefore, such a cornea has a meridian of strongest refractive power and one of weakest refractive power, both of which are usually orientated perpendicular to each other (regular astigmatism). As a result of its toric posterior optic surface, the TrivaT-aA/-aAY also has meridians of strongest and weakest refractive powers. Those are arranged opposite to the corneal effect in order to correct the astigmatism.
Materials or substances in contact with the patient's tissues	The hydrophilic acrylic copolymer consisting of 2-hydroxyethyl methacrylate (HEMA) and methyl methacrylate (MMA) with a UV absorber and with a yellow coloring agent (only for models followed by "Y") is permanently in contact with the ocular tissues.
Method of sterilisation and allowable number of uses	The devices are sterile (steam sterilised) and designed for single use. Reprocessing of an implant is strictly prohibited.
Previous generations of the device(s)	The bifocal Diff-aA/-aAY is a combination of the 1-piece hydrophilic acrylic IOL Aspira-aA/-aAY (IOL platform) and the diffractive optic of the former 3-piece silicone IOL DIFF-s. The Triva-aA/-aAY has been developed based on the optical surface design of the bifocal Diff-aA/-aAY by adding an intermediate focus.

	The toric trifocal TrivaT-aA/-aAY combines the diffractive optic of the Triva-aA on the anterior surface with the toric optic of the Torica-aA/-aAY (and former bifocal ToricaDiff-aA/-aAY) on the posterior surface.
Variants and/or configurations covered by this document	- Optical design - Diff-aA/-aAY (non-toric, bifocal) - Triva-aA/-aAY (non-toric, trifocal) - TrivaT-aA/-aAY (toric, trifocal) - Colour: - Colourless: device names ending with “-aA” - Blue-filtering: device names ending with “-aAY” - Packaging types: - Compact line (CL) - Safeloader (SL) Note: The Safeloader system is a packaging type and not an accessory.
Accessories which are intended to be used in combination with the device	Not applicable Note: Injector systems are independent medical devices which are intended to be used in combination with the device in question (i.e. these devices are not accessories).
Other devices and products which are intended to be used in combination with the device	Implantation of the device can be performed using either an injector system (for Safeloader and Compact line) or forceps (for Compact line only). Tested injector systems for the Triva-aA/-aAY in the Compact line, depending on the power of the IOL, are the Medice AccuJect™ 1.8-1P (LP 604540) for a spherical equivalent (SE) of 10.0 to 25.0 D, the Medice AccuJect™ 2.0-1P (LP 604510) for a SE of 10.0 to 30.0 D and the Medice AccuJect™ 2.2-1P (LP 604530) for a SE of 10.0 to 30.0 D (Medice AG, Switzerland). Tested injector systems for the Triva-aA/-aAY and the TrivaT-aA/-aAY in the Safeloader, depending on the power of the IOL, are the Medice AccuJect™ 2.0-1P (LP 604510, Rev. C (2v7)) for a SE of 10.0 to 30.0 D and a cylinder of 1.0 to 6.0 D and the Medice AccuJect™ 2.2-1P (LP 604530, Rev. C (2v7)) for a SE of 10.0 to 30.0 D and a cylinder of 1.0 to 6.0 D. In addition, the Medice AccuJect™ 1.8-1P (LP 604540, Rev. C (2v7)) is tested for the Triva-aA/-aAY in the Safeloader for a SE of 10.0 to 25.0 D. Tested injector system for the Diff-aA/-aAY in the Safeloader is the Medice AccuJect™ 2.0-1P (LP 604510 Rev. C (2v7)). Note: Injectors designed for use with the HumanOptics Safeloader products feature a yellow plunger. When using a non-tested injector system for implantation of HumanOptics IOLs packaged in the Compact line, the user must refer to the specific instructions for use provided with the injector system to ensure that it is suitable to use with the IOL to be implanted.

4. Risks and warnings

Residual risks and undesirable effects	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>multifocal</u> IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Caution	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection

	• Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal</u> IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Warnings	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Special considerations prior to multifocal IOL implantation	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.
Special considerations prior to toric IOL implantation	• The dioptric power of toric IOLs is indicated as “spherical equivalent (SE) and cylinder (cyl)”. Please check the label carefully. • Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism. • Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation • If possible, posterior corneal astigmatism (tomography) should be considered, especially with astigmatism < 2.0 D.

Other relevant aspects of safety including a summary of any field safety corrective action (FSCA including FSN), if applicable	There have been no recalls since the devices have been launched. In July 2024, HumanOptics issued a Field Safety Notice affecting all devices in a Safeloader packaging manufactured in the period January 1, 2020 to June 28, 2024. During the standardized incoming goods inspection of returns, HumanOptics had observed very rare cases of damaged outer blister packaging after transportation. The FSN informed the users about the observation and reminded them to check the integrity of the sterile barrier system, as requested in the Instructions for Use, along with a visual instruction how to identify potential cracks. This damage is visible in the form of cracking and can be identified by users before use. It could theoretically cause a contamination of the surgical equipment or the intraocular lens. However, the hazard potential due to crack is considered to be very low. Statistical valuation showed that probability of occurrence of the additional risk of endophthalmitis due to the cracks is by several orders of magnitude lower than the general rate of endophthalmitis after cataract surgery. No patient injury has happened to the manufacturer's best knowledge. Meanwhile, the product packaging has been improved to eliminate this issue completely.

5. Summary of clinical evaluation and post-market clinical follow-up (PMCF)

5.1. Summary of clinical data related to equivalent device

Not applicable. The conformity of the device was not assessed on the basis of equivalence.

5.2. Summary of clinical data from conducted investigations of the device before the CE-marking

Not applicable. All clinical investigations were performed after the first CE marking under the European Medical Device Directive 93/42/EEC (see section 5.3).

5.3. Summary of clinical data from other sources

Clinical data confirming the safety and performance of the Multifocal Triva-aA Family have been derived from the following sources:

- Several post-market studies for the devices under consideration and related/parent models
- A comprehensive literature review using PubMed and other databases

Literature search strategies were designed to identify articles relevant to the devices. Both favorable and unfavorable references were identified and summarized.

5.3.1. Post-market studies

The manufacturer HumanOptics has initiated several post-market studies with the Diff-aA/-aAY, Triva-aA/-aAY and TrivaT-aA/-aAY as shown in Table 1. The studies/activities were labeled with acronyms considering: [device in focus]_[city]_[primary end point].

In the two post-market studies conducted with the Diff-aA, a total of 55 eyes of 32 patients were implanted and followed for 1 year. Both studies showed good and stable visual performance and contrast sensitivity with high patient satisfaction. The defocus curve showed the typical two-peaked curve which is expected for a bifocal IOL. Overall, these studies confirmed the long-term safety and clinical performance of the device.

In the case of the Triva-aA/-aAY, final or interim clinical results available to date from 105 eyes of 53 patients were collected in real-world clinical settings and in an ongoing post-market study. The results show that the Triva-aA/-aAY exhibits good visual performance across all ranges – near, intermediate and far distances. Moreover, the enhancement of the intermediate vision was not at the expense of the far and near vision compared to the bifocal Diff-aA. The current data support that the implantation of the Triva-aA/-aAY is a safe and effective surgical treatment for cataract patients.

As for the TrivaT-aA/-aAY, the toric version of the Triva-aA/-aAY, early results from the controlled market release phase (CRMP) of 30 eyes are positive and show that a good visual performance at all three distances can be achieved as well as a high rotational stability of the IOL in the capsular bag. The post-market study in the context of this CMRP is ongoing and data analysis will be completed end of 2024. This new set of data will certainly provide further evidence on the safety and performance of the Multifocal Triva-aA Family, including astigmatism correction and stability in the capsular bag. Further evidence of the rotational stability of the aA IOL platform is provided by three studies with the toric parent models (one ongoing), which currently involve 112 eyes of 70 patients. One of the studies was conducted with the bifocal ToricaDiff-aA and two with the monofocal Torica-aA. In addition, these studies show that the toric optic of HumanOptics IOLs successfully corrects corneal astigmatism. One of the Torica-aA studies also confirmed the good long-term (1.5 years) performance of the lens in terms of visual acuity, astigmatism correction, safety and patient satisfaction. Clinical data for the Torica-aA (Basic UDI-DI 04049154_PC_M2_H2_O3_BX) can be found in the SSCP for the Monofocal Aspira-aA Family [1]. A summary of clinical data for the ToricaDiff-aA is available on request.

Details of the study results with the devices under consideration are provided in the subsequent paragraphs.

Table 1: Overview of clinical studies/activities with the Multifocal Triva-aA Family

Name of activity / Type of activity / Design	Aim	Status / Timelines
Diff-aA/-aAY		
Dif-aA_Gla_Comb Post-market study RCT	Comparative study to confirm the long-term (1-year) performance and safety of the Diff-aA vs. the Diff-sPB after implantation in cataract patients. To assess the patient-related outcome based on patient satisfaction. The add-on lens was implanted in a combined surgery together with the monofocal lens (Aspira-aA)	Completed FU (1 year) completed: 2012
Dif-aA_Gla_Seco Post-market study PNCS	Comparative study to confirm the long-term (1-year) performance and safety of the Diff-aA vs. the Diff-sPB after implantation in cataract patients. To assess the patient-related outcome based on patient satisfaction. Unlike in the previous study, the add-on lens was implanted in the sulcus during a second surgery in pseudophakic eyes.	Completed FU (1 year) completed: 2015
Triva-aA/-aAY		
Tri-aA_CMRP Prospective collection of RWD, CMRP	To assess the clinical performance of the new Triva-aA by 3 selected users, all with an extensive experience with the former models Diff-aA and/or ToricaDiff-aA.	Completed FU (up to 1 year) completed: Q2/2021
Tri-aAY_Pra_RWD Prospective collection of RWD	To assess the clinical performance of the new Triva-aAY by a surgeon with an extensive experience with several models of trifocal IOLs.	Completed FU (4-6 months) completed: Q3/2022
Tri-aAY_Gla_VA Post-market study PSAS	To confirm the results (visual performance and patient-related outcome) obtained by one of the CMRP investigators (above) in the course of a clinical study.	Running Implantation in first patient: 05/2022 Completion of the FU (1 year) planned for Q4/2024 Completed: 1 M Interim Clinical Report
TrivaT-aA/-aAY		
TriT-aAY_Ali_Refr Post-market study IIS, CMRP PSAS	To confirm the performance (incl. residual astigmatism and visual outcomes) and safety (including rotational stability) of the TrivaT-aAY after bilateral implantation in patients with pre-existing corneal astigmatism. In particular, the study focuses on postoperative refractive outcomes. To assess the patient-related outcome based on patient satisfaction. This PSAS was initiated in the course of a CMRP.	Running Implantation in first patient: 05/2023 Completion of FU (4 months) planned for Q4/2024 Some early data with 1M FU available (Interim Clinical Report)

FU=Follow-up; CMRP=controlled market release phase; RWD= real-world data; RCT=randomized controlled trial; PNCS=prospective non-randomized comparative study; IIS= investigator-initiated study; PSAS=prospective single-arm study; M= month

Post-market study Dif-aA_Gla_Comb (Diff-aA vs. Diff-sPB – Glauchau, Germany [2])

The purpose of this prospective post-market study was to compare the long-term (1-year) functional outcomes and patient satisfaction after implantation of either the bifocal hydrophilic acrylic Diff-aA (HumanOptics) or the “add-on group”, where the three-piece bifocal silicone Add-on IOL Diff-sPB was implanted together with the monofocal IOL. The main difference regarding the design of this study (Dif-aA_Gla_Comb) compared to the post-market study Dif-aA_Gla_Seco was that in the “add-on group” the multifocal add-on lens Diff-sPB was implanted together with the monofocal IOL, whereas in the post-market study Dif-aA_Gla_Seco the implantation of the latter had already occurred earlier. 25 eyes (15 patients) were implanted with the Diff-aA and 29 eyes (19 patients) with the Diff-sPB and followed up for 1 year after surgery.

In the Diff-aA group, 88% of eyes from were within ± 0.50 D of emmetropia one year postoperatively. All eyes achieved uncorrected monocular visual acuity (VA) of 0.1 logMAR or better for far vision and 0.3 logMAR or better for intermediate (1 m) and near distance vision. As a result of its bifocal design, the defocus curve showed two peaks at far and near distance foci (0.00 and -2.50 D, respectively) with visual acuities of 0.00 and 0.10 logMAR, respectively. Contrast sensitivity testing showed good outcomes under photopic lighting conditions with values being within or slightly above the normal age-related ranges.

Due to a relatively low number of bilateral cases (10 per group), the patient questionnaire can only provide a limited informative value. As expected with all MIOLs, the majority (70%) of patients experienced disturbing photic phenomena (halos, glare). The degree of impairment was assessed as moderate by one patient and minimal by six patients. 90% or more of the patients were “very satisfied” or “satisfied” with their far and near distance vision without correction at daylight, with the percentages being slightly lower for indoor and night conditions (which is expected for MIOLs). All patients in each group stated to be independent from reading glasses (which supports the statement regarding spectacle-independence made in the indications for use) and all stated that they would choose this lens again.

There were no intraoperative complications. Postoperatively, there was one case of macular oedema (Irvine-Gass syndrome).

The results for the Diff-sPB group were comparable to the Diff-aA. However, minor differences were observed with regards to visual acuity at different distances, the defocus curves, the contrast sensitivity and patient satisfaction. In terms of complications, two cases with mild pigment dispersion were observed (see study Dif-aA_Gla_Seco for explanation).

In conclusion, the results of this clinical investigation confirm the good visual performance and contrast sensitivity for the Diff-aA, with minimal differences were observed between the two tested IOLs. Those might be explained by the fact that the two lens types differ in chemical material (acrylic vs. silicone) and in the optical design (aberration-free vs. aberration-correcting). Taken together, these long-term results indicate a high level of safety and a good performance for Diff-aA combined with a high patient satisfaction.

Post-market study Dif-aA_Gla_Seco (Diff-aA vs. Diff-sPB – Glauchau, Germany) [3]

The purpose of this prospective study was to compare the long-term (1-year) functional outcomes and patient satisfaction after implantation of either the bifocal hydrophilic acrylic Diff-aA (HumanOptics) or a three-piece bifocal silicone Add-on IOL (Diff-sPB, HumanOptics). Both IOLs have the same diffractive optic designs but the add-on is aberration-correcting while the Diff-aA is aberration-free. The add-on lens was implanted in the sulcus during a second

surgery (primary surgery was performed previously with a monofocal IOL in the capsular bag) and the Diff-aA was implanted in the capsular bag. In each group 30 eyes were implanted, corresponding to 17 and 21 patients in the case of the Diff-aA and the Diff-sPB, respectively.

In the Diff-aA group, 97% of eyes from were within ± 0.50 D of emmetropia one year postoperatively. All eyes achieved uncorrected monocular visual acuity (VA) of 0.1 logMAR or better for far vision and 0.3 logMAR or better for intermediate (1 m) and near distance vision. As a result of its bifocal design, the defocus curve showed two peaks at far and near distance foci (0.00 and -2.50 D, respectively) with visual acuities of 0.00 and 0.10 logMAR, respectively. Contrast sensitivity testing showed good outcomes under photopic lighting conditions with values being within or slightly above the normal age-related ranges.

As expected with all MIOLs, the majority of patients experienced halos, and, to a lower extent, also glare but more than one third of the patients experienced no photic disturbances and all patients rated the visual disturbances as minimal or mild. All patients were “very satisfied” (92%) or “satisfied” (8%) with their far and near distance vision without correction at daylight, with the percentages being slightly lower for indoor and night conditions (which is expected for MIOLs). All patients stated that they would choose this lens again.

There were no intraoperative complications. Postoperatively, mild pigment dispersion without clinical consequences was observed in 1 eye from the Diff-aA group. At 1 year, Nd:YAG laser capsulotomy was performed in 7 eyes in the Add-on group and in 5 eyes in the Diff-aA group

In the Diff-sPB group comparable results were obtained but minor differences compared to the Diff-aA were observed in terms of the defocus curves, the contrast sensitivity and patient satisfaction. In terms of complications, the number of cases with mild pigment dispersion was substantially higher (8 cases) in the former group. This was possibly caused by the short distance between the add-on lens and the iris.

In conclusion, the Diff-aA IOL showed good and stable performance in terms of visual acuity and contrast sensitivity over a 1 year follow-up period. The deviation from target refraction was very small for both groups. Interestingly, intermediate distance visual acuity was better than expected for a bifocal lens but the defocus curves showed the typical two-peaked curve shape. These long-term results show the safety and efficacy of the procedures with all patients being highly satisfied with the outcomes.

Real-world data Tri-aA_CMRP (Triva-aA – Controlled market release phase)

Before the official market introduction of the Triva-aA, the IOL was evaluated by 3 ophthalmic surgeons, all with an extensive experience with the bifocal HumanOptics IOLs. Due to the positive experience from this phase, the device was subsequently made available to a broader population of users. In total, more than 50 eyes were implanted with the Triva-aA and followed for a period of at least 6 months.

Two of the sites (Paris and Madrid) just shared their positive impressions with the IOL. The data from Glauchau, which are presented below, were analysed in detail and confirmed the positive feedback from the two first testing sites.

Eight patients were bilaterally implanted with the Triva-aA. Six months after surgery 31% (n = 5) of the eyes were within ± 0.50 D and 88% (n= 14) were within ± 1.00 D of the intended emmetropic target refraction. At this clinical site, monocular and binocular corrected visual acuity (VA) for far, intermediate and near distance were good with mean values ranging between 0.06 logMAR and -0.04 logMAR at the 6-month follow up visit. Without best distance-correction values were not significantly different whatever the distance. All patients achieved uncorrected binocular VA of 0.10 logMAR or better for far and intermediate, and 0.20 logMAR or better for near distance. Between the 6-week and the 6-month follow-up monocular and binocular VA (corrected and

uncorrected) did not significantly change except for monocular UDVA ($p=0.03$) and UNVA ($p=0.04$), which both improved over time.

The aforementioned outcomes are consistent with those of the monocular defocus curve, which is the best objective measure of expected vision at different distances. Indirect comparison with the monocular defocus curve of the bifocal Diff-aA also showed an enhanced intermediate vision with the Triva-aA, while near and far vision were comparable for both lenses, confirming the trifocal design of the lens. When measured binocularly, showed a wide plateau with values better than 0.05 logMAR from -3.0 D to +0.5 D. This indicates good visual performance over various distances.

Decentration greater than 0.4 mm and tilt greater than 7.0 degrees is thought to lower the visual performance of aberration-correcting aspheric IOLs [4], particularly with MIOLs [5]. OCT-based analysis in Glauchau showed, showed that none of Triva-aA IOLs exceeded those values. Mean values for tilt and decentration were $\sim 5.5^\circ$ and ~ 0.15 mm, postoperatively for the IOL as well as preoperatively for the crystalline lens, which is in the range previously described for other IOLs [6]. There was no significant change between the 6-week and the 6-month visit.

All patients were very satisfied or satisfied with their vision and were free from spectacles for far and intermediate vision. Three patients reported to sometimes need glasses for near vision, but only for very small letters. In terms of photic phenomena, one of the patients with large pupils was disturbed by halo and glare. Another patient reported to be sensitive to glare when driving during the night.

No IOL-related complications or serious adverse events were observed. Two eyes of the same patients developed pseudophakic cystoid macular edema (CME), also known as Irvine-Gass syndrome, which resolved upon treatment. This syndrome is considered as a normal side effect of cataract surgery.

In conclusion, the clinical data from Glauchau show that the Triva-aA IOL is safe, stable and effective in providing a comparable visual performance across all ranges – near, intermediate and far distance. The enhancement of the intermediate vision with the trifocal design compared to the bifocal design (Diff-aA) was not at the expense of the far and near vision as shown by the defocus curves of both models. The current data, although from a small series of cases, are valuable because they represent “real-world” data with no pre-established study protocol

Real-world data Tri-aAY_Pra_RWD (Triva-aAY – Prague/Zlín, Czech Republic)

The aim of this prospective, multicentric data collection was to provide an end-user expert opinion on the use and usability of the trifocal Triva-aAY IOL in a real-world setting. Patients with bilateral cataract who had been implanted with the IOL and had attended the routinely scheduled follow-ups at 1-week and 4 to 6 months postoperatively were included in the analysis. Overall, 16 patients (32 eyes) with a median age of 58 years were included in the study.

Measurement of refractive accuracy showed good spherical equivalent (SE) predictability, with 75% of eyes within ± 0.50 D and 94% within ± 1.0 D six months after surgery. The mean postoperative SE was 0.24 ± 0.48 D.

The monocular uncorrected and distance-corrected visual acuity (VA) revealed a good performance of the Triva-aAY at far, intermediate (66 cm) and near (40 cm) distance. At 4-6 months postoperatively, the uncorrected distance VA was 0.3 logMAR or better in 94% of the eyes for distance (UDVA), in 97% of the eyes for intermediate (UIVA) and in 88% of the eyes for near vision (UNVA). With distance-correction, 97% of the eyes achieved 0.3 logMAR or better for distance (CDVA) and intermediate vision (DCIVA) and 100% of the eyes achieved 0.3 logMAR or better for near vision (DCNVA). Binocular performance was even better with all eyes reaching 0.3 logMAR or better for uncorrected and distance-corrected far, intermediate and near vision, except for two eyes (4%) for uncorrected near vision (UNVA). As a result, all patients were spectacle independent at the 4-6-months follow-up.

The defocus curve is the best objective measure of expected vision at different distances. At 4 to 6 months postoperatively, both the monocular and the binocular curve showed the typical progression for a trifocal lens with two small peaks at the expected distance and near focus. A slight drop down was seen at the expected intermediate focus (-1.0 D vergence). Overall, between +1.0 D and -3.0 D of vergence, the Triva-aAY showed a wide range of useful vision, with VA higher than 0.2 logMAR. Binocular visual acuity reached its maximum at 0.0 D of defocus (distance vision) with 0.01 ± 0.05 logMAR.

In terms of photic phenomena, the majority of the patients reported to never or only occasionally experience halos (54%), starburst (64%) or glare (91%). The severity was rated as moderate to severe by 47% of patients for halos, by 37% for starburst and by 10% for glare. However, no patient considered the symptoms as “very bothersome”. Only 23% and 18% of the patients felt “quite bothered” by halos and starburst, respectively, and none of the patients in the case of glare.

No intraoperative or postoperative complications were observed.

In conclusion, this analysis shows that the Triva-aAY provides a good visual performance and a high level of safety in a real-world setting. With an extended plateau from +1.0 D to -3.0 D of defocus in which a visual acuity of 0.2 logMAR or better is maintained, the Triva-aAY is effective in providing a comparable visual performance at all distances. Thus, the data support that the implantation of the Triva-aAY is a safe and effective surgical treatment for cataract patients who also desire spectacle independence.

Post-market study Tri-aAY_Gla_VA (Triva-aAY – Glauchau, Germany)

This prospective, monocentric post-market study was initiated after the same investigator had evaluated the Triva-aA (clear variant) in the context of a controlled market release phase. Thus, the IOL was further assessed by an ophthalmic surgeon who had experience with the Triva-aAY and MIOLs in general. At the time of writing this evaluation report, the interim analysis of the 1 week and 1-month follow-up was available, with 57 eyes of 29 patients implanted with the Triva-aAY. The study will be continued for a long-term evaluation over the course of 1 year. Overall, the purpose of the study is to collect more clinical data for the evaluation of the trifocal Triva-aAY IOL in terms of long-term safety and visual outcomes after bilateral implantation in cataract patients.

One month after surgery, the measurement of refractive accuracy showed a good predictability of the spherical equivalent (SE) with 76.6% of eyes being within ± 0.50 D and 97.9% within ± 1.0 D. The mean postoperative SE was -0.27 ± 0.31 D (range: -1.08 to 0.55 D).

The monocular uncorrected and distance-corrected visual acuity (VA) revealed a good performance of the Triva-aAY at far, intermediate (60 cm) and near (40 cm) distance. At 1 month postoperatively, 97.9% of the eyes achieved an uncorrected monocular VA of 0.3 logMAR or better for far and near distances as well as 93.8% for intermediate distance. A value of 0.0 logMAR or better was reached by 39.6% of eyes for far, 29.2% for intermediate and 43.8% for near distances. With distance-correction, the percentages of patients with a monocular VA of 0.0 logMAR or better increased to 75.0%, 35.4% and 58.3% for far, intermediate and near distance, respectively. The postoperative mean values for monocular UDVA and CDVA were 0.09 ± 0.10 and 0.03 ± 0.06 logMAR, respectively.

The defocus curve is the best objective measure of expected vision at different distances. At 1 month postoperatively, the monocular curve displayed a shape with two peaks at the expected far and near foci and a slight slope in between at the expected intermediate focus. Maximum mean VA was reached at 0.0 D of defocus (far vision) with 0.03 ± 0.06 logMAR. At -2.5 D of defocus (near vision at approximately 40 cm from the eye), the mean VA was 0.07 ± 0.10 logMAR. In the intermediate range, the Triva-aAY showed a continuous performance with lowest values at -1.0 D and -1.5 D (approximately 67 cm from the eye) of defocus, with 0.14 ± 0.10 logMAR and

0.12±0.09 logMAR, respectively. Overall, between +1.00 D and -3.50 D of vergence, the curve showed a wide range of useful vision, with VA values higher than 0.25 logMAR, which indicates a good visual performance over various distances.

Patient satisfaction and subjective reported spectacle-independence will be available upon the 6 month follow-up.

There were no serious adverse events reported during the current course of the study.

In conclusion, the currently available interim data from the Tri-aAY_Gla_VA study show that the Triva-aAY exhibits good visual performance in terms of visual and refractive measurements together with safety. Moreover, the IOL is effective in providing a comparable visual performance across all ranges – near, intermediate and far distance. The current data support that the implantation of the Triva-aAY is a safe and effective surgical treatment for cataract patients.

Post-market study TriT-aAY_Ali_Refr (TrivaT-aAY – controlled market release phase, Valencia and Alicante, Spain)

The goal of this investigator-initiated post-market study is to evaluate the toric TrivaT-aAY in the context of a controlled market release phase (CMRP). Ophthalmic surgeons from 2 centers in Spain (Alicante and Valencia) and with an extensive experience in implanting trifocal IOLs were selected to conduct the study.

The study protocol includes 50 eyes of 25 patients with senile cataract and astigmatism. The standard clinical outcomes for toric MIOLs will be evaluated up to 4 months postoperatively.

Currently, interim results of 30 eyes are available for the 1 month follow-up. Measurement of refractive accuracy showed good spherical equivalent (SE) predictability with mean values of 0.08 ± 0.26 D (ranging from -0.38 to 0.75 D). The mean residual refractive astigmatism was -0.28 ± 0.35 D (ranging from 0.00 to -1.00 D).

The monocular uncorrected and distance-corrected visual acuity (VA) revealed a good performance of the TrivaT-aAY at far, intermediate (60 cm) and near (40 cm) distances 1 month postoperatively. The mean uncorrected distance VA (UDVA) was 0.06 ± 0.07 logMAR; with correction the mean distance VA (CDVA) increased to 0.02 ± 0.04 logMAR. The uncorrected intermediate VA (UIVA) had a mean of 0.13 ± 0.08 logMAR and increased to 0.10 ± 0.08 logMAR with distance correction (DCIVA). The mean uncorrected near VA (UNVA) was 0.14 ± 0.11 logMAR and the mean distance-corrected near VA (DCNVA) was 0.10 ± 0.10 logMAR. Overall, 1 month postoperatively, all eyes achieved good visual acuity at all three distances, which is expected to further improve by the 4-month follow-up.

The rotational stability was assessed with subjective analysis in 15 eyes and revealed a mean absolute IOL rotation of $0.93^\circ \pm 1.28^\circ$, ranging from 0 to 4° , at 1 month postoperatively.

In conclusion, early data with the TrivaT-aAY show that the lens provides a good visual performance with respect to uncorrected as well as distance-corrected VA at far, intermediate and near distances at the 1-month follow-up. Moreover, the rotational stability of the TrivaT-aAY is excellent confirming the suitability of the IOL platform for astigmatism correction.

5.3.2. Literature review

The 7 articles identified through a systematic literature review are of various levels of usefulness to assess the safety and performance of the devices of the Multifocal Triva-aA Family [7-22].

Five eligible publications report the outcomes after implantation with the Diff-aA in a total of 442 eyes [23-27]. Additionally, two publications relating to 36 implanted eyes were selected for the toric parent model Torica-aA [21, 22]. In the case of the Triva-aA/-aAY no eligible articles have been found yet. The TrivaT-aA/-aAY has only recently been launched and is so far only available to selected users in the context of a controlled market release phase.

All articles about the Diff-aA concordantly reported a good visual acuity at near and far distance and (where measured) a high degree of spectacle independence. Glare and halos were observed by several patients but did not impact on the overall patient satisfaction, which was very high. Postoperative complications were not observed¹ and only in one case an intraoperative complication (clinically insignificant decentration due to haptics breakage) was reported (assessed in ref. [23] and [27]). In two of the studies, the Diff-aA was compared to competitor diffractive bifocal IOLs (hydrophobic acrylic) in randomized controlled trials. One [23] showed that the visual performance and safety were similar to the Tecnis ZLB00 IOL. This spherical aberration-correcting IOL was more sensitive to a large angle alpha than the spherical aberration-free Diff-aA, which in turn displayed a stronger image degradation effect of pupil enlargement. The other study [25] showed that the Diff-aA is at least on par with the well-established hydrophobic AcrySof ReSTOR +3.0 MIOL in terms of visual acuity and reading performance and provided functional reading performance even under low illumination and contrast conditions. Finally, one retrospective data analysis [26] found that the rate of clinically significant PCO (Nd:YAG laser capsulotomies) was inferior with the Diff-aA (1.1%) compared to the average of 8 tested multifocal IOLs (17.7%). Overall, these clinical data demonstrate that implantation of the Diff-aA is a safe and effective treatment for cataract patients.

Regarding the toric parent model for the TrivaT-aA/-aAY, the Torica-aA/-aAY, one article [22] reported good rotational stability of the device for 35 eyes. One additional case report [21] confirmed that the device allows to correct astigmatism and provides a good visual acuity.

5.4. Overall summary of the clinical performance and safety

MIOLs are intended to provide vision at multiple distances. Toric MIOLs are further intended to compensate for pre-existing corneal astigmatism. Accordingly, the main criterion for the success of cataract surgery is the long-term visual result at all distances. Therefore, the main performance indicators for the Multifocal Triva-aA Family are the visual acuity at far, (intermediate) and near distance, and the reduction of pre-existing astigmatism for the toric models. Key safety indicators are the rate of complications and, in the case of the toric devices, the rate of IOL rotations.

PERFORMANCE ANALYSIS

Based on the results of the post-market studies and real-world data collections with the Multifocal Triva-aA Family and their multifocal optic parents (chapter 5.3.1), a performance analysis was performed in accordance with the ISO 11979-7:2018-08. A total of 202 eyes of cataract patients were enrolled and implanted with the Diff-aA/-aAY (n=55), the Triva-aA/-aAY (n=105), the TrivaT-aA/-aAY (n=20) or the multifocal toric parent model ToricaDiff-aA/-aAY (n=22). Our analysis included all available results from a period ranging from 4-6 weeks to one year post-operatively. The primary performance endpoint, i.e. achieving monocular CDVA of 0.3 logMAR or better, was met by all eyes, except one at the 3-month follow-up (201/202; 99.5%)

¹ In the presentation by Vo Tan, 2015 [26] some postoperative complications were mentioned but could not be analyzed individually for the Diff-aA.

and by all eyes at the 1-year follow-up (55/55), exceeding the primary safety and performance endpoints (SPE) rates as defined in ISO 11979-7:2018-08 for posterior-chamber lenses (92.5% overall). Similar findings were reported in the literature [24, 27]. **Key performance CDVA outcomes are presented in Table 2.** The results indicate that the visual performance of the Multifocal Triva-aA Family is good and remains stable over the investigated study period of up to one year.

Table 2: Evaluation of performance CDVA outcomes according to DIN EN ISO 11979-7:2018-08

Endpoint	Acceptance criterion	Clinical outcome	Data source
Monocular CDVA (% of eyes reaching 0.3 logMAR or better)	≥ 92.5% of the eyes should achieve CDVA 0.3 logMAR or better [SPE rate]		<u>Post-market studies and RWD collection:</u>
		Diff-aA/-aAY	
		100% [0.00 to 0.10] (n=30) at 4-6-weeks 100% [0.00 to 0.20] (n=30) at 3-month 100% [0.00 to 0.10] (n=30) at 1-year	Diff-aA_Gla_Seco
		100% [0.00 to 0.10] (n=25) at 1-year	Diff-aA_Gla_Comb
		Triva-aA/-aAY	
		100% [0.00 to 0.10] (n=16) at 6-weeks 100% [0.00] (n=16) at 6-month	Tri-aA_CM RP (Glauchau)
		97% [IQR 0.00 to 0.01] (n=32) at 4-6-month	Tri-aA_Pra_RWD
		100% [0.00 to 0.22] (n=57) at 1-month	Tri-aA_Gla_VA*
		TrivaT-aA/-aAY	
		100% [0.00 to 0.16] (n=20) at 1-month	TriT-aAY_Ali_Refr*
		ToricaDiff-aA/-aAY**	
		100% [-0.10 to 0.15] (n=22) at 1-month 100% [-0.10 to 0.15] (n=20) at 3-month	TDI-aA_Hal_VA

* Here, interim data from the 1-month follow-up were used as a longer follow-up period was not available yet

** Related former model (bifocal toric), not available anymore on the market

As MIOs are also designed to also improve near and (for trifocal IOLs) intermediate distance vision, an additional analysis was performed for the **performance endpoint "vision at multiple distances"**. **Key performance outcomes are presented in Table 3** and show that the

acceptance criteria were met in each study. Uncorrected VA values were only minimally worse than the corrected ones for all distances.

Table 3: Evaluation of performance DCIVA and DCNVA outcomes

Endpoint	Acceptance criterion	Clinical outcome	Data source
DCNVA corresponding to the near ADD of the IOL DCIVA corresponding to the intermediate ADD of the IOL	For <u>bifocal</u> and <u>trifocal</u> IOLs: Mean** DCNVA ≤ 0.3 logMAR		<u>Post-market studies and RWD collection:</u>
		Diff-aAI/-aAY	
	For <u>trifocal</u> IOLs: Mean DCIVA ≤ 0.3 logMAR	DCNVA Mean $\pm$ SD: 0.07 ± 0.10 (n=30) at 1-year DCIVA Mean $\pm$ SD: 0.09 ± 0.10 (n=30) at 1-year	Diff-aA_Gla_Seco
		DCNVA Mean $\pm$ SD: 0.08 ± 0.08 (n=25) at 1-year DCIVA Mean $\pm$ SD: 0.09 ± 0.10 (n=25) at 1-year	Diff-aA_Gla_Comb
		Triva-aAI/-aAY	
		DCNVA Mean $\pm$ SD: 0.06 ± 0.06 (n=16) at 6-month DCIVA Mean $\pm$ SD: 0.04 ± 0.06 (n=16) at 6-month	Tri-aA_CM RP (Glauchau)
		DCNVA Median 0.12 [IQR 0.03 to 0.16] (n=32) at 4-6-month DCIVA Median 0.02 [IQR 0.0 to 0.12] (n=32) at 4-6-month	Tri-aA_Pra_RWD
		DCNVA Mean $\pm$ SD: 0.04 ± 0.13 (n=48) at 1-month DCIVA Mean $\pm$ SD: 0.10 ± 0.13 (n=48) at 1-month	Tri-aA_Gla_VA*
		TrivaT-aAI/-aAY	

		DCNVA Mean $\pm$ SD: 0.10 ± 0.10 (n=20) at 1-month DCIVA Mean $\pm$ SD: 0.10 ± 0.08 (n=20) at 1-month	TriT-aAY_Ali_Refr*
		ToricaDiff-aA/-aAY DCNVA Median 0.0 [IQR -0.1 to 0.3] (n=20) at 3-month Median 0.0 [IQR -0.1 to 0.2] (n=20) at 6-month DCIVA Median 0.1 [IQR -0.1 to 0.3] (n=20) at 3-month Median 0.1 [IQR -0.1 to 0.2] (n=20) at 6-month	TDI-aA_Hal_VA

* Here, interim data from the 1-month follow-up were used as a longer follow-up period was not available yet

** DCNVA/DCIVA: median was used if mean was not available

Regarding the performance endpoint “reduction of pre-existing astigmatism”, the results from the TrivaT-aA/-aAY and the toric optic parents Torica-aA and ToricaDiff-aA were evaluated. The interim analysis of the post-market study TriT-aAY_Ali_Refr in the context of CMRP showed a low mean residual refractive astigmatism of -0.17 ± 0.27 D, ranging from 0 to -1.00 D, at 1 month postoperatively. In the post-market study TDI-aA_Hal_VA [28] the ToricaDiff-aA achieved a significant reduction in refractive astigmatism at 6 months postoperatively compared to preoperative values, and residual refractive astigmatism was ≤ 0.25 D in 91% of eyes. In the post-market study Tor-aA_Lin_Rota [29, 30] the mean (absolute) astigmatism was reduced from $1.62 \text{ D} \pm 0.39 \text{ D}$ (range 1.02 to 2.56; topographic measurement) preoperatively to $0.56 \text{ D} \pm 0.34 \text{ D}$ (range 1.25 to 0.00; refractive measurement) 6 months after surgery. Results were similar 1.5 years postoperatively. Vector analysis showed a significant reduction in magnitude of the J0 vector after surgery ($p < 0.0001$) with a mean absolute change in astigmatism of 0.77 ± 0.35 D. The mean J45 vector did not change significantly from preoperative values, which can be explained by very low preoperative values (only 11 out of 40 eyes had oblique astigmatism). The interim analysis of the company-initiated study Tor-aA_Ros_Rota [31] showed a reduction in the magnitude of the J0 and the J45 vectors from pre- to postoperatively. In addition, the correction of astigmatism was similar irrespective of whether the IOL was implanted in eyes with an AL < 25.0 mm or an AL ≥ 25.0 mm. In a case report, Kaindlstorfer et al., 2018 [21] reported a reduction of the refractive astigmatism from 3.75 D to 0.25 D one week after implantation of the Torica-aA. In conclusion, the above findings show that pre-existing regular, stable corneal astigmatism can be corrected with the toric lenses of the aA IOL platform as shown with the TrivaT-aA/-aAY, Torica-aA/-aAY and the ToricaDiff-aA as intended by the manufacturer.

The scientific publications regarding Diff-aA/-aAY add the following information. Two of the published studies are RCTs [23, 25] allowing a direct comparison between the Diff-aA and well-established bifocal aberration-correcting, hydrophobic “gold standards” competitor models (AcrySof ReSTOR +3.0 and Tecnis ZLB00, respectively). Grabner [25] reported the Diff-aA to be at least on par with the AcrySof ReSTOR +3.0 MIOL in terms of visual acuity and reading performance and provided functional reading performance even under low illumination and contrast conditions. The Diff-aA outperformed the competitor model in terms of near VA and reading performance under low luminance and contrast levels. Miháltz and Vécsei-Marlovits [23] showed that distance and near VA achieved with the Diff-aA and with the Tecnis ZLB00 were very similar. In terms of visual quality based on wavefront aberrometry, the spherical aberration-correcting Tecnis MIOL was found to be more sensitive to a large angle alpha than the aberration-free Diff-aA. In turn, image quality deteriorated more strongly with increasing pupil size in the case of the Diff-aA, which had to be expected since the aberration-free design of the Diff-aA does not compensate for the cornea’s spherical aberration [23]. All five publications confirmed a good visual performance of the Diff-aA [23-27]. When tested, the level of spectacle independence and patient satisfaction was found to be high [24-27].

In case of the Torica-aA, a case report [21] confirmed the highly satisfying results in terms of visual performance. Although focusing on a postoperative complication (which was caused by a wrong medical treatment), the early postoperative data demonstrated good visual acuity (UDVA: 0.0 logMAR, CDVA: -0.1 logMAR), good refraction (-2.50 D, as targeted) and a low residual refractive astigmatism (-0.25 D).

Taken together, these data demonstrate that the IOLs of the Multifocal Triva-aA Family achieve the performance intended by the manufacturer as required by the MDR Annex I, when used according to the intended use. Although clinical data for the TrivaT-aA/-aAY are still spare (only a small number of interim results), these data already confirm that the IOL performs as intended without any additional concerns.

SAFETY ANALYSIS

Several complications and undesirable effects have been reported in the literature for foldable posterior-chamber IOLs (including MIOLs) but the incidence is overall very small, particularly, if one takes into account that eligible patients often have several ocular comorbidities. In most cases these side-effects are related to the surgery and not device-related. Regarding the Diff-aA/-aAY, Triva-aA/-aAY and the parent models Torica-aA/-aAY and ToricaDiff-aA/-aAY, the reported outcomes show only a small percentage of complications and undesirable effects which have been already described in the literature.

Based on the results of the post-market studies and real-world data collections with the Multifocal Triva-aA Family and their multifocal optic parents (chapter 5.3.1), a safety analysis was performed. The primary safety endpoint was met as the cumulative and persistent adverse events (complication) rates were within limits compared to the historical SPE rates reported in ISO 11979-7:2018. **Key safety outcomes are presented below in Table 4.**

.

Table 4: Calculation of complication rates for the Multifocal Triva-aA Family and parent models according to DIN EN ISO 11979-7:2018-08 [32]

Endpoint	Acceptance criterion (SPE rate) %	Threshold rate % for n = 100	Threshold rate % for n = 300	Clinical outcomes % for n = 222	Data source
Cumulative					<u>Post-market studies Multifocal Triva-aA Family (n = 160 eyes):</u> <u>Diff-aA/ aAY (n=55)</u> Dif-aA_Gla_Seco Dif-aA_Gla_Comb <u>Triva-aA/-aAY (n=105)</u> Tri-aA_CM RP Tri-aAY_Pra_RWD Tri-aAY_Gla_VA
Cystoid macular oedema (CME)	3.0	8.9	6.0	0.90	
Hypopyon	0.3	3.0	1.8	0.0	
Endophthalmitis	0.1	3.0	1.0	0.0	
Lens dislocated from posterior chamber	0.1	3.0	1.0	0.0	
Pupillary block	0.1	3.0	1.0	0.0	
Retinal detachment	0.3	3.0	1.8	0.45*	
Secondary surgical intervention	0.8	4.2	2.6	0.45	
Persistent					<u>Parent models (n=62):</u> <u>Torica-aA/-aAY (n=40)</u> Tor-aA_Lin_Rota <u>ToricaDiff-aA/-aAY (n=22)</u> TDI-aA_Hal_VA
Corneal stroma oedema	0.3	3.0	1.8	0.0	
Cystoid macular oedema (CME)	0.5	4.2	2.2	0.0	
Iritis	0.3	3.0	1.8	0.0	
Raised IOP requiring treatment	0.4	4.2	1.8	0.0	

* value is not statistically significantly higher than SPE rate (threshold for a study with n=300 is not exceeded)

In terms of rotational stability, a pooled analysis including all members of the HumanOptics aA IOL platform has been performed and shows that the aA IOL platform meets the safety requirements as described in DIN EN ISO 11979-7:2018-08 (Table 5).

Table 5: Calculation of rotational stability according to DIN EN ISO 11979-7:2018-08

Endpoint	Acceptance criterion %	Clinical outcomes % for n = 127 analyzable eyes	Data source
Rotation < 10°	90	95.3	<u>Post-market studies (n = 127):</u> Asp-aA_Wie_Rota Tor-aA_Lin_Rota Tor-aA_Ros_Rota
Rotation < 20°	95	97.6	

Data for mid-term follow-up (1-6 months postoperatively) vs end of surgery (within 24h) are shown.

Undesirable effects commonly reported with Regarding MIOLs include **the perception of positive dysphotopsia** like glare, starbursts and halos, due to the superposition of focused and unfocused multiple images. In rare cases, symptoms may be so disturbing and persistent that the IOL needs to be explanted (and typically replaced with a monofocal IOL). Dysphotopsia symptoms were also observed in the post-market studies with the Diff-aA, the parent/related model ToricaDiff-aA, and in the literature (Dexl et al. [27] and Fan et al. [24]) but more than 90% rated the degree of impairment as low. In addition, the patients' level of satisfaction was high. For the trifocal models, Triva-aA/-aAY and TrivaT-aA/-aAY, symptoms of positive dysphotopsia may be even less pronounced due to the reduced number of diffractive rings compared to the bifocal ones (7 vs 9). In the real-world data collection with the Triva-aAY (Tri-aAY_Pra_RWD) mainly halos and starbursts were reported 6 months postoperatively. However, no patient considered the symptoms as "very bothersome". Glare occurrence was minimal (<10%).

A further prerequisite to achieve successful visual outcomes with MIOLs, is the **accurate alignment of the IOL**, because these lenses are more sensitive to decentration and tilt than monofocal IOLs [5]. However, a certain extent of tilt and decentration is common for IOLs (and the crystalline lens) and remains clinically unnoticed for any IOL design [33]. Regarding the Multifocal Triva-aA Family, the two real-world data collections Tri-aA_CMRP and Tri-aAY_Pra_RWD provide evidence of a good stability of the IOL platform in terms of tilt and decentration.

Regarding complaints related to complications, very few have been reported for the Multifocal Triva-aA Family and the parent models.

Clinical studies with HumanOptics IOLs have been conducted using both the Compact Line as well as the Safeloader-packaged IOLs and confirmed the performance and safety of the tested injectors. However, common injector-associated complications, such as trapped haptics or misconfiguration of the IOL haptics after injector release, have also been observed in a few cases. These were typically associated with handling errors.

Since 2003 HumanOptics has placed several million hydrophilic acrylic IOLs on the market. Over this period of more than 20 years no feedback has been received relating to material degradation other than very rare cases (< 5 ppm/year) of suspected IOL opacification (calcification). The latter is caused by crystallisation of calcium on the IOL surface, usually resulting from environmental circumstances, i.e. it does not constitute a degradation due to aging processes of IOLs manufactured by HumanOptics. The available information confirms that HumanOptics hydrophilic acrylic IOLs are stable over a period of at least 20 years (based on simulated aging tests according to ISO 11979-5). It can be expected that the device is stable indefinitely from the implantation date over the lifetime of the patient.

Taken together, these data demonstrate that the IOLs of the Multifocal Triva-aA Family fulfil the requirements regarding safety as set out by the MDR Annex I. The devices are safe for implantation in the capsular bag for the visual correction of aphakia (and corneal astigmatism), and no unanticipated risks associated with the IOLs were identified.

BENEFIT-RISK ANALYSIS

The main **benefits** of cataract surgery and IOL implantation in general are improved visual function and quality of life. Important aspects involved are:

- Improved vision (visual acuity, contrast sensitivity, depth perception, etc.)
- Reduction in falls, fractures and vehicle crashes

- Reduced eyeglass dependence (through correction of refractive errors and, in case of TIOs, of astigmatism)
- Improvement of quality of life (activity, confidence, mood)
- Reduction of anxiety and depression

In terms of improved vision, the different IOL types currently available on the European market can address single or a combination of vision problems. For instance, **monofocal IOLs** are designed to improve vision at just one set distance (generally far vision); **multifocal IOLs** are designed to additionally improve vision at near and/or intermediate distances and **toric IOLs** specifically correct pre-existing corneal astigmatism (they are available as monofocal toric or multifocal toric IOLs).

Three types of intended clinical benefits can be attributed to the Multifocal Triva-aA Family IOLs. The patient benefit **“restored vision”** is based on achieving good CDVA after IOL implantation, which has been demonstrated in the previous chapters.

When implanted with a MIO, the patient will also expect an improvement in near vision and, in the case of trifocal IOLs, in intermediate vision (intended benefit: **“vision at multiple distances”**). This is confirmed by the good DCNVA and DCIVA values which have been shown in the previous chapters. Importantly, also a high degree of spectacle independence was reached: In the collections of real-world data Tri-aA_CM RP and Tri-aA_Pra_RWD with Triva-aA/-aAY full spectacle independence (i.e. at all distances) was achieved in the majority of patients (range 62.5% to 100%). In those cases where patients reported that they “sometimes” still needed glasses, it was for reading very small letters (near vision). In the case of the bifocal Diff-aA (post-market studies Dif-aA_Gla_[Seco/Comb]), all patients were independent for reading glasses. In line with these results, the majority of patients (78%) reported being completely spectacle free after implantation of the parent model ToricaDiff-aA (TDI-aA_Hal_VA), while the remaining patients needed glasses “sometimes” or “very rarely” (11% each).

This high degree of spectacle independence was also reflected in terms of patient satisfaction: in the post-market studies and collection of real-world data where this was analyzed, all patients reported either to be “very satisfied”/“satisfied” with their vision or that they would choose the same lens model again (see chapter 5.3.1; studies Dif-aA_Gla_Comb, Dif-aA_Gla_Seco, Tri-aA_CM RP, Tor-aA_Lin_Rota and TDI-aA_Hal_VA). It should be highlighted that many of the patients included in the study TDI-aA_Hal_VA were relatively young and underwent refractive lens exchange (RLE). In the case of RLE, the expected benefits are lower than for cataract patients, while patient expectations are certainly higher. Therefore, the positive result (all patients “very satisfied”/“satisfied” and would have chosen the same lens model again) confirms that the HumanOptics MIO design is suitable for RLE patients.

In the case of TIOs (TrivaT-aA/-aAY with their parent models Torica-aA/-aAY and ToricaDiff-aA/-aAY), a reduced dependence on eyeglasses also requires an **“effective reduction of pre-existing astigmatism”**. In line with this, all patients bilaterally implanted with the monofocal Torica-aA achieved complete spectacle independence for the targeted vision distance (far or near) as reported in the post-market study Tor-aA_Lin_Rota. Regarding spectacle independence with the ToricaDiff-aA the majority of patients (78%) reported being completely spectacle free after implantation (see the previous paragraph).

Regarding **risks**, very few complications have been reported for the Multifocal Triva-aA Family, which is in line with the low frequency of complications reported for cataract surgery in general. Nonetheless, complications can always occur but only in extremely rare cases complications

may lead to serious and permanent impairment of the patient. Compared to monofocal lenses, MIOs carry an elevated risk of dysphotopsias / photic phenomena and failed neuroadaptation, which in rare cases may require an IOL exchange. In the case of the Multifocal Triva-aA Family, the reported dysphotopsias such as glare, halos and starbursts have been found to cause mainly mild symptoms. In addition, the relative low number of diffractive rings (7) of HumanOptics trifocal IOLs may help to reduce the occurrence of symptoms. Some of the most relevant risks with IOL implantation in general are ocular inflammations and cystoid macular edema.

As mentioned earlier, postoperative complications are rarely device-related but rather depend on the preoperative ocular comorbidity, which means that IOLs should be implanted into risky patients only after careful assessment of the benefit/risk ratio. Due to the optical effects of MIOs, patient selection and counseling are particularly important. Furthermore, the success of cataract surgery is influenced by the surgical technique, equipment used, and the surgeon's experience. In addition, there is a learning curve when performing implantations with a IOLs with novel designs. Nevertheless, the implantation of IOLs has to be considered as a routine procedure for the target user group, i.e. ophthalmic surgeons. Since the essential aspects of the handling and implantation process are adequately described in the IFU, a special training for the Multifocal Triva-aA Family is not required. However, device handling errors cannot be completely excluded but appear to occur very rarely. A risk that specifically affects TIOs compared to Non-TIOs is postoperative rotation which might sometimes necessitate a secondary procedure to realign the lens. Currently, the respective data for the devices in question compared favourably to the thresholds set by the applicable standards.

In summary, a large body of evidence attests that the Diff-aA/-aAY and Triva-aA/-aAY are reliable products and that the remaining risks are well-controlled. First results with the new TrivaT-aA/-aAY show good agreement with that of its parent models.

A positive benefits-versus-risks profile can be assigned to the medical devices under consideration Diff-aA/-aAY, Triva-aA/-aAY and TrivaT-aA/-aAY, particularly considering the current knowledge and state of the art in the medical fields concerned. There is sufficient evidence that the intended clinical benefits are achieved. The probable benefits outweigh the probable risks for the intended use of the device.

5.5. Ongoing or planned post-market clinical follow-up

Clinical evidence for the safety and performance of the Multifocal Triva-aA Family has been based on the manufacturer's own clinical data and literature results, showing that implantation of the IOLs in the capsular bag is a safe treatment for the correction of aphakia and pre-existing astigmatism for the toric version, providing good and stable visual performance as well as high levels of patient satisfaction.

Furthermore, the devices (except for the recently launched TrivaT-aA/-aAY) are well established in the medical field of cataract surgery, and have been successfully marketed for several years without apparent risks which would have arisen in the course of the PMS activities of the manufacturer. Currently, two post-market studies are running, one with the Triva-aAY and one with the TrivaT-aAY, with the goal to provide supplemental data for the performance of the devices.

6. Possible diagnostic or therapeutic alternatives

Implantation of an IOL is the only well-tolerated treatment available to replace the focusing power of the natural lens after surgical extraction. The preferred placement of an IOL is in the capsular bag. Mild/early cataracts may also be managed by prescribing new eyeglasses to correct a lens-induced change in refractive errors. High-power eye glasses (still in use in

developing countries) and contact lenses (only for high myopic eyes) are alternative treatments for aphakia but do not achieve visual results as good as IOLs.

For presbyopic patients without cataract, a refractive lens exchange (RLE) can be considered, i.e removal of the clear crystalline lens with replacement with a presbyopia-correcting IOL. Other alternatives include progressive eye glasses, multifocal contact lenses, eye drops or to undergo a “PresbyLASIK”, in which an excimer laser is used to sculpt a multi-focal correction into the surface of the eye.

Unlike any of those alternatives, RLE obviates the need for future cataract surgery.

7. Suggested profile and training for users

IOLs must be handled by health professionals and implanted by ophthalmic surgeons. A high level of surgical skill is required for IOL implantation. The surgeon should have observed and/or assisted in a large number of implantations and successfully completed one or more courses on IOL implantation before attempting a surgery. It is also strongly recommended that the surgeon and his/her team receive a training on the specific product and accessories by the respective sales representative before the first implantation of the device. Application guides (PDF) may be available on request. Product application videos can be found on www.humanoptics.com in the “application” section for the respective device.

If you have further questions on the usage of the device, please contact your local sales representative or the manufacturer (application@humanoptics.com).

8. Reference to any harmonised standards and CS applied

Common specification(s) to comply with, if applicable:
N/A
Harmonised standards to apply, if applicable:
DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
DIN EN ISO 14971:2020-07 Medical devices – Application of risk management to medical devices
Relevant adopted monographs of the European Pharmacopoeia
Not applicable

9. References

- [1] Summary of Safety and Clinical Performance (SSCP) for Monofocal Aspira-aA Family (V1.0 Draft / 2023-11)
- [2] Schrecker J, Feith A, Langenbucher A. Comparison of additional pseudophakic multifocal lenses and multifocal intraocular lens in the capsular bag. *Br J Ophthalmol*. 2014;98(7):915-919.
- [3] Schrecker J, Langenbucher A. Visual Performance in the Long Term With Secondary Add-on Versus Primary Capsular Bag Multifocal Intraocular Lenses. *J Refract Surg*. 2016;32(11):742-747.
- [4] Holladay JT, Piers PA, Koranyi G, van der Mooren M, Norrby NE. A new intraocular lens design to reduce spherical aberration of pseudophakic eyes. *J Refract Surg*. 2002;18(6):683-691.
- [5] Sieburth R, Chen M. Intraocular lens correction of presbyopia. *Taiwan J Ophthalmol*. 2019;9(1):4-17.
- [6] Kimura S, Morizane Y, Shiode Y, Hirano M, Doi S, Toshima S, et al. Assessment of tilt and decentration of crystalline lens and intraocular lens relative to the corneal topographic axis using anterior segment optical coherence tomography. *PLoS One*. 2017;12(9):e0184066.
- [7] Schrecker J, Blass S, Langenbucher A. Silicone-diffractive versus acrylic-refractive supplementary IOLs: visual performance and manual handling. *J Refract Surg*. 2014;30(1):41-48.
- [8] Rajabi MT, Korouji S, Farjadnia M, Naderan M, Rajabi MB, Khosravi B, et al. Higher order aberration comparison between two aspherical intraocular lenses: MC6125AS and Akreos advanced optics. *Int J Ophthalmol*. 2015;8(3):565-568.
- [9] Liekfeld A, Ehmer A, Schröter U. Visual function and reading speed after bilateral implantation of 2 types of diffractive multifocal intraocular lenses: Add-on versus capsular bag design. *J Cataract Refract Surg*. 2015;41(10):2107-2114.
- [10] Chang JS, Ng JC, Chan VK, Law AK. Cataract Surgery with a New Fluidics Control Phacoemulsification System in Nanophthalmic Eyes. *Case Rep Ophthalmol*. 2016;7(3):218-226.
- [11] Lasta M, Miháltz K, Kovács I, Vécsei-Marlovits PV. Effect of Spherical Aberration on the Optical Quality after Implantation of Two Different Aspherical Intraocular Lenses. *J Ophthalmol*. 2017;2017:8039719.
- [12] Liu S-Y, Yang X. Visual outcome with the Oculentis Mplus intraocular lens. *Int Eye Sci*. 2017;17(11):2134-2136.
- [13] Naujokaitis T, Scharf D, Baur I, Khoramnia R, Auffarth GU. Bilateral implantation of +56 and +58 diopter custom-made intraocular lenses in patient with extreme nanophthalmos. *Am J Ophthalmol Case Rep*. 2020;20:100963.
- [14] Mayer CS, Hoffmann AM, Prahs P, Reznicek L, Khoramnia R. Functional outcomes after combined iris and intraocular lens implantation in various iris and lens defects. *BMC Ophthalmol*. 2020;20(1):370.
- [15] Sun Y, Zhang Z, Wei Y, Chi W, Zhang S. Intraocular Lens Fixation Technique Without Corneal Incision in Minimally Invasive Vitrectomized Eyes. *Ophthalmol Ther*. 2022;11(2):729-737.
- [16] Yang J, Hai-Bo C, Ming-Bing Z, Jun Y. Error analysis of visual quality in segmented refractive multifocal intraocular lenses implanted eyes. *International Eye Science*. 2022:447-451.
- [17] Mayer CS, Baur ID, Storr J, Khoramnia R. Traumatische Aniridie: Konservativer oder chirurgischer Therapieansatz? *Ophthalmologe*. 2022;119(4):395-399.
- [18] Chuprov AD, Startseva MI, Original, Чупров А, Старцева М. Comparative analysis of clinical and economic efficiency in implantation of various IOL models; Original Сравнительный анализ клинико-экономической эффективности при имплантации различных моделей ИОЛ. *Journal "Medicine"; Original Журнал «Медицина»* №. 2022;2:1.

- [19] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. *J Cataract Refract Surg.* 2022;48(1):75-82.
- [20] Liu X, Shao Y, Lin H, Liu C, Shen J, Zhang L, et al. Corneal densitometry: an innovative method to quantitatively evaluate corneal changes after phacovitrectomy. *BMC Ophthalmol.* 2023;23.
- [21] Kaindlstorfer C, Kneifl M, Reinelt P, Schönherr U. Rotation of a toric intraocular lens from neodymium:YAG laser posterior capsulotomy. *J Cataract Refract Surg.* 2018;44(4):510-511.
- [22] Shukhaiev SV, Ulitina AY, Shukhaeva EN. Comparative evaluation of rotational stability of hydrophobic (AcrySof IQ Toric) and hydrophilic (Torica-aAY) toric IOLs (Russian full-text; English abstract and translation available). *Point of View East – West.* 2019;1:16-19.
- [23] Miháltz K, Vécsei-Marlovits PV. The impact of visual axis position on the optical quality after implantation of multifocal intraocular lenses with different asphericity values. *Graefes Arch Clin Exp Ophthalmol.* 2021;259(3):673-683.
- [24] Fan F, Liu X, Y. F, Luo Y. A pilot study of +3.5 D aspheric diffractive multifocal intraocular lens implantation using real-time intraoperative navigation system. *Chin J Optom Ophthalmol Vis Sci.* 2017;19(2):87-92.
- [25] Grabner G. Best kept secrets Diffractiva Diff-aA. *Cataract & Refractive Surgery Today Europe.* 2017(Jan):52-53.
- [26] Vo Tan P. Multifocal implants and high myopia: my personal experience with the Diffractiva®. Presentation at DIFFRACTIVA symposium, Guangzhou, Mar 31, 2015. 2015.
- [27] Dexl AK, Zaluski S, Rasp M, Grabner G. Visual performance after bilateral implantation of a new diffractive aspheric multifocal intraocular lens with a 3.5 D addition. *Eur J Ophthalmol.* 2014;24(1):35-43.
- [28] Final Clinical Report ToricaDiff-aA (Halle) - Clinical evaluation of a toric multifocal intraocular lens - ToricaDiff-aA (2020-01)
- [29] Gyöngyössy B, Jirak P, Schönherr U. Rotational stability and patient satisfaction after implantation of a new toric IOL. *Eur J Ophthalmol.* 2016;26(4):321-327.
- [30] Gyöngyössy B, Jirak P, Schönherr U. Long-term rotational stability and visual outcomes of a single-piece hydrophilic acrylic toric IOL: a 1.5-year follow-up. *Int J Ophthalmol.* 2017;10(4):573-578.
- [31] Rusch W, Anwar M, Altas C, Steirat M, Fuchsluger T, Walckling M, et al. Influence Of Axial Length On Rotational Stability And Other Clinical Outcomes After Implantation Of A Hydrophilic Toric Intraocular Lens. 41th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; Session: IOLs and Evaluation Testing). 2023.
- [32] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia
- [33] Ashena Z, Maqsood S, Ahmed SN, Nanavaty MA. Effect of Intraocular Lens Tilt and Decentration on Visual Acuity, Dysphotopsia and Wavefront Aberrations. *Vision (Basel).* 2020;4(3).

Part 2 – SSCP intended for patients

This Summary of Safety and Clinical Performance (SSCP) is intended to provide public access to an updated summary of the main aspects of the safety and clinical performance of the device. The information presented below is intended for patients or lay persons. A more extensive summary of its safety and clinical performance prepared for healthcare professionals is found in the first part of this document.

The SSCP is not intended to give general advice on the treatment of a medical condition. Please contact your ophthalmologist in case you have questions about your medical condition or about the use of the device in your situation. This SSCP is not intended to replace an Implant card or the Instructions For Use to provide information on the safe use of the device.

This SSCP is reviewed at least annually and will be updated if needed.


1. Device identification and general information

Device trade name(s)	• Diff-aA (trade name: Diffractiva) • Diff-aAY (trade name: Diffractiva) • Triva-aA • Triva-aAY • TrivaT-aA • TrivaT-aAY (collectively referred to as “Multifocal Triva-aA Family”)
Manufacturer’s name and address	HumanOptics Holding AG Spardorfer Str. 150 91054 Erlangen Germany
Basic UDI-DI(s)	• Diff-aA/-aAY, Triva-aA/-aAY: 04049154_PC_M2_H2_O4_C2 • TrivaT-aA/-aAY: 04049154_PC_M2_H2_O5_C5
Year when the first certificate (CE) was issued covering the device	• Diff-aA/-aAY: 2010 • Triva-aA/-aAY: 2020 • TrivaT-aA/-aAY: 2021

2. Intended use of the device

Intended purpose	The device is a tiny intraocular lens (IOL) intended to be implanted in the capsular bag ² of your eye to replace the natural lens (crystalline lens).




² **Capsular bag:** The capsule which the natural lens sits in.

Indication(s) and intended patient groups	The IOL is used to replace your natural lens after removal by an ophthalmic surgeon. In many cases, the removal of your natural lens is because your lens has become cloudy with the age. This is called “cataract”. Cataract affects your vision.  The device is designed for patients who desire near, intermediate and distance vision and who want to be less dependent on glasses. Therefore, even if your natural lens is clear, the IOL may be used to correct presbyopia. Presbyopia is the gradual loss of your eyes' ability to focus on nearby objects, which affects middle-aged and older adults. Implantation of an IOL is indicated for adult patients only. Note: there are no clinical data available for children, pregnant or breastfeeding women, and immunocompromised persons in connection with the implantation of intraocular lenses from HumanOptics. Depending on your eye anatomy, life-style and preferences, different types of IOLs can be offered to you by your ophthalmologist (see section 3 “mode of action/operation principle”).
Contraindications	There are general contraindications for ocular surgery but no additional contraindications which are specific to the implantation of HumanOptics IOLs.



3. Device description

Description of the device	The Diff-aA/-aAY, Triva-aA/-aAY and TrivaT-aA/-aAY are tiny artificial lenses made of plastic, called intraocular lenses (IOLs), which are implanted into the eye to replace the natural lens. The devices are multifocal³ IOLs (also called presbyopia-correcting IOLs) designed to improve vision at several distances. While the Diff-aA/-aAY is designed to improve near and far vision, the Triva-aA/-aAY can also improve intermediate vision. The TrivaT-aA/-aAY is designed like the Triva-aA/-aAY but can additionally reduce or eliminate corneal astigmatism (see section 3 “mode of action/operation principle”). In many cases, the IOL will be used to replace the clouded lens (cataract) of a person. If lens opacity increases, surgery is the only treatment available to restore good vision. Currently, no alternative or drug therapies are available. However, multifocal IOLs can also be used to decrease spectacle dependence in presbyopic patients with a clear lens (see section 2). The IOLs have a diameter of just a few millimeters and fit perfectly into the eye due to their shape and the flexibility of their material.

³ **Multifocality:** Normal (monofocal) IOLs produce a crisp image on the retina only for objects which have a fixed distance from the eye. Multifocal IOLs are able to image objects from several distances.

	IOLs are composed of an “optic” (the round part of the IOL) and the “haptics”. The optic has a dioptric power that permits the focusing of light onto the retina and the haptics which are attached to the edge of the optic will help maintain the IOL in place in the eye. The IOLs are available as colourless or yellow versions. The yellow variant (with “-Y”) absorbs a portion of the blue light, which may help protect the retina⁴ over the long term.  Schematic frontal and lateral view of the Diff-aA (left) and TrivaT-aA (right) with its toric markings. The IOLs exist also in yellow (Diff-aAY, Triva-aAY and TrivaT-aAY). The dimensions are given in mm.
Mode of action / operating principle	If you have been diagnosed with a cataract, it is because your natural lens has become cloudy and this affects your vision. Most commonly, cataract occurs in older people, and usually develops gradually. Cataract removal is one of the most common operations worldwide. During the procedure, usually an outpatient operation, the natural eye lens is removed and replaced by an IOL. The short procedure is usually painless and causes very little discomfort for the patient (see also chapter 5.1). Step 1 – You have been diagnosed with a cataract  Step 2 – The cataract is break it into tiny pieces using ultrasonic waves or a laser, which are then suctioned out of the eye with a vacuum. The procedure takes usually only a few minutes and is performed under local anesthesia 

⁴ **Retina:** The innermost, light-sensitive layer of tissue of the eye that transmits that information to the brain via the optic nerve

	Step 3 – After removal of the cataract, an IOL is inserted with an injector into the capsular bag of your eye to restore your vision  Step 4 – The IOL is now in place in the capsular bag of your eye  Note: the above surgical procedure is also followed in the case of a clear lens exchange (i.e. in presbyopic patients without a cataract). There is a variety of IOL types to choose from — all with different benefits. Choosing the right IOL can be difficult and it depends on several factors including your eye condition and your lifestyle. The different options should be discussed with your ophthalmologist. Multifocal IOLs are designed to provide clear vision at several distances (far and near for the Diff-aA/-aAY and additionally intermediate distance for Triva-aA/-aAY and TrivaT-aA/-aAY). Therefore, it is expected that you will be less dependent on spectacles. However, glasses or contact lenses may still be needed, especially in the case of the Diff-aA/-aAY for activities at intermediate distances (like computer work). Besides your cataract, you might also suffer from corneal astigmatism. Astigmatism is not a disease and is relatively common: about 20% of cataract patients have moderate to severe astigmatism. Corneal astigmatism happens when your cornea is irregularly shaped. To make it simple, an eye without astigmatism is round like a football ball and with astigmatism, it is more like a rugby ball. The shape of the cornea with astigmatism makes light bend differently as it enters your eye, causing a refractive error, resulting in blurred or distorted vision. A good option to reduce or eliminate your corneal astigmatism is a toric IOL such as the TrivaT-aA/-aAY. Vision without astigmatism

	 Simulated visual impression with astigmatism 
Materials or substances in contact with the patient's tissues	All IOLs manufactured by HumanOptics are made of biocompatible plastic (hydrophilic acrylic material) with an integrated UV absorber. Models with names ending with "Y" additionally contain a yellow coloring agent intended to filter the blue light.
Information about medicinal substances in the device, if any	Not applicable.
Accessories	None which are relevant to patients

4. Risks and warnings

Contact your ophthalmologist if you believe that you are experiencing side effects related to the device or its use or if you are concerned about risks. This document is not intended to replace a consultation with your ophthalmologist if needed.

How potential risks have been controlled or managed	Potential risks have been identified and controlled according to international standards that have been developed by experts from many countries. The corresponding document(s) comprising rules, guidelines and processes are approved and published by a globally recognized organisation.
Remaining risks and undesirable effects	As with any surgical intervention, routine cataract surgery and associated IOL implantation are not entirely risk-free. However, this procedure is among the most common and most successful ones performed nowadays in the world. There are no known risks which are specific to the Diff-aA/-aAY, Triva-aA/-aAY and TrivaT-aA/-aAY and which have not been previously reported for similar IOLs. The following side effects may occur following the IOL implantation but they normally disappear shortly after the surgery: • Foreign body sensation caused by the small incision through which the lens was inserted • Halos, which can form around sources of light, particularly at dusk or in darkness • You may experience increased sensitivity to light in the first few weeks after surgery The following symptoms may indicate complications associated with an IOL implantation. If these occur, you should consult your ophthalmologist: • Sudden or gradual deterioration in vision / blurred vision • Severe redness of the eye • Eye pain Months or years after IOL implantation, your vision might deteriorate and become blurry or hazy. This might be due to posterior capsule opacification (PCO), which occurs when a cloudy layer of scar tissue forms behind the IOL. This condition is not dangerous and can be treated by performing an opening in the capsule (or bag) that contains your IOL. Opening the capsule allows light to enter the eye and improve the clarity of the vision, as long as the rest of the eye is healthy. This procedure is called "Nd:YAG laser capsulotomy". It is fast (a few minutes), painless and is an easy way for restoring your vision. There is no other treatment for PCO and the condition will not improve by itself.
Warnings and precautions	Patients with certain eye diseases and/or malformations have higher-than-normal risks to develop complications during and after IOL implantation. Before the surgery, your ophthalmologist will perform an assessment of your eye and will discuss with you anything to do with your operation. This includes the risks and benefits of surgery and how long it will take to fully recover after the surgery. Difficulties during surgery may also require choosing another lens than planned.
Summary of any field safety corrective action (FSCA including FSN), if applicable	No devices have been recalled since the Diff-aA/-aAY, Triva-aA/-aAY and TrivaT-aA/-aAY have been sold on the market. In July 2024, HumanOptics issued a Field Safety Notice (FSN) ⁵ for healthcare professionals. The notice affected all devices in a certain packaging which had been manufactured in the period

⁵ A field safety notice (FSN) is an important means of communicating safety-related information to users.

	January 1, 2020 to June 28, 2024. The manufacturer had observed very rare cases of damaged outer blister packaging after transportation. The associated hazard potential is considered to be very low and no patient injury has happened to the manufacturer's knowledge. The FSN served to inform the users of this possible risk and to tell them how a damaged packaging can be recognized. No specific action is needed in the case of patients who have been treated with potentially affected products. Meanwhile, the product packaging has been improved.

5. Summary of clinical evaluation and post-market clinical follow-up

5.1. Clinical background of the device

Cataract, i.e. the clouding of the natural lens of the eye, is a very common condition and most commonly affect adults as a result of ageing. Cataract causes deterioration of your vision and usually gets slowly worse over time. Surgery to replace the cloudy lens is the only way to improve your vision. Implantation of an IOL is the only well-tolerated treatment available to replace the cloudy lens. This procedure that usually takes 30 minutes is typically associated with rapid visual recovery, good visual outcomes and minimal complications.

After you and your ophthalmologist have decided that you should have your cataract removed, the length of your eye and the shape of your cornea will be measured. These measurements will help your surgeon choose the best IOL for you. You will also discuss the various lens options available to you. Further information about cataract surgery can be found under:

<https://www.humanoptics.com/en/patienten/help-with-cataracts/>

After implantation, the IOL will normally remain in your eye for the rest of your life, provided there is no medical contraindication. The product lifetime of HumanOptics IOLs is proven for twenty years. Regular ophthalmological check-ups are recommended, as for patients with natural lenses.

5.2. The clinical evidence for the CE-marking

The letters 'CE' which appear on products or packages signify that products sold in the European Economic Area (EEA) have been assessed to meet high safety, health, and environmental protection requirements.

Clinical evidence for CE-marking IOLs should meet the requirements defined by international standards and is based on clinical data with the device, scientific literature and feedback from the market.

The safety and good performance of the Multifocal Triva-aA Family has been shown *via* the results of several clinical studies and data collections of the devices under consideration and their parent models (Torica-aA/-aAY and ToricaDiff-aA/-aAY). This is further supported by several publications in scientific journals. The implantation of the IOLs in the capsular bag is a safe treatment to restore vision after cataract surgery resulting in high levels of patient satisfaction. Furthermore, the HumanOptics multifocal IOLs (Diff-aA/-aAY) are well established in the medical field of cataract surgery, and have been successfully marketed for more than 10 years without safety concerns.

In summary, the IOLs of the Multifocal Triva-aA Family perform well in terms of visual acuity, astigmatism correction for the toric model, patient satisfaction and safety. The clinical outcomes

are good and comparable to other modern IOLs on the market. The lenses show to be a safe and effective option for the surgical treatment of cataract.

5.3. Safety

The implantation of the CE-marked Diff-aA/-aAY, Triva-aA/-aAY and TrivaT-aA/-aAY IOLs considerably improves vision in cataract patients. The CE marking for IOLs is mandatory and means that the manufacturer has verified that the products meet European safety, health or environmental requirements.

The side effects of IOLs are well-known and are similar for all IOLs on the market. They generally occur rarely and are not very harmful (see chapter 4 “Remaining risks and undesirable effects”). Therefore, the benefit-risk ratio of the devices is very positive. The manufacturer continuously analyses complaints, other market feedback and new scientific publications. The results of these activities have been favorable, for the IOLs of the Multifocal Triva-aA Family and their parent models. This further supports the safety of the devices, besides the other data which were evaluated for the CE marking. In addition, further clinical data are currently being collected with the Triva-aAY to obtain additional evidence of the device's performance and demonstrate its long-term safety.

6. Possible diagnostic or therapeutic alternatives

When considering alternative treatments, it is recommended to contact your ophthalmologist who can take into account your individual situation.

Currently, there are no drugs which can stop or reverse the clouding of the natural lens of the eye. In the case of a pronounced cataract, the clouded lens must therefore be removed and replaced by an IOL. The preferred placement of an IOL is in the capsular bag. If you have a mild/early cataract, your ophthalmologist may prescribe only eyeglasses until your cataract becomes more mature.

If you are presbyopic and your natural lens is clear, alternatively to IOL implantation, you may use progressive eye glasses, multifocal contact lenses or eye drops. Another option may be a special type of corneal laser surgery to achieve vision at multiple distances (so-called “PresbyLASIK”). Unlike any of those alternatives, the IOL implantation obviates the need for future cataract surgery.

There are different types of IOLs on the market, each having specific advantages and disadvantages. As mentioned earlier, choosing the right IOL can be difficult and it depends on several factors including your eye condition and your lifestyle. The different options should be discussed with your ophthalmologist.

7. Suggested training for users

Not applicable to patients since the device is to be used only by healthcare professionals.

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Restrisiken und unerwünschte Wirkungen	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>multifocal</u> IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with <u>toric</u> IOLs:

	A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Vorsicht	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal</u> IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Warnhinweise	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Besondere Überlegungen vor der Implantation einer multifokalen IOL	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes.

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Остатъчни рискове и нежелани ефекти	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with multifocal IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time.

	Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Внимание	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal</u> IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Предупреждения	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.

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Zbytková rizika a nežádoucí účinky	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with multifocal IOLs: As with <u>all</u> multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Upozornění	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection

	• Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal</u> IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Upozornění	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Zvláštní úvahy před implantací multifokální nitrooční čočky	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.
Zvláštní úvahy před implantací torické nitrooční čočky	• The dioptric power of toric IOLs is indicated as "spherical equivalent (SE) and cylinder (cyl)". Please check the label carefully. • Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism. • Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation • If possible, posterior corneal astigmatism (tomography) should be considered, especially with astigmatism < 2.0 D.

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Υπολειπόμενοι κίνδυνοι και ανεπιθύμητες ενέργειες	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction

	Further points to note in connection with multifocal IOLs: As with <u>all</u> multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Προσοχή	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal <u>IOLs</u> caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Προειδοποιήσεις	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be

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Riesgos residuales y efectos indeseables	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with multifocal IOLs: As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time.

	Further points to note in connection with toric IOLs: A postoperative rotation_of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Precaución	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal IOLs caution should also <u>be exercised</u> in the following conditions: Amblyopia Optic nerve atrophy Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Advertencias	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Consideraciones especiales previas a la implantación de una LIO multifocal	• It is recommended to target emmetropia.

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Risques résiduels et effets indésirables	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment Cystoid macular edema Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with multifocal <u>IOLs</u>: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with toric IOLs:

	A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Mise en garde	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal IOLs caution <u>should also</u> be exercised in the following conditions: Amblyopia Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Avertissements	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Considérations particulières avant l'implantation d'une lentille intraoculaire multifocale	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes.

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Common specification(s) to comply with, if applicable:
N/A
Harmonised standards to apply, if applicable:
DIN EN ISO 13485:2021-12 Medical devices – Quality management systems – Requirements for regulatory purposes
DIN EN ISO 14971:2022-04 Medical devices – Application of risk management to medical devices
Relevant adopted monographs of the European Pharmacopoeia
Not applicable

9. Références

- [1] Summary of Safety and Clinical Performance (SSCP) for Monofocal Aspira-aA Family (V1.0 Draft / 2023-11)
- [2] Schrecker J, Feith A, Langenbucher A. Comparison of additional pseudophakic multifocal lenses and multifocal intraocular lens in the capsular bag. *Br J Ophthalmol*. 2014;98(7):915-919.
- [3] Schrecker J, Langenbucher A. Visual Performance in the Long Term With Secondary Add-on Versus Primary Capsular Bag Multifocal Intraocular Lenses. *J Refract Surg*. 2016;32(11):742-747.
- [4] Holladay JT, Piers PA, Koranyi G, van der Mooren M, Norrby NE. A new intraocular lens design to reduce spherical aberration of pseudophakic eyes. *J Refract Surg*. 2002;18(6):683-691.
- [5] Sieburth R, Chen M. Intraocular lens correction of presbyopia. *Taiwan J Ophthalmol*. 2019;9(1):4-17.
- [6] Kimura S, Morizane Y, Shiode Y, Hirano M, Doi S, Toshima S, et al. Assessment of tilt and decentration of crystalline lens and intraocular lens relative to the corneal topographic axis using anterior segment optical coherence tomography. *PLoS One*. 2017;12(9):e0184066.
- [7] Schrecker J, Blass S, Langenbucher A. Silicone-diffractive versus acrylic-refractive supplementary IOLs: visual performance and manual handling. *J Refract Surg*. 2014;30(1):41-48.
- [8] Rajabi MT, Korouji S, Farjadnia M, Naderan M, Rajabi MB, Khosravi B, et al. Higher order aberration comparison between two aspherical intraocular lenses: MC6125AS and Akreos advanced optics. *Int J Ophthalmol*. 2015;8(3):565-568.
- [9] Liekfeld A, Ehmer A, Schröter U. Visual function and reading speed after bilateral implantation of 2 types of diffractive multifocal intraocular lenses: Add-on versus capsular bag design. *J Cataract Refract Surg*. 2015;41(10):2107-2114.
- [10] Chang JS, Ng JC, Chan VK, Law AK. Cataract Surgery with a New Fluidics Control Phacoemulsification System in Nanophthalmic Eyes. *Case Rep Ophthalmol*. 2016;7(3):218-226.
- [11] Lasta M, Miháltz K, Kovács I, Vécsei-Marlovits PV. Effect of Spherical Aberration on the Optical Quality after Implantation of Two Different Aspherical Intraocular Lenses. *J Ophthalmol*. 2017;2017:8039719.
- [12] Liu S-Y, Yang X. Visual outcome with the Oculentis Mplus intraocular lens. *Int Eye Sci*. 2017;17(11):2134-2136.
- [13] Naujokaitis T, Scharf D, Baur I, Khoramnia R, Auffarth GU. Bilateral implantation of +56 and +58 diopter custom-made intraocular lenses in patient with extreme nanophthalmos. *Am J Ophthalmol Case Rep*. 2020;20:100963.
- [14] Mayer CS, Hoffmann AM, Prahs P, Reznicek L, Khoramnia R. Functional outcomes after combined iris and intraocular lens implantation in various iris and lens defects. *BMC Ophthalmol*. 2020;20(1):370.
- [15] Sun Y, Zhang Z, Wei Y, Chi W, Zhang S. Intraocular Lens Fixation Technique Without Corneal Incision in Minimally Invasive Vitrectomized Eyes. *Ophthalmol Ther*. 2022;11(2):729-737.
- [16] Yang J, Hai-Bo C, Ming-Bing Z, Jun Y. Error analysis of visual quality in segmented refractive multifocal intraocular lenses implanted eyes. *International Eye Science*. 2022:447-451.
- [17] Mayer CS, Baur ID, Storr J, Khoramnia R. Traumatische Aniridie: Konservativer oder chirurgischer Therapieansatz? *Ophthalmologe*. 2022;119(4):395-399.
- [18] Chuprov AD, Startseva MI, Original, Чупров А, Старцева М. Comparative analysis of clinical and economic efficiency in implantation of various IOL models; Original Сравнительный

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- [19] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. J Cataract Refract Surg. 2022;48(1):75-82.
- [20] Liu X, Shao Y, Lin H, Liu C, Shen J, Zhang L, et al. Corneal densitometry: an innovative method to quantitatively evaluate corneal changes after phacovitrectomy. BMC Ophthalmol. 2023;23.
- [21] Kaindlstorfer C, Kneifl M, Reinelt P, Schönherr U. Rotation of a toric intraocular lens from neodymium:YAG laser posterior capsulotomy. J Cataract Refract Surg. 2018;44(4):510-511.
- [22] Shukhaiev SV, Ulitina AY, Shukhaeva EN. Comparative evaluation of rotational stability of hydrophobic (AcrySof IQ Toric) and hydrophilic (Torica-aAY) toric IOLs (Russian full-text; English abstract and translation available). Point of View East – West. 2019;1:16-19.
- [23] Miháltz K, Vécsei-Marlovits PV. The impact of visual axis position on the optical quality after implantation of multifocal intraocular lenses with different asphericity values. Graefes Arch Clin Exp Ophthalmol. 2021;259(3):673-683.
- [24] Fan F, Liu X, Y. F, Luo Y. A pilot study of +3.5 D aspheric diffractive multifocal intraocular lens implantation using real-time intraoperative navigation system. Chin J Optom Ophthalmol Vis Sci. 2017;19(2):87-92.
- [25] Grabner G. Best kept secrets Diffractiva Diff-aA. Cataract & Refractive Surgery Today Europe. 2017(Jan):52-53.
- [26] Vo Tan P. Multifocal implants and high myopia: my personal experience with the Diffractiva®. Presentation at DIFFRACTIVA symposium, Guangzhou, Mar 31, 2015. 2015.
- [27] Dexl AK, Zaluski S, Rasp M, Grabner G. Visual performance after bilateral implantation of a new diffractive aspheric multifocal intraocular lens with a 3.5 D addition. Eur J Ophthalmol. 2014;24(1):35-43.
- [28] Final Clinical Report ToricaDiff-aA (Halle) - Clinical evaluation of a toric multifocal intraocular lens - ToricaDiff-aA (2020-01)
- [29] Gyöngyössy B, Jirak P, Schönherr U. Rotational stability and patient satisfaction after implantation of a new toric IOL. Eur J Ophthalmol. 2016;26(4):321-327.
- [30] Gyöngyössy B, Jirak P, Schönherr U. Long-term rotational stability and visual outcomes of a single-piece hydrophilic acrylic toric IOL: a 1.5-year follow-up. Int J Ophthalmol. 2017;10(4):573-578.
- [31] Rusch W, Anwar M, Altas C, Steirat M, Fuchsluger T, Walckling M, et al. Influence Of Axial Length On Rotational Stability And Other Clinical Outcomes After Implantation Of A Hydrophilic Toric Intraocular Lens. 41th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; Session: IOLs and Evaluation Testing). 2023.
- [32] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia
- [33] Ashena Z, Maqsood S, Ahmed SN, Nanavaty MA. Effect of Intraocular Lens Tilt and Decentration on Visual Acuity, Dysphotopsia and Wavefront Aberrations. Vision (Basel). 2020;4(3).

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Manufacturer's name and address	HumanOptics Holding AG Spardorfer StraÙe 150 91054 Erlangen Germany
Manufacturer's single registration number (SRN)	DE-MF-000017892
Authorised representative if applicable; name and SRN	N/A
Notified body name and organization number	TÜV NORD CERT GmbH NB 0044


Trgovački naziv(i) uređaja	• Diff-aA (trg. naziv: Diffractiva) • Diff-aAY (trg. naziv: Diffractiva) • Triva-aA • Triva-aAY • TrivaT-aA • TrivaT-aAY (zajednički naziv "Multifokalna obitelj Triva-aA")
Osnovni UDI-DI(ji)	• Diff-aA/-aAY, Triva-aA/-aAY: 04049154_PC_M2_H2_O4_C2 • TrivaT-aA/-aAY: 04049154_PC_M2_H2_O5_C5
CND/EMDN kod(ovi):	• Diff-aA/-aAY, Triva-aA/-aAY: P030102100202 (IOL-ovi, aphakic, multifokalni, asferični, hidrofilni akrilni) • TrivaT-aA/-aAY: P030102100302 (IOLs, APHAKIC, MULTIFOCAL, TORIC, HYDROPHILIC ACRYLIC)
Klasa uređaja	☐ I ☐ IIa ☒ IIb ☐ III
Godina kada je izdan prvi certifikat (CE) za uređaj	• Diff-aA/-aAY: 2010 • Triva-aA/-aAY: 2020 • TrivaT-aA/-aAY: 2021

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4. Rizici i upozorenja

Preostali rizici i neželjeni učinci	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with multifocal IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Oprez	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection

	• Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal</u> IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Upozorenja	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Posebne smjernice prije implantacije multifokalne IOL-a	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.
Posebne smjernice prije implantacije torične IOL-a	• The dioptric power of toric IOLs is indicated as “spherical equivalent (SE) and cylinder (cyl)”. Please check the label carefully. • Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism. • Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation • If possible, posterior corneal astigmatism (tomography) should be considered, especially with astigmatism < 2.0 D.

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Maradék kockázatok és nemkívánatos hatások	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further <u>points</u> to note in connection with <u>multifocal</u> IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further <u>points</u> to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Figyelem	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy

	• Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Figyelmeztetések	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Különleges szempontok a multifokális IOL beültetése előtt	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.
Különleges szempontok a torikus IOL beültetése előtt	• The dioptric power of toric IOLs is indicated as “spherical equivalent (SE) and cylinder (cyl)”. Please check the label carefully. • Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism. • Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation

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- [19] Bonsemeyer MK, Becker E, Liekfeld A. Dysphotopsiae and functional quality of vision after implantation of an intraocular lens with a 7.0 mm optic and plate haptic design. J Cataract Refract Surg. 2022;48(1):75-82.
- [20] Liu X, Shao Y, Lin H, Liu C, Shen J, Zhang L, et al. Corneal densitometry: an innovative method to quantitatively evaluate corneal changes after phacovitrectomy. BMC Ophthalmol. 2023;23.
- [21] Kaindlstorfer C, Kneifl M, Reinelt P, Schonherr U. Rotation of a toric intraocular lens from neodymium:YAG laser posterior capsulotomy. J Cataract Refract Surg. 2018;44(4):510-511.
- [22] Shukhaiev SV, Ulitina AY, Shukhaeva EN. Comparative evaluation of rotational stability of hydrophobic (AcrySof IQ Toric) and hydrophilic (Torica-aAY) toric IOLs (Russian full-text; English abstract and translation available). Point of View East – West. 2019;1:16-19.
- [23] Miháltz K, Vécsei-Marlovits PV. The impact of visual axis position on the optical quality after implantation of multifocal intraocular lenses with different asphericity values. Graefes Arch Clin Exp Ophthalmol. 2021;259(3):673-683.
- [24] Fan F, Liu X, Y. F, Luo Y. A pilot study of +3.5 D aspheric diffractive multifocal intraocular lens implantation using real-time intraoperative navigation system. Chin J Optom Ophthalmol Vis Sci. 2017;19(2):87-92.
- [25] Grabner G. Best kept secrets Diffractiva Diff-aA. Cataract & Refractive Surgery Today Europe. 2017(Jan):52-53.
- [26] Vo Tan P. Multifocal implants and high myopia: my personal experience with the Diffractiva®. Presentation at DIFFRACTIVA symposium, Guangzhou, Mar 31, 2015. 2015.
- [27] Dexl AK, Zaluski S, Rasp M, Grabner G. Visual performance after bilateral implantation of a new diffractive aspheric multifocal intraocular lens with a 3.5 D addition. Eur J Ophthalmol. 2014;24(1):35-43.
- [28] Final Clinical Report ToricaDiff-aA (Halle) - Clinical evaluation of a toric multifocal intraocular lens - ToricaDiff-aA (2020-01)
- [29] Gyöngyössy B, Jirak P, Schönherr U. Rotational stability and patient satisfaction after implantation of a new toric IOL. Eur J Ophthalmol. 2016;26(4):321-327.
- [30] Gyöngyössy B, Jirak P, Schönherr U. Long-term rotational stability and visual outcomes of a single-piece hydrophilic acrylic toric IOL: a 1.5-year follow-up. Int J Ophthalmol. 2017;10(4):573-578.
- [31] Rusch W, Anwar M, Altas C, Steirat M, Fuchsluger T, Walckling M, et al. Influence Of Axial Length On Rotational Stability And Other Clinical Outcomes After Implantation Of A Hydrophilic Toric Intraocular Lens. 41th Congress of the the European Society of Cataract and Refractive Surgery (ESCRS; Session: IOLs and Evaluation Testing). 2023.
- [32] DIN EN ISO 11979-7:2018-08: Ophthalmic implants – Intraocular lenses – Part 7: Clinical investigations of intraocular lenses for the correction of aphakia
- [33] Ashena Z, Maqsood S, Ahmed SN, Nanavaty MA. Effect of Intraocular Lens Tilt and Decentration on Visual Acuity, Dysphotopsia and Wavefront Aberrations. Vision (Basel). 2020;4(3).

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Rischi residui ed effetti indesiderati	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with multifocal IOLs: As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning

	should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Attenzione	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal IOLs <u>caution</u> should also be exercised in the following conditions: Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Avvertenze	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Considerazioni particolari prima dell'impianto di una IOL multifocale	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended.

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Liekančios rizikos ir nepageidaujami reiškiniai	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection <u>with</u> multifocal IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection <u>with</u> toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Atsargiai	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection

	• Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal</u> IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Įspėjimai	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Ypatingi aspektai, į kuriuos reikia atsižvelgti prieš implantuojant daugiažidžius IOL	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.
Ypatingi aspektai, į kuriuos reikia atsižvelgti prieš implantuojant torines IOL	• The dioptric power of toric IOLs is indicated as “spherical equivalent (SE) and cylinder (cyl)”. Please check the label carefully. • Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism. • Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation • If possible, posterior corneal astigmatism (tomography) should be considered, especially with astigmatism < 2.0 D.

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4. Risico's en waarschuwingen

Resterende risico's en ongewenste effecten	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in <u>connection</u> with multifocal IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in <u>connection</u> with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).

Waarschuwing	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal IOLs</u> caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Waarschuwingen	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Speciale overwegingen voorafgaand aan de implantatie van multifocale IOL	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.

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Ryzyko resztkowe i działania niepożądane	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction

	Further points to note in connection with multifocal IOLs: As <u>with</u> all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with toric IOLs: A postoperative <u>rotation</u> of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Ostrzeżenie	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal <u>IOLs</u> caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Ostrzeżenia	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be

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4. Riscos e avisos

Riscos residuais e efeitos indesejáveis	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with multifocal IOLs: As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Precaução	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis • Proliferative diabetic retinopathy

	• Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal IOLs caution <u>should</u> also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Advertências	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Considerações especiais antes da implantação de LIO multifocais	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.
Considerações especiais antes da implantação de LIO tóricas	• The dioptric power of toric IOLs is indicated as “spherical equivalent (SE) and cylinder (cyl)”. Please check the label carefully.

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Riscuri reziduale și efecte nedorite	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: Corneal endothelial decompensation Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with multifocal IOLs: As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with toric IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).

Atenție	Specific caution is required before implanting a lens in a patient with one or more of these conditions: Uveitis Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection • Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For multifocal IOLs caution <u>should also</u> be exercised in the following conditions: Amblyopia Optic nerve atrophy Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Avertismente	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Considerații speciale înainte de implantarea unei lentile intraoculare multifocale	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.

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Kvarstående risker och biverkningar	As with any surgical procedure, cataract surgery with IOL implantation presents risks, which the surgeon must evaluate. Potential complications and undesirable effects in relation to the cataract surgery and lens implantation may include but are not limited to: • Corneal endothelial decompensation • Corneal edema • Retinal detachment • Cystoid macular edema • Increased intraocular pressure • Inflammation (e.g. toxic anterior segment syndrome, endophthalmitis, uveitis, iritis) • Iris trauma • Posterior capsular and zonular rupture • Secondary surgical intervention (e.g. repositioning, removal, or exchange) • Posterior capsule opacification (PCO) • IOL decentration or tilt • Deviation from target refraction Further points to note in connection with <u>multifocal</u> IOLs: • As with all multifocal IOLs, a reduction in contrast sensitivity as compared to a monofocal IOL may occur. This can be more prevalent in low light conditions. Therefore, patients should be cautious when driving at night or in poor visibility conditions. • Some visual effects may be experienced due to the superposition of focused and unfocused multiple images. These may include some perceptions of halos or rings around point sources of light under dark conditions. The perception of the visual phenomena usually decreases with time. Further points to note in connection with <u>toric</u> IOLs: A postoperative rotation of the toric IOL may reduce the astigmatism correction. Misalignment of the toric IOL relative to the intended axis can reduce or negate the refractive benefit, or even worsen the refraction. In such cases, lens repositioning should be considered and should take place within the first two weeks after surgery, prior to IOL encapsulation (shrinking).
Varning	Specific caution is required before implanting a lens in a patient with one or more of these conditions: • Uveitis • Proliferative diabetic retinopathy • Uncontrolled chronic glaucoma • Corneal endothelial dystrophy • Microphthalmos or macrophthalmos • Suspected ocular infection

	• Extreme myopia (ocular axial length > 30 mm) may be a risk factor for IOL capsule inadhesion that may impact negatively visual acuity • Pre-existing ocular conditions which may negatively impact the stability of the implanted IOL (e.g. distorted eye due to previous trauma or developmental disorder, instability of the posterior capsule/zonules) • Surgical difficulties at the time of cataract extraction which might increase the risk for complications (e.g. persistent bleeding, significant iris damage, uncontrolled IOP elevation, significant vitreous prolapse or loss) • Circumstances that would result in damaging the endothelium during implantation For <u>multifocal</u> IOLs caution should also be exercised in the following conditions: • Amblyopia • Optic nerve atrophy • Significant corneal aberrations • Medical conditions that may impair vision Eye anatomies or conditions (such as refractive treatments) that do not allow reliable prediction of postoperative refraction
Varningar	• Due to the hydrophilic properties of the material, substances such as disinfectants, antibiotics or viscoelastics can be absorbed by the lens. This may lead to a toxic lens syndrome. Therefore, at the end of surgery, care should be taken to remove all substances from the eye by using standard irrigation/aspiration techniques. Also be aware that dyes used intraoperatively (e.g. trypan blue) may cause staining of the IOL. • Salt precipitation in/on the lens may occur in very rare cases when using exogenous material such as, but not limited to, add-on lens, air or gas during corneal surgery or vitrectomy. The mechanism and the incidence are to date not clear.
Särskilda överväganden före implantation av multifokala IOL	• It is recommended to target emmetropia. • Patients with significant preoperative astigmatism, determined by keratometry, or expected postoperative astigmatism > 0.50 D may not achieve optimal visual outcomes. • In cases of corneal astigmatism > 1.0 D, implantation of a multifocal toric IOL is recommended. • Lens tilt and decentration may negatively affect the quality of vision.
Särskilda överväganden före implantation av torisk IOL	• The dioptric power of toric IOLs is indicated as "spherical equivalent (SE) and cylinder (cyl)". Please check the label carefully. • Precise biometry, keratometry, topography/tomography and precise lens axis alignment relative to the intended axis are keys to a successful correction of astigmatism. • Eyes with an axial length > 24 mm are more at risk of post-surgery IOL rotation • If possible, posterior corneal astigmatism (tomography) should be considered, especially with astigmatism < 2.0 D.

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