

MODEL (IOL)	DIOPTRER RANGE	INJECTOR SYSTEM/CARTRIDGE		APPROX. EFFECTIVE INCISION SIZE	
ASPIRA-aA(Y)/MC 6125 AS(-Y)	-10.0 – 25.0 D	Medicel AccuJect™ 1.8-1P	LP 604540	2.2 mm	
	-10.0 – 40.0 D	Medicel AccuJect™ 2.0-1P	LP 604510	2.4 mm	
	-10.0 – 60.0 D	Medicel AccuJect™ 2.2-1P	LP 604530	2.5 mm	
	-10.0 – 60.0 D	Medicel ViscoJect™ 2.2-1P BIO	LP 604340C	2.5 mm	
	-10.0 – 30.0 D	R.E.T. Comport	C22R1010	2.7 mm	
DIFF-aA(Y)/MC 6125 Diff (DAY)	10.0 – 30.0 D	See Aspira-aA(Y)			
TORICA-aA(Y)/MC 6125 T(-Y)	sph 10.0 – 30.0 D/cyl 1.0 – 6.0 D	Medicel AccuJect™ 2.0-1P	LP 604510	2.4 mm	
	sph 10.0 – 30.0 D/cyl 1.0 – 6.0 D	Medicel ViscoJect™ 2.2-1P BIO	LP 604340C	2.5 mm	
TORICADIFF-aA(Y)	sph 10.0 – 30.0 D/cyl 1.0 – 6.0 D	Medicel AccuJect™ 2.0-1P	LP 604510 Rev.2v5/6	2.4 mm	
SAFELOADER®	<ul style="list-style-type: none"> <li>▪ ASPIRA-aA(Y)</li> </ul>	Medicel AccuJect™ 1.8-1P	LP 604540 Rev.2v5/6	2.2 mm	
		Medicel AccuJect™ 2.0-1P	LP 604510 Rev.2v5/6	2.4 mm	
		Medicel AccuJect™ 2.2-1P	LP 604530 Rev.2v5/6	2.5 mm	
	<ul style="list-style-type: none"> <li>▪ DIFF-aA(Y)</li> </ul>	See Aspira-aA(Y)			
	<ul style="list-style-type: none"> <li>▪ TORICADIFF-aA(Y)</li> </ul>	Medicel AccuJect™ 2.0-1P	LP 604510 Rev.2v5/6	2.4 mm	
MC X11 ASP	-6.0 – 25.0 D	Medicel AccuJect™ 1.8-1P	LP 604540	2.2 mm	
	-6.0 – 40.0 D	Medicel AccuJect™ 2.0-1P	LP 604510	2.4 mm	
	-6.0 – 40.0 D	Medicel AccuJect™ 2.2-1P	LP 604530	2.5 mm	
	-6.0 – 40.0 D	Medicel ViscoJect™ 2.2-1P BIO	LP 604340 C	2.5 mm	
ASPIRA-aQA/MC 6105	0.0 – 20.0 D	Medicel AccuJect™ 1.8-1P	LP 604540	2.2 mm	
	0.0 – 25.0 D	Medicel AccuJect™ 2.0-1P	LP 604510	2.4 mm	
	0.0 – 30.0 D	Medicel AccuJect™ 2.2-1P	LP 604530	2.5 mm	
	0.0 – 30.0 D	Medicel ViscoJect™ 2.2-1P BIO	LP 604340 C	2.5 mm	
	0.0 – 25.0 D	R.E.T. Comport	C22R1010	2.7 mm	
AS/MC 5812 AS	0.0 – 25.0 D	Medicel AccuJect™ 2.0-1P	LP 604510	2.4 mm	
	0.0 – 30.0 D	Medicel AccuJect™ 2.2-1P	LP 604530	2.5 mm	
	0.0 – 30.0 D	Medicel ViscoJect™ 2.2-1P BIO	LP 604340 C	2.5 mm	
	0.0 – 30.0 D	R.E.T. Comport	C22R1010	2.7 mm	
ASPIRA®3P-aVA	10.0 – 25.0 D	R.E.T. Comport P	P28P0020	3.0 mm	
	10.0 – 30.0 D	Medicel AccuJect™ 3.0-3P	LP 604580	3.2 mm	

MODEL (OTHER IMPLANTS)	SIZE	INJECTOR SYSTEM/CARTRIDGE		APPROX. EFFECTIVE INCISION SIZE
Capsular Tension Ring	CTR 13/11	Geuder Injector Caps. Tension Ring (G-32960)*		1.9 mm
ARTIFICIAL/IRIS with Fiber	10.5 – 9.5 mm	Medicel AccuJect™ 2.0-1P	LP 604510	2.4 mm
	12.0 – 9.5 mm	Medicel AccuJect™ 2.2-1P	LP 604530	2.5 mm
	11.5 – 9.5 mm	Medicel ViscoJect™ 2.2-1P BIO	LP 604340 C	2.5 mm
ARTIFICIAL/IRIS Fiber Free	12.8 – 9.5 mm	Medicel AccuJect™ 1.8-1P	LP 604540	2.2 mm
	12.8 – 9.5 mm	Medicel AccuJect™ 2.0-1P	LP 604510	2.4 mm
	12.8 – 9.5 mm	Medicel AccuJect™ 2.2-1P	LP 604530	2.5 mm
	12.8 – 9.5 mm	Medicel ViscoJect™ 2.2-1P BIO	LP 604340 C	2.5 mm
	12.8 – 9.5 mm	AMO UNFOLDER® Silver T*		2.8 mm

- SL = Safeloader
- SE = Spherical Equivalent=sphere + ½ cylinder
- sph = Sphere
- cyl = Cylinder

\* not available from HumanOptics AG