


ASPIRA-aAY+

TECHNICAL INFORMATION

MONOFOCAL
MONOBLOC
ACRYLIC
ADJUSTED ASPHERIC
PRELOADED



Type

- foldable posterior chamber IOL
- one-piece with C-loop haptics
- also available in preloaded Safeloader autoloading container

Material

- hydrophilic acrylic with UV absorber
- water content 26 %
- blue-light protection
- Abbe number 56
- refractive index 1.46

Specifics

- aspheric surfaces (adjusted aspheric)
- posterior surface with 360° lens epithelial cell barrier

Diopter Range

0.0 to 30.0 D

Sphere in 0.5 D steps
(optic shape)

10.0 to 30.0 D (biconvex)



Sphere in 1.0 D steps
(optic shape)

0.0 to 9.0 D (concave-convex)



Est. A-const

For a complete list of recommended constants,
please visit www.humanoptics.com



SL = Safeloader | CL = Compact Line

CE0044

 HumanOptics Holding AG

Spardorfer Str. 150 · 91054 Erlangen · Germany · Phone +49 (0) 9131 50 66 5-0 · Fax +49 (0) 9131 50 66 5-90
mail@humanoptics.com · www.humanoptics.com